
**SYMBOL GLOSSARY**

Symbol	Description
	Manufacturer/Hersteller/Fabricante/Fabricant/Fabbricante/ Fabricante/Üretici
	Authorized representative in the European Community/European Union/Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/ Representante autorizado en la Comunidad Europea/ Unión Europea/Mandataire dans la Communauté européenne/Union européenne/Mandatario nella Comunità Europea/Unione Europea/Mandatário na Comunidade Europeia/União Europeia/Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi
	Date of manufacture/Herstellungsdatum/Fecha de fabricación/Date de fabrication/Data di fabbricazione/ Data de fabrico/Üretim tarihi
	Use-by date/Haltbarkeitsdatum/Fecha de caducidad/ Date limite d'utilisation/Data di scadenza/Data de validade/Son kullanma tarihi
	Batch code/Chargenbezeichnung/Código de lote/ Code de lot/Codice lotto/Código de lote/Parti kodu
	Catalog number/Bestellnummer/Número de catálogo/ Número de reference/Numero di catalogo/Número de catálogo/Katalog numarası
	Sterilized using ethylene oxide/Sterilisiert mit Ethylenoxid/ Esterilizado con óxido de etileno/Stérilisé à l'oxyde d'éthylène/Sterilizzato con ossido di etilene/Esterilizado por óxido de etileno/Etilen oksit kullanılarak sterilize edilmiştir
	Do not use if the package is damaged and consult instructions for use/Nicht verwenden, wenn Verpackung beschädigt ist, und Gebrauchsanweisung beachten/No utilizar si el paquete está dañado y consultar las instrucciones de uso/Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation/Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso/Não utilizar se existirem danos na embalagem e consultar as instruções de utilização/ Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
	Do not re-sterilize/Nicht reesterilisieren/No reesterilizar/ Ne pas restériliser/Non risterilizzare/Não reesterilizar/ Tekrar sterilize etmeyin
	Do not reuse/Nicht wiederverwenden/No reutilizar/Ne pas réutiliser/Non riutilizzare/Não reutilizar/Tekrar kullanmayın
	Double sterile barrier system/Doppel-Sterilbarrieresystem/ Sistema de barrera estéril doble/Système de double barrière sterile/Sistema a doppia barriera sterile/Sistema de barreira estéril dupla/Çift steril bariyer sistemi
	Fragile, handle with care/Zerbrechlich, sorgfältig handhaben/Frágil, manipular con cuidado/Fragile, manipuler avec soin/Fragile, maneggiare con cura/ Frágil, manusear com cuidado/Kırılğan, dikkatli taşıyın
	Keep dry/Trocken aufbewahren/Mantener seco/Tenir au sec/Mantenere asciutto/Manter seco/Kuru tutun
	Temperature limit Upper = 30 °C, Lower = 15 °C Temperaturgrenze Obere = 30 °C, untere = 15 °C Límites de temperatura Superior = 30 °C, Inferior = 15 °C Limites de température Supérieure = 30 °C, inférieure = 15 °C Limite di temperatura Superiore = 30 °C, inferiore = 15 °C Limite de temperatura Máximo = 30 °C, Mínimo = 15 °C Sıcaklık sınırı Üst = 30 °C, Alt = 15 °C
	Humidity limitation Upper = 60%, Lower = 20% Feuchtigkeitsbeschränkung Obere = 60 %, untere = 20 %

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**Table of Contents**

1. SCOPE.....	2
2. IMPORTANT NOTES.....	2
3. DESCRIPTION.....	2
4. WARNINGS.....	2
5. PRECAUTIONS.....	2
6. DEVICE APPLICATION.....	3
7. PERFORMANCE CHARACTERISTICS.....	3
8. COMPATIBLE DEVICES/ACCESSORIES.....	4
9. INSTRUCTIONS FOR USE.....	4
10. MRI SAFETY INFORMATION.....	5
11. PACKAGING, STERILIZATION, AND REUSE.....	6
12. STORAGE CONDITIONS.....	6
13. DISPOSAL.....	6
14. INCIDENT REPORTING.....	6
15. STIMULATION CHARACTERISTICS .....	43-47
16. ADMINISTRATIVE.....	47

**1 SCOPE**

These instructions for use apply to the following devices:

- AD-TECH® Spencer Depth Electrodes AD Style, S Type Contact
- AD-TECH® Spencer Depth Electrodes LD Style, S Type Contact
- AD-TECH® Spencer Depth Electrodes LD Style, A Type Contact
- AD-TECH® Spencer Depth Electrodes RD Style, S Type Contact
- AD-TECH® Spencer Depth Electrodes SD Style, S Type Contact
- AD-TECH® Spencer Depth Electrodes SD Style, A Type Contact
- AD-TECH® Monopolar Depth Electrode (MD)
- AD-TECH® Stay Flange

**2 IMPORTANT NOTES**

	Before using these products, carefully review these Instructions for Use. Keep it readily available for reference to the Stimulation Parameters.
	Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.


**3 DESCRIPTION**

The Depth Electrodes are slender, cylindrical tubes inserted into the brain. Depth Electrodes have conductive contacts at the distal end to monitor the brain's electrical signals or deliver electrical stimulation.

The contacts connect to internal wires that run throughout the device's body to the tail, where they attach to electrical contacts that allow connection to compatible cable systems.

When connected to compatible equipment, the electrodes serve as the patient-electrode interface for temporary stimulation and recording electrical signals at the subsurface of the brain.

The Depth Electrodes are available in different configurations, including the number of contacts, the contact length, and the contact spacing as specified by the physician. See Figure 1 and Figure 2 for details.

**4 WARNINGS**

	This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories. Incorrect use of this device may result in severe patient harm.
	Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.
	Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.
	The Depth Electrode is for surgical use only. Placement must be in a sterile surgical environment.

	Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.
	The Depth Electrodes shall be surgically removed.
	The Depth Electrode is for short-term use (<30 days) only. Failure to observe the use conditions may result in severe patient harm.
	Do not use the Depth Electrode if the use-by date has been exceeded. Using devices with an expired shelf life may result in severe patient harm.
	Do not use the Depth Electrodes if the packaging is damaged or unintentionally opened before use. Using a device from damaged packing may result in contamination and severe patient harm.
	The Depth Electrodes are for single use only. Do not reuse the device. Reusing the device may result in contamination and severe patient harm.
	Do not re-sterilize the Depth Electrodes. Re-sterilization may result in device failure.
	The Depth Electrode should be introduced easily and NOT FORCED. Forcing the Depth Electrode into position may result in severe patient harm.
	Always use appropriate strain relief to secure the Depth Electrode tail and connecting cables. Failure to use appropriate strain relief may result in dislodgement of the electrode and severe patient harm.
	Do not withdraw the wire carrier (Depth Electrode tail) into the incision. Drawing the carrier tail into the incision may contaminate the wound and result in severe patient harm.
	Failure to comply with the recommended stimulation parameter may result in severe patient harm.
	Disconnect all Depth Electrodes from all monitoring equipment during cardiac defibrillation. Failure to disconnect the Depth Electrodes during cardiac defibrillation may result in severe patient harm.
	Do not use RF surgical equipment on patients undergoing EEG monitoring or stimulation with the Depth Electrodes. Using RF surgical equipment during EEG monitoring may result in severe patient harm.
	The AD and LD Style Spencer Probe Depth Electrodes are MRI unsafe. Using the AD or LD Style Depth Electrodes in an MRI environment may lead to severe patient harm.
	The Monopolar Depth Electrodes are MRI unsafe. Using the Monopolar Depth Electrodes in an MRI environment may lead to severe patient harm.
	The SD and RD Style Spencer Probe Depth Electrodes are MRI unsafe. Using SD or RD Style Depth Electrodes in an MRI environment may lead to severe patient harm.

**5 PRECAUTIONS**

	Handle the Depth Electrode with care to prevent damage. Pulling, stretching, or twisting the Depth Electrode may result in signal loss.
	Ensure that the Depth Electrode contacts are kept free from moisture and contamination.
	Do not remove the stylet before placement of the Depth Electrode. Removing and replacing the stylet before placement may damage the device. The stylet ensures that the Depth Electrode remains rigid during the delivery procedure.
	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.

	<p>Upon completion of placement and tunneling, connect the electrode contacts outside the body to the electrode cable.</p> <p>At no time should the electrode contacts outside the body be exposed or unprotected.</p> <p>Connection to an electrode cable protects the electrode contacts outside the body.</p> <p>If the patient becomes mobile, the electrode contacts must remain connected to the electrode cable.</p>
	<p>Observe the recommended connection instructions to reduce Electromagnetic Interference (EMI).</p> <p>Failure to observe the recommended connection instructions may result in unusable EEG signals.</p>
	<p>Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.</p>
	<p>Dispose of the Depth Electrodes following the healthcare facility's biohazardous and pathological waste procedures.</p>
	<p>Dispose of any accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharp's disposal container.</p>
	<p>Federal Law (USA) restricts this device to sale by or on the order of a physician.</p>

## 6 DEVICE APPLICATION

### INTENDED PURPOSE/USE

The Depth Electrodes are for temporary (<30 days) use with recording, monitoring, and stimulation equipment for the recording, monitoring, and stimulation of electrical signals at the subsurface level of the brain. The recording of electrical activity supports the definition of the location of epileptogenic foci and brain mapping.

### INDICATIONS FOR USE (EU)

The Depth Electrodes are indicated for patients undergoing diagnostic electrical brain mapping as part of neurological investigations and treatment planning.

### CONTRAINDICATIONS

	<p>The Depth Electrodes are not for use on any patient whom the physician/surgeon considers at risk of infection.</p>
	<p>The Depth Electrodes are not intended for continuous stimulation. Stimulation should only be applied to support brain mapping.</p>

### PATIENT TARGET GROUP

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

### INTENDED USER

The Depth Electrodes are to be used by physicians/surgeons experienced in Depth Electrodes and the associated Placement Accessories.

### ENVIRONMENT AND FACILITIES

Depth Electrodes are for surgical use only.

Placement must be in a sterile surgical environment.

The Depth Electrodes shall be surgically removed.

Depth Electrodes are placed at the target location using an appropriate stereotactic guidance system and supporting equipment.

### EXPECTED CLINICAL BENEFIT

Brain mapping allows the physician to gather detailed information about the electrical characteristics of the brain to identify areas related to physiological symptoms, enabling them to develop appropriate patient-specific treatment plans.

### UNDESIRABLE SIDE EFFECTS

Use of the Depth Electrodes carries the following clinical risks:

- Abscess
- Bleeding
- CSF leakage/drainage/pseudomeningocele
- Death
- Epidural hemorrhage
- Intracerebral hemorrhage
- Medical complications including deep vein thrombosis, pulmonary embolisms, allergic reactions, and psychiatric changes.
- Meningitis
- Permanent neurologic deficit
- Subarachnoid hemorrhage
- Subdural hemorrhage
- Superficial wound infection
- Transient neurologic deficit
- Worsening seizures

For further details, refer to the Summary of Safety and Clinical Performance (SSCP).

### SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) Unique ID: **SSCP-00002**

The SSCP is available from:

- EUDAMED under SRN: US-MF-00004795 and the Unique ID
- Email request to [sales@adtechmedical.com](mailto:sales@adtechmedical.com), quoting the SSCP Unique ID

### ELECTRONIC INSTRUCTIONS FOR USE

An electronic copy of these instructions for use is available from:

- Ad-Tech's website [www.adtechmedical.com](http://www.adtechmedical.com)
- Email request to [sales@adtechmedical.com](mailto:sales@adtechmedical.com), quoting the IFU document number AT10019-B

## 7 PERFORMANCE CHARACTERISTICS

### DEVICE ATTRIBUTES

Spencer Probe Depth Electrodes				
Style:	AD	LD	RD	SD
Attribute				
Diameter	1.96 mm	1.3 mm	0.86 mm	1.12 mm
Strain Relief (Stay Flange included)	No	No	No	Yes
Contact Material	Platinum	Platinum	Platinum	Platinum
No. of Contacts	2 – 16	2 – 16	2 – 16	2 – 16
Contact Length Standard (S)	1.27 mm	1.57 mm	2.29 mm	2.41 mm
Contact Diameter Standard (S)	1.91 mm	1.28 mm	0.86 mm	1.12 mm
Contact Length Alternative (A)	Not Available	5.33 mm	Not Available	1.32 mm
Contact Diameter Alternative (A)	Not Available	1.30 mm	Not Available	1.10 mm
Recording Area	26.27 – 36.27 mm	42.07 – 66.57 mm	27.29 – 77.29 mm	3.52 – 112.41 mm
Stylet Standard	Stainless Steel	Stainless Steel	Tungsten	Stainless Steel
Stylet Alternative	Not Available	Not Available	Not Available	Tungsten

Monopolar Depth Electrodes	
Diameter	1.12 mm
Strain Relief	No
Stylet Cover	No
Contact Material	Platinum
No. of Contacts	1
Contact Length Standard (S)	2.41 mm
Contact Diameter Standard (S)	1.10 mm
Recording Area	2.41 mm
Stylet Standard	Stainless Steel

### MONITORING CHARACTERISTICS

The Depth Electrodes connect to EEG monitoring equipment capable of processing signals greater than 32 Hz with a minimum input voltage of 50  $\mu$ V.

### STIMULATION CHARACTERISTICS

	<p>Failure to comply with the recommended stimulation parameter may result in severe patient harm.</p>
---	--

The Depth Electrodes connect to EEG stimulation equipment with a constant-current generator of less than or equal to 34 mA and a constant-voltage stimulator of less than or equal to 10 V.

## Instructions for Use Spencer Probe and Monopolar Depth Electrodes

Charge density is the most valuable and accurate measure of the intensity of tissue stimulation by an electrode, considering current strength, pulse duration, and electrode contact size. AD-TECH® has calculated the charge density using the formula:

$$\text{Charge Density} = \frac{\mu\text{C}}{\text{cm}^2} = \frac{\text{Stimulation Intensity (mA)} \times \text{Pulse Duration (mS)}}{\text{Electrode Surface Area (cm}^2\text{)}}$$

Refer to the Stimulation Parameters tables for the safe tissue stimulation parameters on pages 43 - 47.

### 8 COMPATIBLE DEVICES/ACCESSORIES

#### STAY FLANGE

The AD-TECH® Stay Flange is an optional accessory for the Depth Electrodes. AD-TECH® supplies each SD Style Spencer Probe Depth Electrode with a Stay Flange. The Stay Flange provides strain relief to the Depth Electrode tail, anchoring it via sutures to the patient's skull.

AD-TECH® does not supply the Stay Flange with the AD, LD, or RD Style Spencer Probe Depth Electrodes or the Monopolar Depth Electrode. However, they are available separately using part number FLANGE-X.

#### ANCHOR BOLTS

The AD-TECH® Anchor Bolts are optional accessories for use with LD, RD, and SD Style Spencer Probe Depth Electrodes. Anchor Bolts are secured in the skull to provide an access point. The Anchor Bolts minimize concerns about potential cerebrospinal fluid (CSF) leakage and infection of the subdural space while stabilizing the electrode. The LD, RD, and SD Style Depth Electrodes are compatible with the following Anchor Bolts:

AD Style Spencer Probe Depth Electrodes and Monopolar Depth Electrodes are **NOT compatible** with the Anchor Bolts.

#### SLOTTED CANNULAS

The AD-TECH® Slotted Cannulas are optional accessories for use with the LD, RD, and SD Style Spencer Probe Depth Electrodes. Slotted Cannulas create a path within the brain to facilitate the Depth Electrode placement. A Slotted Cannula is applied to place a Depth Electrode through a craniotomy or a burr hole without using an Anchor Bolt. The following slotted cannulas are available:

AD Style Spencer Probe Depth Electrodes and Monopolar Depth Electrodes are **NOT compatible** with the Slotted Cannulas.

#### OBTURATORS

The AD-TECH® Obturators are optional accessories for use with RD, SD, and LD Style Spencer Probe Depth Electrodes. Obturators create a path within the brain to facilitate the placement of Depth Electrodes. The Obturator is applied when placing Depth Electrodes through a burr hole with an applied Anchor Bolt. The following Obturators are available:

#### DISPOSABLE DRILL KITS

The AD-TECH® Disposable Drill Kits are an optional accessory to the Spencer Probe Depth Electrodes and Monopolar Depth Electrodes. When used with a compatible Surgical Drill, the Drill Kit makes a burr hole to accept the Anchor Bolt and Depth Electrode. The following Drill Kits are available:

#### RULER

The AD-TECH® Ruler is an optional accessory to the Spencer Probe Depth Electrodes and Monopolar Depth Electrodes. The Ruler is used to make critical measurements when placing the Depth Electrodes.

#### CONNECTION SYSTEMS

AD-TECH® provides cables in several standard configurations to suit the user's needs.

Contact AD-TECH® sales support to specify the compatible cable configuration fully.

The following connection cable models are available:

Catalog Number	Description	Compatible Electrodes	Cable Configuration (y)
L-DC-yDIN-aa*	TECH-ATTACH Connection System Cable	Spencer Probe Depth Electrodes, Style AD, LD, RD, SD	2-64
L-DCL-yDIN-aa*	TECH-ATTACH Connection System Cable	Spencer Probe Depth Electrodes, Style AD, LD, RD, SD	2-64
L-SRL-yDIN-aa*	Light Weight CABRIO Cables	Spencer Probe Depth Electrodes, Style AD, LD, RD, SD	2-64

Catalog Number	Description	Compatible Electrodes	Cable Configuration (y)
711-LWKDINX	Lead Wire Set	Monopolar Depth Electrode	N/A
711-LWKDINX-09	Lead Wire Set	Monopolar Depth Electrode	N/A

\*Note: 'aa' denotes the standard configurations available from AD-TECH®. Contact an AD-TECH® sales representative for support.

The following connection blocks are available for use with the TECH-ATTACH Connection System Cable:

Catalog Number	Description	Compatible Cables
DC-yX	y-contact connection block	Lightweight TECH-ATTACH Connection System Cable
DC-ySX	y-contact slotted connection block	Lightweight TECH-ATTACH Connection System Cable

### 9 INSTRUCTIONS FOR USE

	<i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories. Incorrect use of this device may result in severe patient harm.</i>
	<i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i>
	<i>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</i>
	<i>The Depth Electrode is for surgical use only. Placement must be in a sterile surgical environment. Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i>
	<i>The Depth Electrodes shall be surgically removed.</i>
	<i>The Depth Electrode is for short-term use (&lt;30 days) only. Failure to observe the use conditions may result in severe patient harm.</i>
	<i>The Depth Electrode should be introduced easily and NOT FORCED. Forcing the Depth Electrode into position may result in severe patient harm.</i>
	<i>Always use appropriate strain relief to secure the Depth Electrode tail and cables. Failure to use appropriate strain relief may result in dislodgement of the electrode and severe patient harm.</i>
	<i>Do not withdraw the wire carrier (Depth Electrode tail) into the incision. Drawing the carrier tail into the incision may contaminate the wound and result in severe patient harm.</i>
	<i>Disconnect all Depth Electrodes from all monitoring equipment during cardiac defibrillation. Failure to disconnect the Depth Electrodes during cardiac defibrillation may result in severe patient harm.</i>

	Do not use RF surgical equipment on patients undergoing EEG monitoring or stimulation with the Depth Electrodes. Using RF surgical equipment during EEG monitoring may result in severe patient harm.
	Handle the Depth Electrode with care to prevent damage. Pulling, stretching, or twisting the Depth Electrode may result in signal loss.
	Ensure that the Depth Electrode contacts are kept free from moisture and contamination.
	Do not remove the stylet before placement of the Depth Electrode. Removing and replacing the stylet before placement may damage the device. The stylet ensures that the Depth Electrode remains rigid during the delivery procedure.
	Upon completion of placement and tunneling, connect the electrode contacts outside the body to the electrode cable. At no time should the electrode contacts outside the body be exposed or unprotected. Connection to an electrode cable protects the electrode contacts outside the body. If the patient becomes mobile, the electrode contacts must remain connected to the electrode cable.
	Observe the recommended connection instructions to reduce Electromagnetic Interference (EMI). Failure to observe the recommended connection instructions may result in unusable EEG signals.
	Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.

**ELECTRODE PLACEMENT**

The physician can place the Depth Electrode through a standard twist drill or burr hole, or craniotomy site.

Prepare the path through the brain tissue using the chosen method, Obturator or Slotted Cannula.

Regarding Figure 3, the steps for placing the Depth Electrode are as follows:

**Steps 1 and 2:** Slide the Depth Electrode along the pre-planned/prepared trajectory to the target location. Withdraw the stylet while maintaining the Depth Electrode position.

**Step 3:** Create a strain relief S-shape in the tail of the Depth Electrode. Open the Stay Flange, c.

**Step 4:** Clamp the Stay Flange around the body of the Depth Electrode in the center of the strain relief S-shape.

**Step 5:** Suture the Stay Flange in positions a and b to the patient's scalp, securing the tail of the Depth Electrode.

When using the Depth Electrode with the Anchor Bolt, follow the placement instructions in the Anchor Bolt instructions for use.

When using the Depth Electrode in combination with the Slotted Cannula, follow the placement instructions in the Slotted Cannula instructions for use.

When using the Depth Electrode with the Obturator, follow the instructions in the Obturator instructions for use.

**CONNECTING THE ELECTRODES**

Once the Depth Electrode is secured, immediately connect the exposed contacts to the compatible AD-TECH® Connection System.

To connect the EEG cable, see the CABRIO®/TECH-ATTACH® Connection Systems instructions for use. Also, consult the electrode code chart provided with the electrodes.

When connecting the electrodes to the EEG monitoring/stimulation equipment, consider the following best practices to reduce the likelihood of Electromagnetic Interference (EMI) and unwanted EEG signal interference/artifacts:

- Ensure that the third-party EEG equipment meets international safety standards, including those for electrical safety and Electromagnetic Compatibility (EMC).
- Connect the EEG equipment following the original equipment manufacturer's instructions for use, including observing any electrical grounding requirements.
- Ensure that the EEG amplifier provides EMI filtering and has a Common Mode Rejection Ratio (CMRR) of at least 80 dB at 50/60 Hz.
- Reduce any sources of potential EMI by switching off any unnecessary equipment, disconnecting it from the mains network, and removing mobile phones or other telecommunication equipment from the vicinity.
- Utilize a patient ground reference connection.
- Ensure the Depth Electrode external (proximal) contacts are kept clean and dry to avoid high impedance connections or crosstalk.
- Ensure the Connector Blocks are closed tightly around the Depth Electrode proximal contacts to avoid high impedance or intermittent connections.

- Bundle/braid all recording electrode cables together to help reduce 50 Hz/60 Hz noise generated by the mains network.
- Avoid running the connection cables near or over potential sources of EMI, such as other electrical equipment or connecting cables. If it is necessary to cross cables, ensure they cross perpendicular to each other to limit cross-coupling of EMI.
- Monitor the EEG signals for signs of noise and other signal artifacts that are not being generated by the patient's brain activity and take appropriate action to address the root causes.

Ensure the EEG monitor displays the expected output signals on all electrode contacts to confirm Depth Electrode operation.

During the monitoring period, ensure that the Depth Electrodes are not subject to excessive tensile force, particularly during a patient's seizure.

Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.

**TROUBLESHOOTING**


- In the event of signal loss, perform the following troubleshooting steps:
1. Confirm the correct connection of the electrode tails with the cable connection system.
  2. Verify proper DIN end connection to the third-party EEG system.
  3. Confirm the Depth Electrode tail is not under tension and that there are no bends smaller than 1.5 mm in diameter.
  4. Confirm there is a proper patient ground reference.
  5. For signal artifact issues, consult the EEG equipment manufacturer's troubleshooting guide.

**DEVICE REMOVAL**

Remove the devices within the specified lifetime of <30 days.

Remove the Depth Electrodes surgically within a sterile surgical environment.

**Figure 3 – Using the Depth Electrode**

**10 MRI SAFETY INFORMATION**

	The AD and LD Style Spencer Probe Depth Electrodes are MRI unsafe. Using the AD or LD Style Depth Electrodes in an MRI environment may lead to severe patient harm.
	The Monopolar Depth Electrodes are MRI unsafe. Using the Monopolar Depth Electrodes in an MRI environment may lead to severe patient harm.
	The RD and SD Style Spencer Probe Depth Electrodes are MRI unsafe. Using the RD or SD Style Depth Electrodes in an MRI environment may lead to severe patient harm.

**MRI SAFETY, AD, AND LD STYLE SPENCER PROBE DEPTH ELECTRODES**

AD-TECH® has not verified the safety and performance of the AD and LD Style Spencer Probe Depth Electrodes in the MRI environment.

AD-TECH® has classified these Depth Electrodes as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice for Marking Medical Devices and Other Items for Safety in The Magnetic Resonance Environment defines MRI Unsafe as:

*A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.*

AD-TECH® has labeled the AD and LD Style Spencer Probe Depth Electrodes and the Monopolar Depth Electrodes as MRI Unsafe using the symbol:

**MRI SAFETY, MONOPOLAR DEPTH ELECTRODES**

AD-TECH® has not verified the safety and performance of the Monopolar Depth Electrodes in the MRI environment. AD-TECH® has classified these Depth Electrodes as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice for Marking Medical Devices and Other Items for Safety in The Magnetic Resonance Environment defines MRI Unsafe as:

*A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.*

AD-TECH® has labeled the Monopolar Depth Electrodes as MRI Unsafe using the symbol:

**MRI SAFETY, RD, AND SD-STYLE SPENCER PROBE DEPTH ELECTRODES**

In the EU, AD-TECH® has not verified the conditional use of the RD and SD Style Depth Electrodes and Anchor Bolts in 1.5 T and 3 T MRI scanners. AD-TECH® has classified these Depth Electrodes as MRI Unsafe and labelled them appropriately.

ASTM F2503-20 Standard Practice for Marking Medical Devices and Other Items for Safety in The Magnetic Resonance Environment defines MRI Unsafe as:

*A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.*

AD-TECH® has labeled the RD and SD Style Depth Electrodes as MRI Unsafe using the symbol:

**11 PACKAGING, STERILIZATION, AND REUSE**

Observe the following packaging, sterilization, and reuse instructions:

	AD-TECH® sterilizes the Depth Electrodes using ethylene oxide.
	Do not use the Depth Electrode if the use-by date has been exceeded. Using devices with an expired shelf life may result in severe patient harm.
	Do not use the Depth Electrodes if the packaging is damaged or unintentionally opened before use. Using a device from damaged packing may result in contamination and severe patient harm.
	The Depth Electrode is for single use only. Do not reuse the device. Reusing the device may result in contamination and severe harm to the patient.
	Do not re-sterilize the Depth Electrodes. Re-sterilization may result in device failure.
	The Depth Electrodes are packaged in a double sterile barrier to support surgical asepsis procedures.

**12 STORAGE CONDITIONS**

	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
--	--

Store the Depth Electrodes in their outer packaging, observing the handling and storage conditions detailed below.

	Fragile, handle with care.
--	----------------------------

	Keep away from sunlight.
	Keep dry
	Temperature limit Upper = 30 °C, Lower = 15 °C
	Humidity limitation Upper = 60%, Lower = 20%

**13 DISPOSAL**

	Dispose of the Depth Electrodes following the healthcare facility's biohazardous and pathological waste procedures.
	Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharp's disposal container.

**14 INCIDENT REPORTING**

**EUROPEAN UNION**

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or another person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

**CANADA**

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the following:

- Death of a patient, user, or other person,
- Serious deterioration in the state of health of a patient, user, or other person,
- Potential for death or serious deterioration in the health of a patient, user, or other person if the incident were to recur.

**UNITED STATES**

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening,
- Results in permanent impairment of a body function or permanent damage to a body structure,
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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