DIRECTIONS FOR USE

SPINAL ELECTRODE

STERILE	EO
---------	----

USA
Rx ONLY

CE
2797

EC	REP
----	-----

E.C. REP Ltd
5 Fitzwilliam Square East
Dublin 2, D02 R744, Ireland
Tel: +353 1 2 544 944

TABLE OF CONTENTS

English	Page 4-5
Čeština	Page 6-7
Dansk	Page 8-9
Deutsch	Page 10-11
Español	Page 12-13
Ελληνικά	Page 14-15
Français	Page 16-17
Italiano	Page 18-19
Nederlands	Page 20-21
Turkish	Page 22-23
Slovenčina	Page 24-25
Portuguese	Page 26-27

Please refer to Ad-Tech's Symbols Glossary @ www.adtechmedical.com

DIRECTIONS FOR USE

SPINAL ELECTRODE

Purpose of Electrodes: AD-TECH's Spinal Electrodes are intended for nerve monitoring & recording during spinal surgery. The product should only be used by a physician/surgeon trained in the use of spinal electrodes.

Indications: AD-TECH's Spinal Electrodes are designed for monitoring & recording Motor Evoked Potentials (MEP) of the descending pathway during spinal surgery.

Contraindications: AD-TECH's Spinal Electrodes should not be used on any patient who the physician/surgeon considers at risk for infection or on whom MEP studies cannot safely be performed.

Use: The Spinal Electrode kit is supplied STERILE. Sterility is guaranteed unless package is damaged or seal is broken.

The Spinal Electrode has not been evaluated for safety in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Spinal electrode in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

This product should only be used by a physician/surgeon trained in the use of spinal electrodes. The selection of the style of electrode is made by the physician/surgeon.

The spinal electrode should be handled with care to prevent damage (a direct pull on the electrode may cause a loss of contact recordings).

Placement is selected by physician/surgeon. The spinal electrode is supplied with an introducing stylet or introducing needle depending on the model selected. The spinal electrode should introduce easily and **SHOULD NOT BE FORCED**. After placement, the introducing stylet (if used) can be withdrawn from electrode and discarded. With the stylet removed, the electrode is flexible.

To connect the spinal electrode to an MEP amplifier, use the connection cable integrated with the electrode. Also, consult electrode code chart.

Charge density is the most useful and accurate measure of the intensity of stimulation of tissue by an electrode, because it takes into account not only current strength and pulse duration, but also electrode contact size. Charge density increases as current and/or pulse duration increase and/or contact size decreases. In order to take pulse duration and current strength into consideration for the charge density delivered by the electrode contact supplied, the Stimulation Parameters Table is provided as a guide to their safe use. Contact Ad-Tech if you have any questions.

For Single Use Only. Do Not Re-Sterilize or Reuse.

Not Intended for Implantation (21 CFR 860.3(d): > 30 days). For Surgical Use Only. Do not use if packaging is damaged.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

CAUTION: Reuse of this device is prohibited as it may malfunction and cause contamination and risk to the patient.

CAUTION: Disconnect from monitoring equipment during cardiac defibrillation.

STORAGE: SPINAL electrodes need to be stored and used within typical hospital/office room ambient temperature and humidity conditions.

STIMULATION PARAMETERS

Contact Length: 1.32 mm
Contact Diameter: 1.1 mm

Charge Density ($\mu\text{C}/\text{cm}^2$)	Current (mA)									
	0.2	0.4	0.6	0.8	1	1.5	2	3	4	5
Pulse Duration (ms)	0.1	0.44	0.88	1.31	1.75	2.19	3.29	4.38	8.77	10.96
	0.2	0.88	1.75	2.63	3.51	4.38	6.57	8.77	17.53	21.91
	0.3	1.31	2.63	3.94	5.26	6.57	9.86	13.15	19.72	26.30
	0.4	1.75	3.51	5.26	7.01	8.77	13.15	17.53	26.30	
	0.5	2.19	4.38	6.57	8.77	10.96	16.43	21.91		
	0.75	3.29	6.57	9.86	13.15	16.43	24.65			
	1	4.38	8.77	13.15	17.53	21.91				
	1.25	5.48	10.96	16.43	21.91	27.39				
	1.5	6.57	13.15	19.72	26.30					
	1.75	7.67	15.34	23.01						
	2	8.77	17.53	26.30						
	2.5	10.96	21.91							
	3	13.15	26.30							

$$Charge\ Density\ \frac{\mu C}{cm^2} = \frac{Stimulation\ Intensity\ (mA) * Pulse\ Duration\ (ms)}{Electrode\ Surface\ Area\ (cm^2)}$$

The electrode surface area is 0.046 cm². The maximum charge density limit for safe tissue stimulation is 30 $\mu\text{C}/\text{cm}^2$.

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400 WEST OAKVIEW PARKWAY, OAK CREEK, WISCONSIN, 53154 U.S.A
Toll Free (U.S.A): 800.776.1555