**INSTRUCTIONS FOR USE**

**FORAMEN OVALE  
DEPTH ELECTRODES**

Table of Contents **en**

**1 SYMBOL GLOSSARY** .....2

**2 SCOPE** .....3

**3 IMPORTANT NOTES** .....3

**4 DESCRIPTION** .....3

**5 WARNINGS** .....3

**6 PRECAUTIONS** .....4

**7 DEVICE APPLICATION** .....4

7.1 INTENDED PURPOSE/USE .....4

7.2 INDICATIONS FOR USE (EU) .....4

7.3 CONTRAINDICATIONS .....4

7.4 PATIENT TARGET GROUP .....4

7.5 INTENDED USER .....4

7.6 ENVIRONMENT AND FACILITIES .....4

7.7 EXPECTED CLINICAL BENEFIT .....4

7.8 UNDESIRABLE SIDE EFFECTS .....5

7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE .....5

7.10 ELECTRONIC INSTRUCTIONS FOR USE .....5

**8 PERFORMANCE CHARACTERISTICS** .....5

8.1 DEVICE ATTRIBUTES .....5

8.2 MONITORING CHARACTERISTICS .....5

8.3 STIMULATION CHARACTERISTICS .....5

**9 COMPATIBLE DEVICES/ACCESSORIES** .....6

9.1 INTRODUCER NEEDLE .....6

9.2 CONNECTION SYSTEMS .....6

**10 INSTRUCTIONS FOR USE** .....6

10.1 ELECTRODE PLACEMENT .....7

10.2 CONNECTING THE ELECTRODES .....7

10.3 TROUBLE SHOOTING .....7

10.4 DEVICE REMOVAL .....7

**11 MRI SAFETY INFORMATION** .....7

**12 PACKAGING, STERILIZATION, AND REUSE** .....8

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### 2 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Foramen Ovale Depth Electrodes (FO)

### 3 IMPORTANT NOTES

	<p>Before using these products, carefully review these Instructions for Use.</p> <p>Keep it readily available for reference to the Stimulation Parameters.</p>
--	--

	<p>Carefully read all the symbol-marked cautions and warnings.</p> <p>Incorrect use of the products may lead to severe patient harm.</p>
--	--

### 4 DESCRIPTION



The Depth Electrodes are slender, cylindrical tubes inserted into the brain.

Depth Electrodes have conductive contacts at the distal end to monitor the brain's electrical signals or deliver electrical stimulation.

The contacts connect to internal wires that run throughout the device's body to the tail, where they attach to electrical contacts that allow connection to compatible cable systems.

When connected to compatible equipment, the electrodes serve as the patient-electrode interface for temporary stimulation and recording electrical signals at the subsurface of the brain.



The Depth Electrodes are available in different configurations, including the number of contacts, the contact length, and the contact spacing as specified by the physician. See Figure 1 and Figure 2 for details.

### 5 WARNINGS



	<p>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes.</p> <p>Incorrect use of this device may result in severe patient harm.</p>
	<p>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</p>
	<p>The Depth Electrode is for surgical use only.</p> <p>Placement must be in a sterile surgical environment.</p> <p>Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</p>
	<p>The Depth Electrodes shall be surgically removed.</p>
	<p>The Depth Electrode is for short-term use (&lt;30 days) only.</p> <p>Failure to observe the use conditions may result in severe patient harm.</p>
	<p>Do not use the Depth Electrode if the use-by date has been exceeded.</p> <p>Using devices with an expired shelf life may result in severe patient harm.</p>

**Instructions for Use Foramen Ovale Depth Electrode**

	<p>Do not use the Depth Electrodes if the packaging is damaged or unintentionally opened before use.</p> <p>Using a device from damaged packing may result in contamination and severe patient harm.</p>
	<p>The Depth Electrodes are for single use only. Do not reuse the device.</p> <p>Reusing the device may result in contamination and severe patient harm.</p>
	<p>Do not re-sterilize the Depth Electrodes.</p> <p>Re-sterilization may result in device failure.</p>
	<p>The Depth Electrode should be introduced easily and NOT FORCED.</p> <p>Forcing the Depth Electrode into position may result in severe patient harm.</p>
	<p>Always use appropriate strain relief to secure the Depth Electrode tail and cables.</p> <p>Failure to use appropriate strain relief may result in dislodgement of the electrode and severe patient harm.</p>
	<p>Failure to comply with the recommended stimulation parameter may result in severe patient harm.</p>
	<p>Do not use RF surgical equipment on patients undergoing EEG monitoring or stimulation with the Depth Electrodes.</p> <p>Using RF surgical equipment during EEG monitoring may result in severe patient harm.</p>
	<p>Do not use RF surgical equipment on patients undergoing EEG monitoring with the Depth Electrodes.</p> <p>Using RF surgical equipment during EEG monitoring may result in severe patient harm.</p>
	<p>The Foramen Ovale Depth Electrodes are MRI unsafe.</p> <p>Using the Foramen Ovale Depth Electrodes in an MRI environment may lead to severe patient harm.</p>

**6 PRECAUTIONS**

	<p>Handle the Depth Electrode with care to prevent damage.</p> <p>Pulling, stretching, or twisting the Depth Electrode may result in signal loss.</p>
	<p>Ensure that the Depth Electrode contacts are kept free from moisture and contamination.</p>
	<p>Observe the recommended storage conditions.</p> <p>Failure to observe storage conditions may result in device failure.</p>
	<p>Upon completion of placement and tunneling, connect the electrode contacts outside the body to the electrode cable.</p> <p>At no time should the electrode contacts outside the body be exposed or unprotected.</p> <p>Connection to an electrode cable protects the electrode contacts outside the body.</p> <p>If the patient becomes mobile, the electrode contacts must remain connected to the electrode cable.</p>
	<p>Observe the recommended connection instructions to reduce Electromagnetic Interference (EMI).</p>

	<p>Failure to observe the recommended connection instructions may result in unusable EEG signals.</p>
	<p>Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.</p>
	<p>Dispose of the Depth Electrodes following the healthcare facility's biohazardous and pathological waste procedures.</p>
	<p>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</p>
	<p>Federal Law (USA) restricts this device to sale by or on the order of a physician.</p>

**7 DEVICE APPLICATION**

**7.1 INTENDED PURPOSE/USE**

The Depth Electrodes are for temporary (<30 days) use with recording, monitoring, and stimulation equipment for the recording, monitoring, and stimulation of electrical signals at the subsurface level of the brain. The recording of electrical activity supports the definition of the location of epileptogenic foci and brain mapping.

**7.2 INDICATIONS FOR USE (EU)**

The Depth Electrodes are indicated for patients undergoing diagnostic electrical brain mapping as part of neurological investigations and treatment planning.

**7.3 CONTRAINDICATIONS**

	<p>The Depth Electrodes are not for use on any patient whom the physician/surgeon considers at risk of infection.</p>
	<p>The Depth Electrodes are not intended for continuous stimulation. Stimulation should only be applied to support brain mapping.</p>

**7.4 PATIENT TARGET GROUP**

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

**7.5 INTENDED USER**

The Depth Electrodes are to be used by physicians/surgeons experienced in Depth Electrodes and the associated Placement Accessories.

**7.6 ENVIRONMENT AND FACILITIES**

Depth Electrodes are for surgical use only.

Placement and must be in a sterile surgical environment.

The Depth Electrodes shall be surgically removed.

Foramen Ovale Depth Electrodes are placed at the target location using high-resolution Fluoroscopy.

**7.7 EXPECTED CLINICAL BENEFIT**

Brain mapping allows the physician to gather detailed information about the electrical characteristics of the brain to identify areas related to

**Instructions for Use Foramen Ovale Depth Electrode**

physiological symptoms enabling them to develop appropriate patient-specific treatment plans.

**7.8 UNDESIRABLE SIDE EFFECTS**

Use of the Depth Electrodes carries the following clinical risks:

- Abscess
- Bleeding
- CSF leakage/drainage/pseudomeningocele
- Death
- Epidural hemorrhage
- Intracerebral hemorrhage
- Medical complications including deep vein thrombosis, pulmonary embolisms, allergic reactions, and psychiatric changes
- Meningitis
- Permanent neurologic deficit
- Subarachnoid hemorrhage
- Subdural hemorrhage
- Superficial wound infection
- Transient neurologic deficit
- Worsening seizures

For further details, refer to the Summary of Safety and Clinical Performance (SSCP).

**7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**

The Summary of Safety and Clinical Performance (SSCP) Unique ID: **SSCP-00002**

The SSCP is available from:

- EUDAMED under SRN: US-MF-000004795 and the Unique ID
- AD-TECH® website: [www.adtechmedical.com/SSCP](http://www.adtechmedical.com/SSCP)
- Email request to [sales@adtechmedical.com](mailto:sales@adtechmedical.com), quoting the SSCP Unique ID

**7.10 ELECTRONIC INSTRUCTIONS FOR USE**

An electronic copy of this instructions for use is available from:

- Ad-Tech's website [www.adtechmedical.com](http://www.adtechmedical.com)
- Email request to [sales@adtechmedical.com](mailto:sales@adtechmedical.com), quoting the IFU document number AT10035-1-B

**8 PERFORMANCE CHARACTERISTICS**

**8.1 DEVICE ATTRIBUTES**

**Foramen Ovale Depth Electrode**

Diameter	1.02 mm
Contact Material	Platinum
No. of Contacts	1 - 6
Contact Length Standard (S)	2.41 mm
Contact Diameter Standard (S)	1.10 mm
Recording Area	12 mm – 52 mm
Introducer Needle	17 Gauge (1.016 mm)

**8.2 MONITORING CHARACTERISTICS**

The Depth Electrodes connect to EEG monitoring equipment capable of processing signals greater than 32 Hz with a minimum input voltage of 50 µV.

**8.3 STIMULATION CHARACTERISTICS**

	<p><i>Failure to comply with the recommended stimulation parameter may result in severe patient harm.</i></p>
---	---

The AD-TECH® Depth Electrodes connect to EEG stimulation equipment with a constant-current generator of less than or equal to 34 mA and a constant-voltage stimulator of less than or equal to 10 V.

Charge density is the most valuable and accurate measure of the intensity of tissue stimulation by an electrode, taking into account current strength, pulse duration, and electrode contact size. AD-TECH® has calculated the charge density using the formula:

$$\begin{array}{c} \text{Charge} \\ \text{Density} \end{array} = \frac{\begin{array}{c} \mu\text{C} \\ \text{Stimulation Intensity (mA)} * \text{Pulse} \\ \text{Duration (mS)} \end{array}}{\text{cm}^2} = \frac{\text{Electrode Surface Area (cm}^2\text{)}}{\text{cm}^2}$$

Refer to the Stimulation Parameters tables below for the safe tissue stimulation parameters.

**Instructions for Use Foramen Ovale Depth Electrode**

Foramen Ovale Depth Electrode, S Type Contact														
Contact Length 2.41 mm, Contact Diameter 1.10 mm														
The electrode surface area is 0.081 cm <sup>2</sup> . The maximum charge density limit for safe tissue stimulation is 30 µC/cm <sup>2</sup> .														
Charge Density (µC/cm <sup>2</sup> )	Current (mA)													
	0.2	0.4	0.6	0.8	1	1.5	2	3	4	5	7.5	10	15	
Pulse Duration (ms)	0.1	0.25	0.49	0.74	0.99	1.23	1.85	2.47	3.70	4.94	6.17	9.25	12.34	18.51
	0.2	0.49	0.99	1.48	1.97	2.47	3.70	4.94	7.40	9.87	12.34	18.51	24.68	
	0.3	0.74	1.48	2.22	2.96	3.70	5.55	7.40	11.10	14.81	18.51	27.76		
	0.4	0.99	1.97	2.96	3.95	4.94	7.40	9.87	14.81	19.74	24.68			
	0.5	1.23	2.47	3.70	4.94	6.17	9.25	12.34	18.51	24.68				
	0.75	1.85	3.70	5.55	7.40	9.25	13.88	18.51	27.76					
	1	2.47	4.94	7.40	9.87	12.34	18.51	24.68						
	1.25	3.08	6.17	9.25	12.34	15.42	23.13							
	1.5	3.70	7.40	11.10	14.81	18.51	27.76							
	1.75	4.32	8.64	12.95	17.27	21.59								
	2	4.94	9.87	14.81	19.74	24.68								
	2.5	6.17	12.34	18.51	24.68									
	3	7.40	14.81	22.21	29.61									

**9 COMPATIBLE DEVICES/ACCESSORIES**

**9.1 INTRODUCER NEEDLE**

AD-TECH® provides each Foramen Ovale Depth Electrode with a 17 Gauge Introducer Needle.

**9.2 CONNECTION SYSTEMS**

AD-TECH® provides cables in several standard configurations to suit the user's needs.

Contact AD-TECH® sales support to specify the compatible cable configuration fully.

The following connection cable models are available:

Catalog Number	Description	Compatible Electrodes	Cable Configuration (y)
FO-LDC-yDIN	TECH-ATTACH Connection System Cable	Foramen Ovale Depth Electrode	3 - 6

The following connection blocks are available for use with the TECH-ATTACH Connection System Cable:

Catalog Number	Description	Compatible Cables	Number of Contacts (yy)
2DC-FO-yyX	FO style TECH-ATTACH Disposable Connector Blocks	FO style TECH-ATTACH Connection System Cable, FO-LDC-4DIN, FO-LDC-6DIN (5mm), FO-LDC-6DIN (10mm)	04

**10 INSTRUCTIONS FOR USE**

	<i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes. Incorrect use of this device may result in severe patient harm.</i>
	<i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i>
	<i>The Depth Electrode is for surgical use only. Placement and extraction must be in a sterile surgical environment. Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i>
	<i>The Depth Electrodes shall be surgically removed.</i>
	<i>The Depth Electrode is for short-term use (&lt;30 days) only. Failure to observe the use conditions may result in severe patient harm.</i>
	<i>The Depth Electrode should be introduced easily and NOT FORCED. Forcing the Depth Electrode into position may result in severe patient harm.</i>
	<i>Always use appropriate strain relief to secure the Depth Electrode tail and cables. Failure to use appropriate strain relief may result in dislodgement of the electrode and severe patient harm.</i>
	<i>Disconnect all Depth Electrodes from all monitoring equipment during cardiac defibrillation. Failure to disconnect the Depth Electrodes during cardiac defibrillation may result in severe patient harm.</i>

## Instructions for Use Foramen Ovale Depth Electrode

	<p>Do not use RF surgical equipment on patients undergoing EEG monitoring or stimulation with the Depth Electrodes.</p> <p>Using RF surgical equipment during EEG monitoring may result in severe patient harm.</p>
	<p>Handle the Depth Electrode with care to prevent damage.</p> <p>Pulling, stretching, or twisting the Depth Electrode may result in signal loss.</p>
	<p>Ensure that the Depth Electrode contacts are kept free from moisture and contamination.</p>
	<p>Upon completion of placement and tunneling, connect the electrode contacts outside the body to the electrode cable.</p> <p>At no time should the electrode contacts outside the body be exposed or unprotected.</p> <p>Connection to an electrode cable protects the electrode contacts outside the body.</p> <p>If the patient becomes mobile, the electrode contacts must remain connected to the electrode cable.</p>
	<p>Observe the recommended connection instructions to reduce Electromagnetic Interference (EMI).</p> <p>Failure to observe the recommended connection instructions may result in unusable EEG signals.</p>
	<p>Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.</p>

### 10.1 ELECTRODE PLACEMENT

Foramen Ovale Depth Electrodes are placed bilaterally through the Foramen Ovale into the Ambient Cistern. In this position, the electrodes come to be located intracranially but extracerebrally.

The steps for placing the Depth Electrode are as follows:

**Step 1:** Mark the entry point on the patient's cheek as per pre-operative planning.

**Step 2:** Prepare the patient's cheek with an iodine solution starting at the incision site and circling outwards. Drape the area surrounding the incision site.

**Step 3:** Puncture the skin with the introducer needle, put a finger behind the last molar, and use it to guide the needle tip into position.

**Step 4:** Advance the needle towards the region of the Foramen Ovale under fluoroscopic guidance and use the lateral views provided to determine the position of the needle tip.

**Step 5:** When the needle passes the Foramen Ovale, remove the stylet and replace it with the Foramen Ovale Depth Electrode.

**Step 6:** The electrode should pass through the Foramen Ovale close to the mandibular nerve entering the cranial cavity. Finally, position the electrode in the Ambient Cistern to allow recording from the Mesial temporal lobe area.

**Step 7:** Assess the correct placement of the Foramen Ovale Depth Electrode, ensuring there is no penetration into the foramina of the skull base, such as the inferior orbital fissure and the jugular foramen.

**Step 8:** Remove the introducer needle and secure the electrode's tail to the patient's cheek using surgical tape and a dressing.

### 10.2 CONNECTING THE ELECTRODES

Once the Depth Electrode is secured, immediately connect the exposed contacts to the compatible AD-TECH® Connection System.

To connect to EEG cable, see TECH-ATTACH® Connection instructions for use. Also, consult the electrode code chart provided with the electrodes.

When connecting the electrodes to the EEG monitoring/stimulation equipment, consider the following best practices to reduce the likelihood of Electromagnetic Interference (EMI) and unwanted EEG signal interference/artifacts:

- Ensure that the third-party EEG equipment meets international safety standards, including those for electrical safety and Electromagnetic Compatibility (EMC).
- Connect the EEG equipment following the original equipment manufacturer's instructions for use, including observing any electrical grounding requirements.
- Ensure that the EEG amplifier provides EMI filtering and has a Common Mode Rejection Ratio (CMRR) of at least 80dB at 50/60Hz.
- Reduce any sources of potential EMI by switching off any unnecessary equipment, disconnecting it from the mains network, and removing mobile phones or other telecommunication equipment from the vicinity.
- Utilize a patient ground reference connection.
- Ensure the Depth Electrode external (proximal) contacts are kept clean and dry to avoid high impedance connections or cross-talk.
- Ensure the Connector Blocks are closed tightly around the Depth Electrode proximal contacts to avoid high impedance or intermittent connections.
- Bundle/braid all recording electrode cables together to help reduce 50Hz/60Hz noise generated by the mains network.
- Avoid running the connection cables near or over potential sources of EMI, such as other electrical equipment or connecting cables. If it is necessary to cross cables, ensure they cross perpendicular to each other to limit cross-coupling of EMI.
- Monitor the EEG signals for signs of noise and other signal artifacts that are not being generated by the patient's brain activity and take appropriate action to address the root causes.
- Refer to Section 11, MRI Safety Information, for the safe use of MRI conditional devices in the MRI environment.

Ensure the EEG monitor displays the expected output signals on all electrode contacts to confirm Depth Electrode operation.

During the monitoring period, ensure that the Depth Electrodes are not subject to excessive tensile force, particularly during a patient's seizure.

Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.

### 10.3 TROUBLE SHOOTING

In the event of signal loss, perform the following troubleshooting steps:

1. Confirm the correct connection of the electrode tails with the cable connection system.
2. Verify proper DIN end connection to the third-party EEG system.
3. Confirm the Depth Electrode tail is not under tension and that there are no bends smaller than 1.5 mm in diameter.
4. Confirm there is a proper patient ground reference.
5. For signal artifact issues, consult the EEG equipment manufacturer's troubleshooting guide.

### 10.4 DEVICE REMOVAL

Remove the devices within the specified lifetime of <30 days.

Remove the Depth Electrodes surgically within a sterile surgical environment.

## 11 MRI SAFETY INFORMATION

	<p><i>The Foramen Ovale Depth Electrodes are MRI unsafe.</i></p> <p><i>Using the Foramen Ovale Depth Electrodes in an MRI environment may lead to severe patient harm.</i></p>
---	--

AD-TECH® has not verified the safety and performance of the Foramen Ovale Depth Electrodes in the MRI environment.

AD-TECH® has classified Foramen Ovale Depth Electrodes as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices And Other Items For Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

## Instructions for Use Foramen Ovale Depth Electrode

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

AD-TECH® has labeled the Foramen Ovale Depth Electrodes as MRI Unsafe using the symbol:

### 12 PACKAGING, STERILIZATION, AND REUSE

Observe the following packaging, sterilization, and reuse instructions:

	AD-TECH® sterilizes the Depth Electrodes using ethylene oxide.
	Do not use the Depth Electrode if the use-by date has been exceeded. Using devices with an expired shelf life may result in severe patient harm.
	Do not use the Depth Electrodes if the packaging is damaged or unintentionally opened before use. Using a device from damaged packing may result in contamination and severe patient harm.
	The Depth Electrode is for single use only. Do not reuse the device. Reusing the device may result in contamination and severe harm to the patient.
	Do not re-sterilize the Depth Electrodes. Re-sterilization may result in device failure.
	The Depth Electrodes are packaged in a double sterile barrier to support surgical asepsis procedures.

### 13 STORAGE CONDITIONS

	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
--	--

Store the Depth Electrodes in their outer packaging, observing the handling and storage conditions detailed below.

	Fragile, handle with care.
	Keep away from sunlight.
	Keep dry.
	Temperature limit: Upper = 30°C, Lower = 15°C
	Humidity limitation: Upper = 60%, Lower = 20%

### 14 DISPOSAL

	Dispose of the Depth Electrodes following the healthcare facility's biohazardous and pathological waste procedures.
	Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.

### 15 INCIDENT REPORTING

#### 15.1 EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or another person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

#### 15.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the following:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in the health of a patient, user or other person if the incident were to recur

#### 15.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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