INSTRUCTIONS FOR USE

**BEHNKE FRIED AND WIRE BUNDLE
MACRO-MICRO
DEPTH ELECTRODES**

Table of Contents **en**

1 SYMBOL GLOSSARY..... 2

2 SCOPE..... 3

3 IMPORTANT NOTES 3

4 DESCRIPTION..... 3

5 WARNINGS 4

6 PRECAUTIONS 5

7 DEVICE APPLICATION 5

7.1 INTENDED PURPOSE/USE 5

7.2 INDICATIONS FOR USE (EU) 5

7.3 CONTRAINDICATIONS 5

7.4 PATIENT TARGET GROUP 5

7.5 INTENDED USER..... 5

7.6 ENVIRONMENT AND FACILITIES 5

7.7 EXPECTED CLINICAL BENEFIT 5

7.8 UNDESIRABLE SIDE EFFECTS 5

7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE 5

7.10 ELECTRONIC INSTRUCTIONS FOR USE 6

8 PERFORMANCE CHARACTERISTICS 6

8.1 DEVICE ATTRIBUTES 6

8.2 MONITORING CHARACTERISTICS 6

8.3 STIMULATION CHARACTERISTICS 6

9 COMPATIBLE DEVICES/ACCESSORIES 7

9.1 STAY FLANGE 7

9.2 ANCHOR BOLTS..... 7

9.3 SLOTTED CANNULAS..... 7

9.4 OBTURATORS 7

9.5 DISPOSABLE DRILL KITS..... 7

9.6 RULER 8

9.7 CONNECTION SYSTEMS 8

10 INSTRUCTIONS FOR USE 8

10.1 ELECTRODE PLACEMENT 9

10.2 CONNECTING THE ELECTRODES 9

10.3 TROUBLESHOOTING 10

10.4 DEVICE REMOVAL 10

11 MRI SAFETY INFORMATION..... 10

12 PACKAGING, STERILIZATION, AND REUSE 10

13 STORAGE CONDITIONS..... 10

14 DISPOSAL 10

15 INCIDENT REPORTING..... 11


15.1 EUROPEAN UNION 11

15.2 CANADA 11

15.3 UNITED STATES..... 11

16 ADMINISTRATIVE..... 63

1 SYMBOL GLOSSARY

Symbol	Description
	Manufacturer/Hersteller/Fabricante/Fabricant/ Fabbricante/Fabricante/Üretici
	Authorized representative in the European Community/European Union/ Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/ Representante autorizado en la Comunidad Europea/Unión Europea/Mandataire dans la Communauté européenne/Union européenne/ Mandatario nella Comunità Europea/Unione Europea/Mandatário na Comunidade Europeia/União Europeia/Avrupa Topluluğu/ Avrupa Birliği yetkili temsilcisi
	Date of manufacture/Herstellungsdatum/ Fecha de fabricación/Date de fabrication/ Data di fabbricazione/Data de fabrico/ Üretim tarihi
	Use-by date/Haltbarkeitsdatum/Fecha de caducidad/Date limite d'utilisation/Data di scadenza/Data de validade/Son kullanma tarihi
	Batch code/Chargenbezeichnung/Código de lote/Code de lot/Codice lotto/Código de lote/Parti kodu
	Catalog number/Bestellnummer/Número de catálogo/Numéro de référence/Numero di catalogo/Número de catálogo/Katalog numarası
	Sterilized using ethylene oxide/Sterilisiert mit Ethylenoxid/Esterilizado con óxido de etileno/Stérilisé à l'oxyde d'éthylène/ Sterilizzato con ossido di etilene/Esterilizado por óxido de etileno/Etilen oksit kullanılarak sterilize edilmiştir
	Do not use if the package is damaged and consult instructions for use/Nicht verwenden, wenn Verpackung beschädigt ist und Gebrauchsanweisung beachten/No utilizar si el paquete está dañado y consultar las instrucciones de uso/Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation/Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso/Não utilizar se existirem danos na embalagem e consultar as instruções de utilização/Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
	Do not re-sterilize/Nicht reesterilisieren/ No reesterilizar/Ne pas restériliser/ Non risterilizzare/Não reesterilizar/ Tekrar sterilize etmeyin
	Do not reuse/Nicht wiederverwenden/ No reutilizar/Ne pas réutiliser/Non riutilizzare/ Não reutilizar/Tekrar kullanmayın
	Double sterile barrier system/Doppel-Sterilbarrieresystem/Sistema de barrera estéril doble/Système de double barrière stérile/Sistema a doppia barriera sterile/ Sistema de barreira estéril dupla/Çift steril bariyer sistemi

Instructions for Use Behnke Fried and Wire Bundle Macro-Micro Depth Electrodes

Symbol	Description
	Fragile, handle with care/Zerbrechlich, sorgfältig handhaben/Frágil, manipular con cuidado/Fragile, manipuler avec soin/Fragile, maneggiare con cura/Frágil, manusear com cuidado/Kırılğan, dikkatli taşıyın
	Keep dry/Trocken aufbewahren/Mantener seco/Tenir au sec/Mantenere asciutto/Manter seco/Kuru tutun
	Temperature limit Upper = 30°C, Lower = 15°C Temperaturgrenze obere = 30 °C, untere = 15 °C Limites de temperatura Superior = 30 °C, Inferior = 15 °C Limites de température Supérieure = 30 °C, inférieure = 15 °C Limite di temperatura Superiore = 30 °C, inferiore = 15 °C Limite de temperatura Máximo = 30 °C, Mínimo = 15 °C Sıcaklık sınırı Üst=30 °C, Alt=15 °C
	Humidity limitation Upper = 60%, Lower = 20% Feuchtigkeitsbeschränkung obere = 60 %, untere = 20 % Limites de humedat Superior = 60 %, Inferior = 20 % Limites d'humidité Supérieure = 60 %, inférieure = 20 % Limitazione di umidità Superiore = 60%, inferiore = 20% Limite de humidade Máximo = 60%, Mínimo = 20% Nem sınırları Üst=%60, Alt=%20
	Consult instructions for use/ Gebrauchsanweisung beachten/Consultar las instrucciones de uso/Consulter la notice d'utilisation/Consultare le istruzioni per l'uso/ Consultar as instruções de utilização/ Kullanma talimatına başvurun
	Caution/Vorsicht/Atención/Mise en garde/ Attenzione/Atenção/Dikkat
	Medical device/Medizinprodukt/Producto sanitario/Dispositif médical/Dispositivo medico/Dispositivo médico/Tıbbi cihaz
	Unique device identifier/Einmalige Produktkennung (UDI)/Identificador único de dispositivo/Identifiant unique du dispositif/ Identificativo univoco del dispositivo/ Identificação única do dispositivo/Benzersiz cihaz tanımlayıcı
	A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment./Ein Medizinprodukt, das inakzeptable Risiken für den Patienten, medizinisches Personal oder andere Personen in der MR-Umgebung darstellt./Un producto sanitario que plantea riesgos inaceptables para el paciente, el personal médico u otras personas dentro del entorno de RM./Un dispositif médical présentant des risques inacceptables pour le patient, le personnel médical, ou toute autre personne à l'intérieur de l'environnement de RM./Un dispositivo medico che comporta rischi inaccettabili per il paziente, gli operatori o qualsiasi altro individuo in ambiente RM./Um dispositivo médico que represente

Symbol	Description
	riscos inaceitáveis para o doente, equipa médica ou outras pessoas dentro do ambiente de RM./MR ortamında hasta, sağlık çalışanları veya diğer kişiler için kabul edilemeyecek riskler taşıyan bir tıbbi cihaz. Federal Law (USA) restricts this device to sale by or on the order of a physician./In den USA darf dieses Produkt laut Bundesgesetz nur durch einen Arzt oder auf ärztliche Verordnung verkauft werden./La legislación federal (EE. UU.) restringe la venta de este dispositivo a médicos o por prescripción facultativa./La loi fédérale (États-Unis) limite la vente de ce dispositif par un médecin ou sur ordonnance médicale./Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo a un medico o su prescrizione di un medico./A lei federal (EUA) restringe a venda deste dispositivo a médicos ou mediante prescrição médica./Federal Yasalar (ABD), bu cihazın satışını hekim tarafından veya talimatıyla yapılacak şekilde sınırlamaktadır.

2 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Behnke Fried Depth Electrodes (BF)
- AD-TECH® Wire Bundle Depth Electrodes (WB)

3 IMPORTANT NOTES

	<i>Before using these products, carefully review these Instructions for Use. Keep it readily available for reference to the Stimulation Parameters.</i>
	<i>Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.</i>

4 DESCRIPTION

The Behnke Fried and Wire Bundle Depth Electrodes provide macro and micro monitoring and stimulation.




The Behnke Fried Macro Depth Electrode is a cylindrical tube with conductive contacts at the distal end to monitor the brain's electrical signals or deliver electrical stimulation at the macro level. The contacts connect to internal wires that run throughout the device's body to the tail, where they attach to electrical contacts that allow connection to compatible cable systems. An inner lumen runs the electrode length to accommodate the Wire Bundle Depth Electrode. The Behnke Fried Macro Depth Electrode incorporates a stylet to keep it ridged during placement.

The Wire Bundle Depth Electrode consists of an outer tube with a wire bundle that passes through the inner lumen and protrudes from the distal end. The surgeon trims the wires to the required length to monitor the brain's electrical signals at a micro level. The wire bundle terminates to proximal contacts allowing connection to a compatible cable system.

The Wire Bundle Depth Electrode fits inside the Behnke Fried Depth Electrode, creating a seal to prevent CSF leakage. The surgeon places the complete assembly in the desired location to provide macro monitoring and stimulation.

Instructions for Use Behnke Fried and Wire Bundle Macro-Micro Depth Electrodes

The Behnke Fried Macro Depth Electrodes provides seven to nine macro contacts. The Wire Bundle Depth Electrode provides eight or nine micro-contacts. Figures 1, 2, and 3 for details.

Note #3: Micro wires extend 20 mm past the tip; trim to the desired length before use.

5 WARNINGS

	<p>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes.</p> <p>Incorrect use of this device may result in severe patient harm.</p>
	<p>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</p>
	<p>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</p>
	<p>The Depth Electrode is for surgical use only.</p> <p>Placement must be in a sterile surgical environment.</p> <p>Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</p>
	<p>The Depth Electrodes shall be surgically removed.</p>
	<p>Do not use the Behnke Fried or Wire Bundle Depth Electrodes independently.</p> <p>Failure to use the Behnke Fried or Wire Bundle Depth Electrodes as directed may result in severe patient harm.</p>

	<p>The Depth Electrode is for short-term use (<30 days) only.</p> <p>Failure to observe the use conditions may result in severe patient harm.</p>
	<p>Do not use the Depth Electrode if the use-by date has been exceeded.</p> <p>Using devices with an expired shelf life may result in severe patient harm.</p>
	<p>Do not use the Depth Electrodes if the packaging is damaged or unintentionally opened before use.</p> <p>Using a device from damaged packing may result in contamination and severe patient harm.</p>
	<p>The Depth Electrodes are for single use only. Do not reuse the device.</p> <p>Reusing the device may result in contamination and severe patient harm.</p>
	<p>Do not re-sterilize the Depth Electrodes.</p> <p>Re-sterilization may result in device failure.</p>
	<p>The Depth Electrode should be introduced easily and NOT FORCED.</p> <p>Forcing the Depth Electrode into position may result in severe patient harm.</p>
	<p>Always use appropriate strain relief to secure the Depth Electrode tail and cables.</p> <p>Failure to use appropriate strain relief may result in dislodgement of the electrode and severe patient harm.</p>
	<p>Do not withdraw the wire carrier (Depth Electrode tail) into the incision.</p> <p>Drawing the carrier tail into the incision may contaminate the wound and result in severe patient harm.</p>
	<p>Failure to comply with the recommended stimulation parameter may result in severe patient harm.</p>
	<p>Disconnect all Depth Electrodes from all monitoring equipment during cardiac defibrillation.</p> <p>Failure to disconnect the Depth Electrodes during cardiac defibrillation may result in severe patient harm.</p>
	<p>Do not use RF surgical equipment on patients undergoing EEG monitoring or stimulation with the Depth Electrodes.</p> <p>Using RF surgical equipment during EEG monitoring may result in severe patient harm.</p>
	<p>The Behnke Fried and Wire Bundle Depth Electrodes MRI unsafe.</p> <p>Using the Behnke Fried and Wire Bundle Depth Electrodes in an MRI environment may severely harm the user and patient.</p>

6 PRECAUTIONS

	<p>Handle the Depth Electrode with care to prevent damage.</p> <p>Pulling, stretching, or twisting the Depth Electrode may result in signal loss.</p>
	<p>Ensure that the Depth Electrode contacts are kept free from moisture and contamination.</p>
	<p>Do not remove the stylet before placement of the Depth Electrode. Removing and replacing the stylet before placement may damage the device. The stylet ensures that the Depth Electrode remains rigid during the delivery procedure.</p>
	<p>Observe the recommended storage conditions.</p> <p>Failure to observe storage conditions may result in device failure.</p>
	<p>Upon completion of placement and tunneling, connect the electrode contacts outside the body to the electrode cable.</p> <p>At no time should the electrode contacts outside the body be exposed or unprotected.</p> <p>Connection to an electrode cable protects the electrode contacts outside the body.</p> <p>If the patient becomes mobile, the electrode contacts must remain connected to the electrode cable.</p>
	<p>Observe the recommended connection instructions to reduce Electromagnetic Interference (EMI).</p> <p>Failure to observe the recommended connection instructions may result in unusable EEG signals.</p>
	<p>It is essential to place and secure a strain relief using the EEG cables so there will not be a direct pull on the electrodes.</p> <p>A direct pull on the electrode may cause a loss of contact recordings or other damage.</p>
	<p>Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.</p>
	<p>Dispose of the Depth Electrodes following the healthcare facility's biohazardous and pathological waste procedures.</p>
	<p>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</p>
	<p>Federal Law (USA) restricts this device to sale by or on the order of a physician.</p>

7 DEVICE APPLICATION

7.1 INTENDED PURPOSE/USE

The Depth Electrodes are for temporary (<30 days) use with recording, monitoring, and stimulation equipment for the recording, monitoring, and stimulation of electrical signals at the subsurface level of the brain. The recording of electrical activity supports the definition of the location of epileptogenic foci and brain mapping.

7.2 INDICATIONS FOR USE (EU)

The Depth Electrodes are indicated for patients undergoing diagnostic electrical brain mapping as part of neurological investigations and treatment planning.

7.3 CONTRAINDICATIONS

	<p>The Depth Electrodes are not for use on any patient whom the physician/surgeon considers at risk of infection.</p>
	<p>The Depth Electrodes are not intended for continuous stimulation. Stimulation should only be applied to support brain mapping.</p>

7.4 PATIENT TARGET GROUP

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

7.5 INTENDED USER

The Depth Electrodes are to be used by physicians/surgeons experienced in Depth Electrodes and the associated Placement Accessories.

7.6 ENVIRONMENT AND FACILITIES

Depth Electrodes are for surgical use only. Placement must be in a sterile surgical environment. The Depth Electrodes shall be surgically removed. Depth Electrodes are placed at the target location using an appropriate stereotactic guidance system and supporting equipment.

7.7 EXPECTED CLINICAL BENEFIT

Brain mapping allows the physician to gather detailed information about the electrical characteristics of the brain to identify areas related to physiological symptoms enabling them to develop appropriate patient-specific treatment plans.

7.8 UNDESIRABLE SIDE EFFECTS

Use of the Depth Electrodes carries the following clinical risks:

- Abscess
- Bleeding
- CSF leakage/drainage/pseudomeningocele
- Death
- Epidural hemorrhage
- Intracerebral hemorrhage
- Medical complications including deep vein thrombosis, pulmonary embolisms, allergic reactions, and psychiatric changes
- Meningitis
- Permanent neurologic deficit
- Subarachnoid hemorrhage
- Subdural hemorrhage
- Superficial wound infection
- Transient neurologic deficit
- Worsening seizures

For further details, refer to the Summary of Safety and Clinical Performance (SSCP).

7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) Unique ID: **SSCP-00002**

The SSCP is available from:

Instructions for Use Behnke Fried and Wire Bundle Macro-Micro Depth Electrodes

- EUDAMED under SRN: US-MF-000004795 and the Unique ID
- Email request to sales@adtechmedical.com, quoting the SSCP Unique ID

7.10 ELECTRONIC INSTRUCTIONS FOR USE

An electronic copy of this instructions for use is available from:

- Ad-Tech's website www.adtechmedical.com
- Email request to sales@adtechmedical.com, quoting the IFU document number AT10036-1-B

8 PERFORMANCE CHARACTERISTICS

8.1 DEVICE ATTRIBUTES

Behnke Fried Depth Electrode	
Diameter	1.28 mm
Contact Material	Platinum
No. of Contacts	7-12
Contact Length	1.57 mm
Contact Diameter	1.28 mm
Contact Spacing	3.0 mm - 10.0 mm
Recording Area	31.57 mm – 111.57 mm
Stylet	Stainless Steel

Wire Bundle Depth Electrode	
Diameter	0.038 mm
Contact Material	Platinum
No. of Contacts	Up to 9
Contact Length	20 mm (Maximum)
Contact Diameter	38 Micron / 0.44 mm

8.2 MONITORING CHARACTERISTICS

The Depth Electrodes connect to EEG monitoring equipment capable of processing signals greater than 32 Hz with a minimum input voltage of 50 µV.

8.3 STIMULATION CHARACTERISTICS

	<i>Failure to comply with the recommended stimulation parameter may result in severe patient harm.</i>
---	--

The Depth Electrodes connect to EEG stimulation equipment with a constant-current generator of less than or equal to 34 mA and a constant-voltage stimulator of less than or equal to 10 V.

Charge density is the most valuable and accurate measure of the intensity of tissue stimulation by an electrode, taking into account current strength, pulse duration, and electrode contact size. AD-TECH® has calculated the charge density using the formula:

$$\text{Charge Density} \frac{\mu\text{C}}{\text{cm}^2} = \frac{\text{Stimulation Intensity (mA)} * \text{Pulse Duration (mS)}}{\text{Electrode Surface Area (cm}^2\text{)}}$$

Instructions for Use Behnke Fried and Wire Bundle Macro-Micro Depth Electrodes

Refer to the Stimulation Parameters tables below for the safe tissue stimulation parameters.

Note: The micro contacts provided by the Wire Bundle are not used for stimulation.

Behnke Fried Depth Electrode, S Type Contact														
Contact Length 1.57 mm, Contact Diameter 1.28 mm														
The electrode surface area is 0.061 cm ² . The maximum charge density limit for safe tissue stimulation is 30 µC/cm ² .														
Charge Density (µC/cm ²)	Current (mA)													
	0.2	0.4	0.6	0.8	1	1.5	2	3	4	5	7.5	10	15	
Pulse Duration (ms)	0.1	0.33	0.66	0.99	1.32	1.64	2.47	3.29	4.93	6.58	8.22	12.33	16.45	24.67
	0.2	0.66	1.32	1.97	2.63	3.29	4.93	6.58	9.87	13.16	16.45	24.67		
	0.3	0.99	1.97	2.96	3.95	4.93	7.40	9.87	14.80	19.74	24.67			
	0.4	1.32	2.63	3.95	5.26	6.58	9.87	13.16	19.74	26.31				
	0.5	1.64	3.29	4.93	6.58	8.22	12.33	16.45	24.67					
	0.75	2.47	4.93	7.40	9.87	12.33	18.50	24.67						
	1	3.29	6.58	9.87	13.16	16.45	24.67							
	1.25	4.11	8.22	12.33	16.45	20.56								
	1.5	4.93	9.87	14.80	19.74	24.67								
	1.75	5.76	11.51	17.27	23.02	28.78								
	2	6.58	13.16	19.74	26.31									
	2.5	8.22	16.45	24.67										
	3	9.87	19.74	29.60										

9 COMPATIBLE DEVICES/ACCESSORIES

9.1 STAY FLANGE

The AD-TECH® Stay Flange is an optional accessory for the Depth Electrodes. The Stay Flange provides strain relief to the Depth Electrode tail, anchoring it via sutures to the patient's skull. Stay Flanges are available separately using part number FLANGE-X.

9.2 ANCHOR BOLTS

The AD-TECH® Anchor Bolts are optional accessories for use with the Behnke Fried Depth Electrode. Anchor Bolts are secured in the skull to provide an access point. The Anchor Bolts minimize concerns about potential cerebrospinal fluid (CSF) leakage and infection of the subdural space while stabilizing the electrode. The Behnke Fried Depth Electrodes are compatible with the following Anchor Bolts:

Catalog #	Description	Material
LSBK1-AX-06 (Clear)	Use with BF, LD, or MM Depth Electrodes	Titanium
LSBK1-BX-06 (Clear)	Use with BF, LD, or MM Depth Electrodes	Titanium
LSBK1-CX-06 (Clear)	Use with BF, LD, or MM Depth Electrodes	Titanium

9.3 SLOTTED CANNULAS

The AD-TECH® Slotted Cannulas are optional accessories for use with the Behnke Fried Depth Electrode. Slotted Cannulas create a path within the brain to facilitate the Depth Electrode placement. A Slotted Cannula is applied to place a Depth Electrode through a craniotomy or a burr hole without using an Anchor Bolt. The following slotted cannulas are available:

Catalog #	Description
2SCK1-190X	Two-Piece Slotted Cannula - 190mm (2.11 diameter) 2/Box Sterile, Reusable.
2SCK1-240X	Two-Piece Slotted Cannula - 240mm (2.11 diameter) 2/Box Sterile, Reusable.

9.4 OBTURATORS

The AD-TECH® Obturators are optional accessories for use with the Behnke Fried Depth Electrode. Obturators create a path within the brain to facilitate the placement of Depth Electrodes. The Obturator is applied when placing Depth Electrodes through a burr hole with an applied Anchor Bolt. The following Obturators are available:

Catalog #	Description
OB-20-190X	An obturator is used to make a track inside the brain. It is used with RD, SD, MM, BF/WB, and LD style electrodes; it has a 0.86mm diameter and 190mm length. Sterile, Reusable
OB-17-240X	An obturator is used to make a track inside the brain. It is used with SD, MM, BF/WB, and LD style electrodes; it has a 1.17mm diameter and 240mm length. Sterile, Reusable

9.5 DISPOSABLE DRILL KITS

The AD-TECH® Disposable Drill Kits are an optional accessory with the Behnke Fried Depth Electrode. When used with a compatible Surgical Drill, the Drill Kit makes a burr hole to accept the Anchor Bolt or Depth Electrode. The following Drill Kits are available:

Instructions for Use Behnke Fried and Wire Bundle Macro-Micro Depth Electrodes

Catalog Number	Description	Drill			Drill Stop Quantity	Drill Wrench Quantity	Sterile Yes/No
		Quantity	Length (YY)	Diameter (a.a)			
DDK2-2.4-30N	Disposable Drill Kit – Non-sterile	2	30 cm	2.4 mm	2	2	No
DRL-2.4-22N	Disposable Drill Bit – Non-sterile	1	22 cm	2.4 mm	0	0	No
DDK2-2.4-16X	Disposable Drill Kit – Sterile	2	16 cm	2.4 mm	2	2	Yes
DDK2-2.4-30X	Disposable Drill Kit – Sterile	2	30 cm	2.4 mm	2	2	Yes

9.6 RULER

The AD-TECH® Ruler is an optional accessory to the Behnke Fried and Wire Bundle Depth Electrodes. The Ruler is used to make critical measurements when placing the Depth Electrodes.

Catalog #	Description
RULER-20N	20cm Channeled Ruler for Drill Stop & Depth Electrode measurement. Non-sterile, Reusable.

9.7 CONNECTION SYSTEMS

AD-TECH® provides cables in several standard configurations to suit the user's needs.

Contact AD-TECH® sales support to specify the compatible cable configuration fully.

The following connection cable models are available:

Catalog Number	Description	Compatible Electrodes	Cable Configuration (y)
L-DC-yDIN-aa*	TECH-ATTACH Connection System Cable	Behnke Fried Depth Electrodes Wire Bundle Depth Electrode	7-12 8 or 9
L-DCL-yDIN-aa*	TECH-ATTACH Connection System Cable	Behnke Fried Depth Electrodes Wire Bundle Depth Electrode	7-12 8 or 9
L-SRL-yDIN-aa*	Light Weight CABRIO Cables	Wire Bundle Depth Electrode	8 or 9

*Note: 'aa' denotes the standard configurations available from ADTECH®. Contact an AD-TECH® sales representative for support.

Catalog Number	Description	Compatible Cables	Number of Contacts (y)
DC-yX-aa*	y-contact connection block**	TECH-ATTACH Connection System Cable	7-10
DC-ySC	y-contact slotted connection block	TECH-ATTACH Connection System Cable	7-10

** A special slotted connector block that connects the macro contacts to the TECH-ATTACH cable is included with each BF-series electrode.

10 INSTRUCTIONS FOR USE

	<i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes. Incorrect use of this device may result in severe patient harm.</i>
	<i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i>
	<i>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</i>
	<i>The Depth Electrode is for surgical use only. Placement must be in a sterile surgical environment. Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i>
	<i>The Depth Electrodes shall be surgically removed.</i>
	<i>Do not use the Behnke Fried or Wire Bundle Depth Electrodes independently. Failure to use the Behnke Fried or Wire Bundle Depth Electrodes as directed may result in severe patient harm.</i>
	<i>The Depth Electrode is for short-term use (<30 days) only. Failure to observe the use conditions may result in severe patient harm.</i>
	<i>The Depth Electrode should be introduced easily and NOT FORCED. Forcing the Depth Electrode into position may result in severe patient harm.</i>
	<i>Always use appropriate strain relief to secure the Depth Electrode tail and cables. Failure to use appropriate strain relief may result in dislodgement of the electrode and severe patient harm.</i>

The following connection blocks are available for use with the TECH-ATTACH Connection System Cable:

Instructions for Use Behnke Fried and Wire Bundle Macro-Micro Depth Electrodes

	<p><i>Do not withdraw the wire carrier (Depth Electrode tail) into the incision.</i></p> <p><i>Drawing the carrier tail into the incision may contaminate the wound and result in severe patient harm.</i></p>
	<p><i>Disconnect all Depth Electrodes from all monitoring equipment during cardiac defibrillation.</i></p> <p><i>Failure to disconnect the Depth Electrodes during cardiac defibrillation may result in severe patient harm.</i></p>
	<p><i>Do not use RF surgical equipment on patients undergoing EEG monitoring or stimulation with the Depth Electrodes.</i></p> <p><i>Using RF surgical equipment during EEG monitoring may result in severe patient harm.</i></p>
	<p><i>Handle the Depth Electrode with care to prevent damage.</i></p> <p><i>Pulling, stretching, or twisting the Depth Electrode may result in signal loss.</i></p>
	<p><i>Do not remove the stylet before placement of the Depth Electrode. Removing and replacing the stylet before placement may damage the device. The stylet ensures that the Depth Electrode remains rigid during the delivery procedure.</i></p>
	<p><i>Upon completion of placement and tunneling, connect the electrode contacts outside the body to the electrode cable.</i></p> <p><i>At no time should the electrode contacts outside the body be exposed or unprotected.</i></p> <p><i>Connection to an electrode cable protects the electrode contacts outside the body.</i></p> <p><i>If the patient becomes mobile, the electrode contacts must remain connected to the electrode cable.</i></p>
	<p><i>Ensure that the Depth Electrode contacts are kept free from moisture and contamination.</i></p>
	<p><i>Observe the recommended connection instructions to reduce Electromagnetic Interference (EMI).</i></p> <p><i>Failure to observe the recommended connection instructions may result in unusable EEG signals.</i></p>
	<p><i>Ensure that the Depth Electrode implant sites are kept clean and dry.</i></p> <p><i>Monitor for any signs of infection.</i></p>

10.1 ELECTRODE PLACEMENT

The physician can place the Depth Electrode through a standard twist drill or burr hole, or craniotomy site.

Prepare the path through the brain tissue using the chosen method, Obturator or Slotted Cannula.

Regarding Figure 4, the steps for placing the Depth Electrode are as follows:

Step 1: Cutting the micro-wires to the desired length

- a) Slide the micro-wire Cover Tube (Green) over the micro-wires to cover them.
- b) Insert the covered micro-wires 4 – 5 cm into the lumen of the Behnke Fried Depth Electrode.
- c) Pull the Cover Tube to the proximal end of the Wire Bundle Depth Electrode until secure.

- d) Complete the insertion of the Wire Bundle Depth Electrode into the Behnke Fried Depth Electrode until the taper fittings tighten.
- e) Ensure that the Cover Tube does not enter completely into the inner lumen of the Behnke Fried Depth Electrode.
- f) The micro-wires will now be protruding from the end of the assembled Depth Electrode.
- g) Measure and trim the micro-wires to the desired length.
- h) Safely discard the trimmings to avoid contamination of the patient.
- i) Separate the wire bundle from the Behnke Fried and reinsert the stylet into the Behnke Fried.

Step 2:

- a) Slide the Depth Electrode along the pre-planned/prepared trajectory to the target location. Withdraw the stylet while maintaining the Depth Electrode Position.

When using the Depth Electrode with the Anchor Bolt, follow the placement instructions in the Anchor Bolt instructions for use.

When using the Depth Electrode in combination with the Slotted Cannula, follow the placement instructions in the Slotted Cannula instructions for use.

When using the Depth Electrode with the Obturator, follow the instructions in the Obturator instructions for use.

10.2 CONNECTING THE ELECTRODES

Once the Depth Electrode is secured, immediately connect the exposed contacts to the compatible AD-TECH® Connection System.

To connect the EEG cable, see the CABRIO®/TECH-ATTACH® Connection Systems instructions for use. Also, consult the electrode code chart provided with the electrodes.

When connecting the electrodes to the EEG monitoring/stimulation equipment, consider the following best practices to reduce the likelihood of Electromagnetic Interference (EMI) and unwanted EEG signal interference/artifacts:

- Ensure that the third-party EEG equipment meets international safety standards, including those for electrical safety and Electromagnetic Compatibility (EMC).
- Connect the EEG equipment following the original equipment manufacturer's instructions for use, including observing any electrical grounding requirements.
- Ensure that the EEG amplifier provides EMI filtering and has a Common Mode Rejection Ratio (CMRR) of at least 80dB at 50/60Hz.
- Reduce any sources of potential EMI by switching off any unnecessary equipment, disconnecting it from the mains network, and removing mobile phones or other telecommunication equipment from the vicinity.
- Utilize a patient ground reference connection.
- Ensure the Depth Electrode external (proximal) contacts are kept clean and dry to avoid high impedance connections or cross-talk.
- Ensure the Connector Blocks are closed tightly around the Depth Electrode proximal contacts to avoid high impedance or intermittent connections.
- Bundle/braid all recording electrode cables together to help reduce 50Hz/60Hz noise generated by the mains network.
- Avoid running the connection cables near or over potential sources of EMI, such as other electrical equipment or connecting cables. If it is necessary to cross cables, ensure they cross perpendicular to each other to limit cross-coupling of EMI.
- Monitor the EEG signals for signs of noise and other signal artifacts that are not being generated by the patient's brain activity and take appropriate action to address the root causes.
- Refer to Section 11, MRI Safety Information, for the safe use of MRI conditional devices in the MRI environment.

Ensure the EEG monitor displays the expected output signals on all electrode contacts to confirm Depth Electrode operation.

During the monitoring period, ensure that the Depth Electrodes are not subject to excessive tensile force, particularly during a patient's seizure.

Ensure that the Depth Electrode implant sites are kept clean and dry. Monitor for any signs of infection.

Instructions for Use Behnke Fried and Wire Bundle Macro-Micro Depth Electrodes

10.3 TROUBLESHOOTING

In the event of signal loss, perform the following troubleshooting steps:

1. Confirm the correct connection of the electrode tails with the cable connection system.
2. Verify proper DIN end connection to the third-party EEG system.
3. Confirm the Depth Electrode tail is not under tension and that there are no bends smaller than 1.5 mm in diameter.
4. Confirm there is a proper patient ground reference.
5. For signal artifact issues, consult the EEG equipment manufacturer's troubleshooting guide.

10.4 DEVICE REMOVAL

Remove the devices within the specified lifetime of <30 days.

Remove the Depth Electrodes surgically within a sterile surgical environment.

11 MRI SAFETY INFORMATION

	<p><i>The Behnke Fried and Wire Bundle Depth Electrodes MRI unsafe.</i></p> <p><i>Using the Behnke Fried and Wire Bundle Depth Electrodes in an MRI environment may severely harm the user and patient.</i></p>
---	---

AD-TECH® has not verified the safety and performance of the Behnke Fried and Wire Bundle Depth Electrodes in the MRI environment.

AD-TECH® has classified the Behnke Fried and Wire Bundle Depth Electrodes as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices And Other Items For Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

AD-TECH® has labeled the Behnke Fried and Wire Bundle Depth Electrodes as MRI Unsafe using the symbol:

12 PACKAGING, STERILIZATION, AND REUSE

Observe the following packaging, sterilization, and reuse instructions:

	<p><i>AD-TECH® sterilizes the Depth Electrodes using ethylene oxide.</i></p>
	<p><i>Do not use the Depth Electrode if the use-by date has been exceeded.</i></p> <p><i>Using devices with an expired shelf life may result in severe patient harm.</i></p>
	<p><i>Do not use the Depth Electrodes if the packaging is damaged or unintentionally opened before use.</i></p> <p><i>Using a device from damaged packing may result in contamination and severe patient harm.</i></p>

	<p><i>The Depth Electrode is for single use only. Do not reuse the device.</i></p> <p><i>Reusing the device may result in contamination and severe harm to the patient.</i></p>
	<p><i>Do not re-sterilize the Depth Electrodes.</i></p> <p><i>Re-sterilization may result in device failure.</i></p>
	<p><i>The Depth Electrodes are packaged in a double sterile barrier to support surgical asepsis procedures.</i></p>

13 STORAGE CONDITIONS

	<p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p>
---	---

Store the Depth Electrodes in their outer packaging, observing the handling and storage conditions detailed below.

	<p><i>Fragile, handle with care.</i></p>
	<p><i>Keep away from sunlight.</i></p>
	<p><i>Keep dry.</i></p>
	<p><i>Temperature limit:</i></p> <p><i>Upper = 30°C, Lower = 15°C</i></p>
	<p><i>Humidity limitation:</i></p> <p><i>Upper = 60%, Lower = 20%</i></p>

14 DISPOSAL

	<p><i>Dispose of the Depth Electrodes following the healthcare facility's biohazardous and pathological waste procedures.</i></p>
	<p><i>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</i></p>

15 INCIDENT REPORTING

15.1 EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or another person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

15.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the following:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in the health of a patient, user or other person if the incident were to recur

15.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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