**INSTRUCTIONS FOR USE**  
**ANCHOR BOLT**

**1 SYMBOL GLOSSARY ..... 2**

**2 SCOPE..... 3**

**3 IMPORTANT NOTES ..... 3**

**4 DESCRIPTION ..... 3**

**5 WARNINGS ..... 4**

**6 PRECAUTIONS ..... 4**

**7 DEVICE APPLICATION ..... 4**

7.1 INTENDED USE/PURPOSE ..... 4

7.2 INDICATIONS FOR USE ..... 4

7.3 CONTRAINDICATIONS ..... 4

7.4 PATIENT TARGET GROUP ..... 4

7.5 INTENDED USER ..... 4

7.6 ENVIRONMENT AND FACILITIES ..... 5

7.7 EXPECTED CLINICAL BENEFIT ..... 5

7.8 UNDESIRABLE SIDE EFFECTS ..... 5

7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE ..... 5

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### 2 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® LSB Series Anchor Bolts

### 3 IMPORTANT NOTES

	<p>Before using these products, carefully review these Instructions for Use. Keep it readily available for reference.</p>
	<p>Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.</p>

### 4 DESCRIPTION

The Anchor Bolts are optional accessories to support the placement of the following AD-TECH® Depth Electrodes:

- AD-TECH® Spencer Probe Depth Electrodes, Styles LD, RD, and SD
- AD-TECH® Behnke Fried Depth Electrode (BF)
- AD-TECH® Macro-Micro Depth Electrode (MM)

Anchor Bolts are secured in the skull to provide an access point and minimize concerns about potential cerebrospinal fluid (CSF) leakage and infection of the subdural space while stabilizing the Electrode.

The Anchor Bolt consists of a tube through which the Depth Electrode can pass. The proximal end is threaded, allowing the surgeon to screw it into the patient's skull via an appropriately sized burr hole. The distal end is threaded to accept the Anchor Bolt Cap.

The surgeon passes the Anchor Bolt Cap over the distal end of the Depth Electrode and screws it tight, compressing the internal Gasket against its outer wall. The compression between the Gasket and the Electrode Body is the primary seal to prevent CSF leakage.

A Silicone Cap protects the Anchor Bolt Cap and acts as additional strain relief for the Depth Electrode tail.

The Anchor Bolts are available in Titanium Grade 2 and Anodized Titanium Grade 2.

The color of the Anchor Bolt gasket denotes the compatible Depth Electrodes as follows:

## Instructions for Use Anchor Bolt

- Spencer Probe Depth Electrode, SD Style = Green
- Spencer Probe Depth Electrode, RD Style = Blue
- Spencer Probe Depth Electrode, LD Style = Clear
- Behnke Fried (BF) and Macro-Micro Depth Electrodes (MM) = Clear

### 5 WARNINGS

	<p><i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories.</i></p> <p><i>Incorrect use of this device may result in severe patient harm.</i></p>
	<p><i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i></p>
	<p><i>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</i></p>
	<p><i>The Anchor Bolt is for surgical use only.</i></p> <p><i>Placement must be in a sterile surgical environment.</i></p> <p><i>Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i></p>
	<p><i>The Anchor Bolts shall be surgically removed.</i></p>
	<p><i>The Anchor Bolt is for short-term use (&lt;30 days) only.</i></p> <p><i>Failure to observe the use conditions may result in severe patient harm.</i></p>
	<p><i>Do not use the Anchor Bolt if the packaging is damaged or unintentionally opened before use.</i></p> <p><i>Using a device from damaged packing may result in contamination and severe patient harm.</i></p>
	<p><i>The Anchor Bolt is for single use only. Do not reuse the device.</i></p> <p><i>Reuse of the device may result in contamination and severe patient harm.</i></p>
	<p><i>Do not re-sterilize the Anchor Bolts.</i></p> <p><i>Re-sterilization may result in device failure.</i></p>
	<p><i>The Anchor Bolts are MRI unsafe.</i></p> <p><i>Using the Anchor Bolts in an MRI environment may lead to severe patient harm.</i></p>

### 6 PRECAUTIONS

	<p><i>Handle the Anchor Bolts with care to prevent damage.</i></p> <p><i>Do not apply excessive torque to the Anchor Bolt during placement. When using the recommended burr hole diameter, the Anchor Bolt will screw easily into the skull.</i></p> <p><i>Excess torque may result in a bent or broken anchor bolt.</i></p>
---	--

	<p><i>Do not use the Anchor Bolts if the third-party interface (e.g., stereotactic/head frame) is damaged or flawed.</i></p>
	<p><i>Discontinue Anchor Bolt placement if you experience any instrument (e.g., the Placement/Removal Wrench, Cranial Drill, etc.) bending during the placement procedure.</i></p>
	<p><i>Protect the Anchor Bolt from the blunt impact that may occur if the patient has a seizure. Significant force may fracture the Anchor Bolt.</i></p>
	<p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p>
	<p><i>Dispose of the Anchor Bolts following the healthcare facility's biohazardous and pathological waste procedures.</i></p>
	<p><i>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</i></p>
	<p><i>Federal Law (USA) restricts this device to sale by or on the order of a physician.</i></p>

### 7 DEVICE APPLICATION

#### 7.1 INTENDED USE/PURPOSE

The Anchor Bolts are optional accessories for temporary use (<30 days) with compatible Depth Electrodes. The Anchor Bolts may be applied when desired to minimize concerns about potential cerebrospinal fluid (CSF) leakage and infection of the subdural space while stabilizing the Depth Electrode. The Anchor Bolts are secured in the skull to provide an access point for Depth Electrodes.

#### 7.2 INDICATIONS FOR USE

The Ad-Tech Anchor Bolts are optional accessories for use with Depth Electrodes. The Anchor Bolts may be applied when it is desired to minimize concerns about potential cerebrospinal fluid (CSF) leakage and infection of the subdural space while stabilizing the Electrode. Anchor Bolts are secured in the skull to provide an access point for Depth Electrodes.

#### 7.3 CONTRAINDICATIONS

	<p><i>The Anchor Bolts are not for use on any patient whom the physician/surgeon considers at risk of infection.</i></p>
	<p><i>The Anchor Bolts are not for use with patients with low bone density or softening of the skull.</i></p>

#### 7.4 PATIENT TARGET GROUP

The Anchor Bolts are used in pediatric and adult populations, regardless of gender.

#### 7.5 INTENDED USER

The Anchor Bolts are to be used by physicians/surgeons experienced in Depth Electrodes and the associated Placement Accessories.

## Instructions for Use Anchor Bolt

### 7.6 ENVIRONMENT AND FACILITIES

Anchor Bolts are for surgical use only.

Placement and removal must be in a sterile surgical environment.

Anchor Bolts are placed at the target location using an appropriate stereotactic guidance system and supporting equipment.

The surgeon must use the AD-TECH® LSB Anchor Bolt Wrench to place and remove the Anchor Bolts.

### 7.7 EXPECTED CLINICAL BENEFIT

When used, the Anchor Bolts minimize concerns about potential cerebrospinal fluid (CSF) leakage and infection of the subdural space while stabilizing the Depth Electrode.

### 7.8 UNDESIRABLE SIDE EFFECTS

Use of the Depth Electrodes and the Anchor Bolts carries the following clinical risks:

- Abscess
- Bleeding
- CSF leakage/drainage/pseudomeningocele
- Death
- Epidural hemorrhage
- Intracerebral hemorrhage
- Medical complications including deep vein thrombosis, pulmonary embolisms, allergic reactions, and psychiatric changes
- Meningitis
- Permanent neurologic deficit
- Subarachnoid hemorrhage
- Subdural hemorrhage
- Superficial wound infection
- Transient neurologic deficit
- Worsening seizures

For further details, refer to the Summary of Safety and Clinical Performance (SSCP)

### 7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) Unique ID: **SSCP-00002**

The SSCP is available from:

- EUDAMED under SRN: US-MF-000004795 and the Unique ID
- Email request to [sales@adtechmedical.com](mailto:sales@adtechmedical.com), quoting the SSCP Unique ID

### 7.10 ELECTRONIC INSTRUCTIONS FOR USE

An electronic copy of this instructions for use is available from:


- Ad-Tech's website [www.adtechmedical.com](http://www.adtechmedical.com)
- Email request to [sales@adtechmedical.com](mailto:sales@adtechmedical.com), quoting the IFU document number AT10042-LSB-B

## 8 PERFORMANCE CHARACTERISTICS

### 8.1 DEVICE ATTRIBUTES

AD-TECH® designed the Anchor Bolts to fracture if subjected to excessive side loads such as a blunt impact during a patient's seizure. Therefore, protecting the patient's skull from traumatic injury.

The following AD-TECH® Anchor Bolts are available:

Catalog Number (Gasket Color)	Description	Distal Opening Diameter (D)	Post Length (L)
<b>Titanium</b>			
LSBK1-AX-05 (Green)	Use with SD Depth Electrodes	1.5 mm	21 mm
LSBK1-AX-06 (Clear)	Use with BF, LD, or MM Depth Electrodes	1.5 mm	21 mm
LSBK1-BX-05 (Green)	Use with SD Depth Electrodes	1.5 mm	13 mm
LSBK1-BX-06 (Clear)	Use with BF, LD, or MM Depth Electrodes	1.5 mm	13 mm
LSBK1-CX-05 (Green)	Use with SD Depth Electrodes	1.5 mm	26 mm
LSBK1-CX-06 (Clear)	Use with BF, LD, or MM Depth Electrodes	1.5 mm	26 mm
LSBK1-CX-04 (Blue)	Use with RD Depth Electrodes	1.5 mm	26 mm
<b>Anodized Titanium</b>			
LSBK2-AX-04 (Blue)	Use with RD Depth Electrodes	0.99 mm	21 mm
LSBK2-BX-04 (Blue)	Use with RD Depth Electrodes	0.99 mm	13 mm
			

## 9 COMPATIBLE DEVICES/ACCESSORIES

### 9.1 LSB ANCHOR BOLT WRENCH

Use the AD-TECH® LSB Anchor Bolt Wrench to screw the Anchor Bolt into the patient's skull and remove it at the end of the monitoring period.

Catalog Number	Shaft Diameter (inches)	Length (inches)	Supplied Sterile Yes/No
LSB-PWL-2.4-6X	0.094 (2.4 mm)	6	Yes
LSB-PWL-2.4-10X	0.094 (2.4 mm)	10	Yes
LSB-PWL-2.4-6N	0.094 (2.4 mm)	6	No
LSB-PWL-2.4-10N	0.094 (2.4 mm)	10	No
LSB-PWL-2.4-11N	0.094 (2.4 mm)	11	No

## Instructions for Use Anchor Bolt

### 9.2 DISPOSABLE DRILL KIT

Use the AD-TECH® Disposable Drill Kit to drill a pilot hole in the patient's skull to accept the Anchor Bolt. The following Drill Kits are available:

Catalog Number	Description	Drill			Drill Stop Quantity	Drill Wrench Quantity	Sterile Yes/No
		Quantity	Length (YY)	Diameter (a.a)			
DDK2-2.4-30N	Disposable Drill Kit – Non-sterile	2	30cm	2.4mm	2	2	No
DRL-2.4-22N	Disposable Drill Bit – Non-sterile	1	22cm	2.4mm	0	0	No
DDK2-2.4-16X	Disposable Drill Kit – Sterile	2	16cm	2.4mm	2	2	Yes
DDK2-2.4-30X	Disposable Drill Kit – Sterile	2	30cm	2.4mm	2	2	Yes

### 9.3 RULER

The AD-TECH® Ruler is an optional accessory to the Depth Electrodes and Anchor Bolt. The Ruler is used to make critical measurements when placing the Depth Electrodes.

Catalog #	Description
RULER-20N	20cm Channeled Ruler for drill stop & depth electrode measurement. Non-Sterile, Reusable.

## 10 INSTRUCTIONS FOR USE

	<p><i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories.</i></p> <p><i>Incorrect use of this device may result in severe patient harm.</i></p>
	<p><i>The Anchor Bolt is for surgical use only.</i></p> <p><i>Placement and removal must be in a sterile surgical environment.</i></p> <p><i>Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i></p>

	<p><i>The Anchor Bolt is for short-term use (&lt;30 days) only.</i></p> <p><i>Failure to observe the use conditions may result in severe patient harm.</i></p>
	<p><i>Handle the Anchor Bolts with care to prevent damage.</i></p> <p><i>Do not apply excessive torque to the Anchor Bolt during placement. When using the recommended burr hole diameter, the Anchor Bolt will screw easily into the skull.</i></p> <p><i>Excess torque may result in a bent or broken anchor bolt.</i></p>
	<p><i>Do not use the Anchor Bolts if the third-party interface (e.g., stereotactic/head frame) is damaged or flawed.</i></p>
	<p><i>Discontinue Anchor Bolt placement if you experience any instrument (e.g., the LSB Anchor Bolt Wrench, Cranial Drill, etc.) bending during the placement procedure.</i></p>
	<p><i>Protect the Anchor Bolt from the blunt impact that may occur if the patient has a seizure. Significant force may fracture the Anchor Bolt.</i></p>

Anchor Bolts are placed at the target location using an appropriate stereotactic guidance system and supporting equipment.

When used with the compatible stereotactic system, the AD-TECH® Drill Sleeve Guide supports the alignment of the Cranial Drill Bit and LSB Anchor Bolt Wrench. Refer to the Disposable Cranial Drill Bit Kit IFU for details.

Regarding Figure 1, the steps for placing the Anchor Bolt and associated Depth Electrode are as follows:

**Step 1:** Using the Cranial Drill Kit and a Drill Bit of the appropriate diameter, drill a hole in the patient's skull in the position and trajectory determined by pre-operative planning.

**Step 2:** Remove the Anchor Bolt Cap. Place the Anchor Bolt in the pre-drilled hole and use the LSB Anchor Bolt Wrench in a clockwise direction to screw it into place. The Anchor Bolt has a tapered thread and will naturally tighten when rotated. Do not over-tighten.

**Step 3:** Measure the height of the Anchor Bolt protruding from the skull. Assuming that the surface of the patient's head is the zero reference point, measure from the distal end of the Depth Electrode a distance equal to the desired insertion depth (from the surface of the head) plus the height of the Anchor Bolt, marking the Depth Electrode with indelible ink.

Note: The physician may use the optional AD-TECH® Obturator to create a path through the brain tissue for the Depth Electrode. Refer to the Obturator instructions for use for details.

**Step 4:** Insert the Depth Electrode into the Anchor Bolt Cap. Feed the tip of the Depth Electrode into the Anchor Bolt and push it to the measured depth, aligning the mark on the Depth Electrode body with the top of the Anchor Bolt Cap. Partially tighten the Anchor Bolt Cap and remove the Stylet from the Depth Electrode. Tighten the Anchor Bolt Cap by rotating it clockwise five times.

**Step 5:** Bending the tail end of the Depth Electrode to one side, position the Silicone Cap over the Depth Electrode and Anchor Bolt Cap. The Depth Electrode is now secured.






**Step 6:** Removal of the Anchor Bolt. Remove the Silicone Cap from the Anchor Bolt Cap and Depth Electrode. Untwist the Anchor Bolt Cap anticlockwise and remove it. Slowly withdraw the Depth Electrode from the Anchor Bolt. Use the LSB Anchor Bolt Wrench to unscrew the Anchor Bolt in an anticlockwise direction. Dispose of the Anchor Bolt as directed.

In the event of Anchor Bolt breakage:

1. Stop use immediately and monitor the patient.
2. Attempt to remove the remaining Anchor Bolt using the LSB Anchor Bolt Wrench.
3. Use other surgical equipment, e.g., a clamp, to grip and unscrew the retained portion.
4. Consider other surgical interventions as necessary.

Instructions for Use Anchor Bolt

Figure 1: Electrode Placement with the Anchor Bolt

<p>(1)</p>	
<p>(2)</p>	
<p>(4)</p>	 <p>Place the Depth Electrode to the desired position.</p> <p>(Tighten Partially)</p> <p>Tighten the Anchor Bolt Cap 5 turns.</p>
<p>(5)</p>	
<p>(6)</p>	

## Instructions for Use Anchor Bolt

### 11 MRI SAFETY INFORMATION

	<p><i>The Anchor Bolts are MRI unsafe.</i></p> <p><i>Using the Anchor Bolts in an MRI environment may lead to severe patient harm.</i></p>
---	--

AD-TECH® has not verified the safety and performance of the Anchor Bolts in the MRI environment. AD-TECH® has classified these Anchor Bolts as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice for Marking Medical Devices and Other Items for Safety in The Magnetic Resonance Environment defines MRI Unsafe as:

*A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.*

AD-TECH® has labeled the Anchor Bolts as MRI Unsafe using the symbol:

### 12 PACKAGING, STERILISATION, AND REUSE

Observe the following packaging, sterilization, and reuse instructions:

	<p><i>AD-TECH® sterilizes the Anchor Bolts using ethylene oxide.</i></p>
	<p><i>Do not use the Anchor Bolts if the packaging is damaged or unintentionally opened before use.</i></p> <p><i>Using a device from damaged packing may result in contamination and severe patient harm.</i></p>
	<p><i>Do not re-sterilize the device.</i></p> <p><i>Re-sterilization may result in device failure.</i></p>
	<p><i>The Anchor Bolt is for single use only. Do not reuse the device.</i></p> <p><i>Reuse of the device may result in contamination and severe patient harm.</i></p>
	<p><i>The Anchor Bolt is packaged in a double sterile barrier to support surgical asepsis procedures.</i></p>

### 13 STORAGE CONDITIONS

	<p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p>
---	---

Store the Anchor Bolts in their outer packaging, observing the handling and storage conditions detailed below.

	<p>Fragile, handle with care.</p>
	<p>Keep away from sunlight.</p>

	<p>Keep dry</p>
	<p>Temperature limit</p> <p>Upper = 30°C, Lower = 15°C</p>
	<p>Humidity limitation</p> <p>Upper = 60%, Lower = 20%</p>

### 14 DISPOSAL

	<p><i>Dispose of the Anchor Bolts following the healthcare facility's biohazardous and pathological waste procedures.</i></p>
	<p><i>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</i></p>

### 15 INCIDENT REPORTING

#### 15.1 EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

#### 15.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

#### 15.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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