**INSTRUCTIONS FOR USE
SLOTTED CANNULA**

1 SYMBOL GLOSSARY 2

2 SCOPE 3

3 IMPORTANT NOTES 3

4 DESCRIPTION 3

5 WARNINGS 4

6 PRECAUTIONS 4

7 DEVICE APPLICATION 4

7.1 INTENDED USE/PURPOSE 4

7.2 INDICATIONS FOR USE 4

7.3 CONTRAINDICATIONS 4

7.4 PATIENT TARGET GROUP 4

7.5 INTENDED USER 4

7.6 ENVIRONMENT AND FACILITIES 5

7.7 EXPECTED CLINICAL BENEFIT 5

7.8 UNDESIRABLE SIDE EFFECTS 5

7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE 5

7.10 ELECTRONIC INSTRUCTIONS FOR USE 5

8 PERFORMANCE CHARACTERISTICS 5

8.1 DEVICE ATTRIBUTES 5

9 COMPATIBLE DEVICES/ACCESSORIES 5

9.1 SLOTTED CANNULA SLEEVE GUIDE (SCSG) 5

9.2 DISPOSABLE DRILL KIT 5

9.3 RULER 6

10 INSTRUCTIONS FOR USE 6

11 MRI SAFETY INFORMATION 6

12 PACKAGING, STERILIZATION, AND REUSE 6

12.1 STERILE/NON-STERILE SLOTTED CANNULAS 7

12.2 REUSE 7

13 STORAGE CONDITIONS 7

14 DISPOSAL 7

15 INCIDENT REPORTING 7

15.1 EUROPEAN UNION 7

15.2 CANADA 7



15.3 UNITED STATES 7

16 ADMINISTRATIVE 37

1 SYMBOL GLOSSARY

Symbol	Description
	Manufacturer/Hersteller/Fabricante/Fabricant/ Fabbicante/Fabricante/Üretici
	Authorized representative in the European Community/European Union/Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/Representante autorizado en la Comunidad Europea/Unión Europea/Mandataire dans la Communauté européenne/Union européenne/Mandatario nella Comunità Europea/Unione Europea/Mandatário na Comunidade Europeia/União Europeia/Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi
	Date of manufacture/Herstellungsdatum/Fecha de fabricación/Date de fabrication/Data di fabbricazione/Data de fabrico/Üretim tarihi
	Use-by date/Haltbarkeitsdatum/Fecha de caducidad/Date limite d'utilisation/Data di scadenza/Data de validade/Son kullanma tarihi
	Batch code/Chargenbezeichnung/Código de lote/ Code de lot/Codice lotto/Código de lote/Parti kodu
	Catalog number/Bestellnummer/Número de catálogo/Numéro de référence/Numero di catalogo/Número de catálogo/Katalog numarası
	Non-sterile. Sterilize before use./Nicht verwenden, wenn Verpackung beschädigt ist, und Gebrauchsanweisung beachten./No estéril. Esterilizar antes del uso./Non stérile. Stériliser avant utilisation./Non sterile. Sterilizzare prima dell'uso./Não estéril. Esterilizar antes de utilizar./ Steril değildir. Kullanmadan önce sterilize edin.
	Sterilized using ethylene oxide/Sterilisiert mit Ethylenoxid/Esterilizado con óxido de etileno/ Stérilisé à l'oxyde d'éthylène/Sterilizzato con ossido di etilene/Esterilizado por óxido de etileno/Etilen oksit kullanılarak sterilize edilmiştir
	Do not use if the package is damaged and consult instructions for use/Nicht verwenden, wenn Verpackung beschädigt ist, und Gebrauchsanweisung beachten/No utilizar si el paquete está dañado y consultar las instrucciones de uso/Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation/Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso/ Não utilizar se existirem danos na embalagem e consultar as instruções de utilização/Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
	Double sterile barrier system/Doppel-Sterilbarriere- system/Sistema de barrera estéril doble/Système de double barrière stérile/Sistema a doppia barriera sterile/Sistema de barreira estéril dupla/ Çift steril bariyer sistemi
	Fragile, handle with care/Zerbrechlich, sorgfältig handhaben/Frágil, manipular con cuidado/Fragile, manipuler avec soin/Fragile, maneggiare con cura/Frágil, manusear com cuidado/Kırılgan, dikkatli taşıyın
	Keep dry/Trocken aufbewahren/Mantener seco/ Tenir au sec/Mantenere asciutto/Manter seco/ Kuru tutun

Instructions for Use Slotted Cannula

Symbol	Description
	<p>Temperature limit Upper = 30 °C, Lower = 15 °C</p> <p>Temperaturgrenze obere = 30 °C, untere = 15 °C</p> <p>Límites de temperatura Superior = 30 °C, Inferior = 15 °C</p> <p>Limites de température Supérieure = 30 °C, inférieure = 15 °C</p> <p>Limite di temperatura Superiore = 30 °C, inferiore = 15 °C</p> <p>Limite de temperatura Máximo = 30 °C, Mínimo = 15 °C</p> <p>Sıcaklık sınırı Üst=30 °C, Alt=15 °C</p>
	<p>Humidity limitation Upper = 60%, Lower = 20%</p> <p>Feuchtigkeitsbeschränkung obere = 60 %, untere = 20 %</p> <p>Límites de humedad Superior = 60 %, Inferior = 20 %</p> <p>Limites d'humidité Supérieure = 60 %, inférieure = 20 %</p> <p>Limitazione di umidità Superiore = 60%, inferiore = 20%</p> <p>Limite de humidade Máximo = 60%, Mínimo = 20%</p> <p>Nem sınırları Üst=%60, Alt=%20</p>
	<p>Consult instructions for use/Gebrauchsanweisung beachten/Consultar las instrucciones de uso/ Consulter la notice d'utilisation/Consultare le istruzioni per l'uso/Consultar as instruções de utilização/Kullanma talimatına başvurun</p>
	<p>Caution/Vorsicht/Atención/Mise en garde/ Attenzione/Atenção/Dikkat</p>
	<p>Medical device/Medizinprodukt/Producto sanitario/ Dispositif médical/Dispositivo medico/Dispositivo médico/Tıbbi cihaz</p>
	<p>Unique device identifier/Einmalige Produktkennung (UDI)/Identificador único de dispositivo/Identifiant unique du dispositif/Identificativo unívoco del dispositivo/Identificação única do dispositivo/ Benzersiz cihaz tanımlayıcı</p>
	<p>A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment./Ein Medizinprodukt, das inakzeptable Risiken für den Patienten, medizinisches Personal oder andere Personen in der MR-Umgebung darstellt./Un producto sanitario que plantea riesgos inaceptables para el paciente, el personal médico u otras personas dentro del entorno de RM./Un dispositif médical présentant des risques inacceptables pour le patient, le personnel médical, ou toute autre personne à l'intérieur de l'environnement de RM./ Un dispositivo medico che comporta rischi inaccettabili per il paziente, gli operatori o qualsiasi altro individuo in ambiente RM./ Um dispositivo médico que represente riscos inaceitáveis para o doente, equipa médica ou outras pessoas dentro do ambiente de RM./ MR ortamında hasta, sağlık çalışanları veya diğer kişiler için kabul edilemeyecek riskler taşıyan bir tıbbi cihaz.</p>
	<p>Federal Law (USA) restricts this device to sale by or on the order of a physician./In den USA darf dieses Produkt laut Bundesgesetz nur durch einen Arzt oder auf ärztliche Verordnung verkauft werden./La legislación federal (EE. UU.) restringe</p>

Symbol	Description
	<p>la venta de este dispositivo a médicos o por prescripción facultativa./La loi fédérale (États-Unis) limite la vente de ce dispositif par un médecin ou sur ordonnance médicale./Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo a un medico o su prescrizione di un medico./A lei federal (EUA) restringe a venda deste dispositivo a médicos ou mediante prescrição médica./Federal Yasalar (ABD), bu cihazın satışını hekim tarafından veya talimatıyla yapılacak şekilde sınırlamaktadır.</p>

2 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Slotted Cannula – Non-sterile
- AD-TECH® Slotted Cannula – Sterile

3 IMPORTANT NOTES

	<p><i>Before using these products, carefully review these Instructions for Use.</i></p> <p><i>Keep it readily available for reference.</i></p>
	<p><i>Carefully read all the symbol-marked cautions and warnings.</i></p> <p><i>Incorrect use of the products may lead to severe patient harm.</i></p>

4 DESCRIPTION

The Slotted Cannulas are optional accessories to support the placement of the following Depth Electrodes:

- AD-TECH® Spencer Probe Depth Electrodes, Styles RD, and SD
- AD-TECH® Behnke Fried Depth Electrode (BF)
- AD-TECH® Macro-Micro Depth Electrode (MM)

The Slotted Cannula may be applied when placing Depth Electrodes through a craniotomy or burr hole without using an Anchor Bolt. The Slotted Cannula is applied to create a path within the brain to facilitate the placement of the Depth Electrodes.

The Slotted Cannula is a two-part device consisting of a thin split tube known as the Outer Sheath and an Inner Obturator.

The Slotted Cannula, including the Inner Obturator, is advanced along the planned placement trajectory until reaching the target location.

The surgeon removes the Inner Obturator and inserts the Depth Electrode through the Outer Sheath.



After placement of the Depth Electrode, the surgeon removes the Outer Sheath.

The Slotted Cannulas are available either sterile or non-sterile.

The Slotted Cannulas are reusable up to twenty-five times.

See Figure 1 for details of the available Slotted Cannulas.

Instructions for Use Slotted Cannula

Figure 1 – Available Slotted Cannulas	
REF	2SCK1-190X Contains: 2 Slotted Cannulas (2SC-190X)
REF	2SCK1-240X Contains: 2 Slotted Cannulas (2SC-240X)
REF	2SC-190X Contains: 1 Slotted Cannula
REF	2SC-240X Contains: 1 Slotted Cannula
	
	

5 WARNINGS

	<p><i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories.</i></p> <p><i>Incorrect use of this device may result in severe patient harm.</i></p>
	<p><i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i></p>
	<p><i>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</i></p>
	<p><i>The Slotted Cannula is for surgical use only.</i></p> <p><i>Placement and extraction must be in a sterile surgical environment.</i></p> <p><i>Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i></p>
	<p><i>Do not use the Slotted Cannula if the packaging is damaged or unintentionally opened before use.</i></p> <p><i>Using a device from damaged packing may result in contamination and severe patient harm.</i></p>
	<p><i>Clean and sterilize the Slotted Cannula as directed.</i></p> <p><i>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</i></p>
	<p><i>The Depth Electrode should be introduced easily and NOT FORCED.</i></p> <p><i>Forcing the Depth Electrode into position may result in severe patient harm.</i></p>
	<p><i>The Slotted Cannulas are MRI unsafe.</i></p> <p><i>Using the Slotted Cannulas in an MRI environment may severely harm the user and patient</i></p>

6 PRECAUTIONS

	<p><i>Handle the Slotted Cannulas with care to prevent damage.</i></p>
	<p><i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i></p>
	<p><i>Do not use the Slotted Cannulas if the third-party interface (e.g., stereotactic/head frame) is damaged or flawed.</i></p>
	<p><i>Discontinue using the Slotted Cannula if you experience bending during the placement procedure.</i></p>
	<p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p>
	<p><i>Dispose of the Slotted Cannulas following the healthcare facility's biohazardous and pathological waste procedures.</i></p>
	<p><i>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</i></p>
	<p><i>Federal Law (USA) restricts this device to sale by or on the order of a physician.</i></p>

7 DEVICE APPLICATION

7.1 INTENDED USE/PURPOSE

The Slotted Cannulas are optional accessories for use with Depth Electrodes. The Slotted Cannula may be applied when placing Depth Electrodes through a craniotomy or burr hole without using an Anchor Bolt. The Slotted Cannula is applied to create a path within the brain to facilitate the placement of the Depth Electrode.

7.2 INDICATIONS FOR USE

The Depth Electrodes are indicated for patients undergoing diagnostic electrical brain mapping as part of neurological investigations and treatment planning.

7.3 CONTRAINDICATIONS

	<p><i>The Slotted Cannulas are not for use on any patient whom the physician/surgeon considers at risk of infection.</i></p>
---	--

7.4 PATIENT TARGET GROUP

The Slotted Cannulas are used in pediatric and adult populations, regardless of gender.

7.5 INTENDED USER

The Slotted Cannulas are to be used by physicians/surgeons experienced in Depth Electrodes and the associated Placement Accessories.

Instructions for Use Slotted Cannula

7.6 ENVIRONMENT AND FACILITIES

Slotted Cannulas are for surgical use only.

Placement and removal must be in a sterile surgical environment.

Slotted Cannulas are placed at the target location using an appropriate stereotactic guidance system and supporting equipment.

7.7 EXPECTED CLINICAL BENEFIT

The Slotted Cannula creates a path through the brain to receive the Depth Electrode, reducing tissue trauma and supporting accurate placement.

7.8 UNDESIRABLE SIDE EFFECTS

Use of the Depth Electrodes and the Slotted Cannulas carries the following clinical risks:

- Abscess
- Bleeding
- CSF leakage/drainage/pseudomeningocele
- Death
- Epidural hemorrhage
- Intracerebral hemorrhage
- Medical complications including deep vein thrombosis, pulmonary embolisms, allergic reactions, and psychiatric changes
- Meningitis
- Permanent neurologic deficit
- Subarachnoid hemorrhage
- Subdural hemorrhage
- Superficial wound infection
- Transient neurologic deficit
- Worsening seizures

For further details, refer to the Summary of Safety and Clinical Performance (SSCP)

7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) Unique ID: **SSCP-00002**

The SSCP is available from:

- EUDAMED under SRN: US-MF-000004795 and the Unique ID
- Email request to sales@adtechmedical.com, quoting the SSCP Unique ID

7.10 ELECTRONIC INSTRUCTIONS FOR USE

An electronic copy of this instructions for use is available from:

- Ad-Tech's website www.adtechmedical.com
- Email request to sales@adtechmedical.com, quoting the IFU document number AT10068-B

8 PERFORMANCE CHARACTERISTICS

8.1 DEVICE ATTRIBUTES

The following Slotted Cannulas are available:

Catalog Number	Description	Quantity	Length	Diameter	Supplied Sterile
2SC-190N	One-piece Slotted Cannula	1	190 mm	2.4 mm	No
2SC-190X	One-piece Slotted Cannula	1	190 mm	2.4 mm	Yes

Catalog Number	Description	Quantity	Length	Diameter	Supplied Sterile
2SC-240N	One-piece Slotted Cannula	1	240 mm	2.4 mm	No
2SC-240X	One-piece Slotted Cannula	1	240 mm	2.4 mm	Yes
2SCK1-190X	Two-piece Slotted Cannula	2	190 mm	2.4 mm	Yes
2SCK1-240X	Two-piece Slotted Cannula	2	240 mm	2.4 mm	Yes

9 COMPATIBLE DEVICES/ACCESSORIES

9.1 SLOTTED CANNULA SLEEVE GUIDE (SCSG)

The AD-TECH® Slotted Cannula Sleeve Guide (SCSG) is an accessory to the Slotted Cannula. The SCSG fits compatible stereotactic head frames to provide trajectory alignment and set the penetration depth for the Slotted Cannula. The following SCSGs are available:

Catalog Number	Name	Length	Inner Diameter	Slot Width	Compatible Electrodes	Supplied Sterile Yes/No
SCSG-90-2.18N	Slotted Cannula Sleeve Guide	90 mm	2.18 mm	1.19 mm	SD, RD	No
SCSG2-90-2.18N	Slotted Cannula Sleeve Guide	90 mm	2.18 mm	1.73 mm	LD, BF, MM	No

9.2 DISPOSABLE DRILL KIT

The AD-TECH® Disposable Drill Kits are an optional accessory to the Depth Electrodes. When used with a compatible Surgical Drill, the Drill Kit makes a burr hole to accept the Slotted Cannula and the Depth Electrode. The following Drill Kits are available:

Catalog Number	Description	Drill			Drill Stop Quantity	Drill Wrench Quantity	Sterile Yes/No
		Quantity	Length (YY)	Diameter (a.a)			
DDK2-2.4-30N	Disposable Drill Kit – Non-sterile	2	30 cm	2.4 mm	2	2	No
DRL-2.4-22N	Disposable Drill Bit – Non-sterile	1	22 cm	2.4 mm	0	0	No

Instructions for Use Slotted Cannula

Catalog Number	Description	Drill			Drill Stop Quantity	Drill Wrench Quantity	Sterile Yes/No
		Quantity	Length (YY)	Diameter (a.a)			
DDK2-2.4-16X	Disposable Drill Kit – Sterile	2	16 cm	2.4 mm	2	2	Yes
DDK2-2.4-30X	Disposable Drill Kit – Sterile	2	30 cm	2.4 mm	2	2	Yes
DDK2-2.8-30X	Disposable Drill Kit – Sterile	2	30 cm	2.8 mm	2	2	Yes
DDK2-3.2-30X	Disposable Drill Kit – Sterile	2	30 cm	3.2 mm	2	2	Yes

9.3 RULER

The AD-TECH® Ruler is an optional accessory to the Depth Electrodes and Slotted Cannula. The Ruler is used to make critical measurements when placing the Depth Electrodes.

Catalog #	Description
RULER-20N	20 cm Channeled Ruler for Drill Stop & Depth Electrode measurement. Non-Sterile, Reusable.

10 INSTRUCTIONS FOR USE

	<i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories.</i> <i>Incorrect use of this device may result in severe patient harm.</i>
	<i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i>
	<i>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</i>
	<i>The Slotted Cannula is for surgical use only.</i> <i>Placement and extraction must be in a sterile surgical environment.</i> <i>Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i>
	<i>The Depth Electrode should be introduced easily and NOT FORCED.</i> <i>Forcing the Depth Electrode into position may result in severe patient harm.</i>

	<i>Handle the Slotted Cannulas with care to prevent damage.</i>
	<i>Do not use the Slotted Cannulas if the third-party interface (e.g., stereotactic/head frame) is damaged or flawed.</i>
	<i>Discontinue using the Slotted Cannula if you experience bending during the placement procedure.</i>

Slotted Cannulas are placed at the target location using an appropriate stereotactic guidance system and supporting equipment, including the Slotted Cannula Sleeve Guide (SCSG).

Regarding Figure 2, the steps for placing the Depth Electrode using the SCSG are as follows:

Step 1: Following the Disposable Drill Kit instructions, drill a burr hole in the patient's head to accept the Slotted Cannula.

Step 2: Fit the SCSG to the placement platform (head frame) and confirm the trajectory and depth.

Step 3: Slide the assembled Slotted Cannula into the SCSG, through the burr hole, and into the patient's brain.

Step 4: Remove the Slotted Cannula's Inner Obturator, keeping the Outer Sheath in position.

Step 5: Prepare the Depth Electrode for insertion, measuring and marking the insertion depth on the body of the electrode using the top of the Slotted Cannula as the reference point.

Step 6: Insert the Depth Electrode into the top of the Outer Sheath and push it into the patient's brain, stopping at the calculated depth. The Depth Electrode should be introduced easily and NOT FORCED.

Step 7: Withdraw the Outer Sheath of the Slotted Cannula, guiding the Depth Electrode tail through the slot in the Outer Sheath body.

Step 8: Secure the Depth Electrode to the patient's scalp using the Stay Flange or an alternative method of anchoring.

11 MRI SAFETY INFORMATION

	<i>The Slotted Cannulas are MRI unsafe.</i> <i>Using the Slotted Cannulas in an MRI environment may lead to severe patient harm.</i>
---	---

AD-TECH® has not verified the safety and performance of the Slotted Cannulas in the MRI environment. AD-TECH® has classified the Slotted Cannulas as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices And Other Items For Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

AD-TECH® has labeled the Slotted Cannulas as MRI Unsafe using the symbol:

12 PACKAGING, STERILIZATION, AND REUSE

	<i>Clean and sterilize the Slotted Cannula as directed.</i> <i>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</i>
---	---

Instructions for Use Slotted Cannula

12.1 STERILE/NON-STERILE SLOTTED CANNULAS

Observe the following packaging, sterilization, and reuse instructions:

	<i>AD-TECH® sterilizes the Slotted Cannulas using ethylene oxide.</i>
	<i>Do not use the Slotted Cannula if the packaging is damaged or unintentionally opened before use.</i> <i>Using a device from damaged packing may result in contamination and severe patient harm.</i>
	<i>The Slotted Cannula is packaged in a double sterile barrier to support surgical asepsis procedures.</i>

	Non-sterile. Sterilize before use.
---	------------------------------------

Before and after use, clean and sterilize the Slotted Cannula as per AT10054-B AD-TECH® CLEANING AND STERILIZATION INSTRUCTIONS PLACEMENT TOOLS AND ACCESSORIES.

NOTE: If the device was provided sterile by the manufacturer, the product has already undergone 1 cycle.

12.2 REUSE

	<i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i>
--	--

The Slotted Cannulas can be cleaned, re-sterilized, and reused twenty-five times.

13 STORAGE CONDITIONS

	<i>Observe the recommended storage conditions.</i> <i>Failure to observe storage conditions may result in device failure.</i>
---	--

Store the Slotted Cannulas in their outer packaging, observing the handling and storage conditions detailed below.

	Fragile, handle with care.
	Keep away from sunlight.
	Keep dry
	Temperature limit Upper = 30°C, Lower = 15°C

	Humidity limitation Upper = 60%, Lower = 20%
---	---

14 DISPOSAL

	<i>Dispose of the Slotted Cannulas following the healthcare facility's biohazardous and pathological waste procedures.</i>
	<i>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</i>

15 INCIDENT REPORTING

15.1 EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

15.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

15.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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EC REP Limited 5 Fitzwilliam Square East Dublin 2, D02 R744, Ireland Tel: +353 1 2 544 944