
**INSTRUCTIONS FOR USE
OBTURATORS**


| | | |
|-----------|--|-----------|
| 1 | SYMBOL GLOSSARY | 2 |
| 2 | SCOPE | 3 |
| 3 | IMPORTANT NOTES | 3 |
| 4 | DESCRIPTION | 3 |
| 5 | WARNINGS | 4 |
| 6 | PRECAUTIONS | 4 |
| 7 | DEVICE APPLICATION | 4 |
| 7.1 | INTENDED USE/PURPOSE | 4 |
| 7.2 | INDICATIONS FOR USE..... | 4 |
| 7.3 | CONTRAINDICATIONS | 4 |
| 7.4 | PATIENT TARGET GROUP | 4 |
| 7.5 | INTENDED USER..... | 4 |
| 7.6 | ENVIRONMENT AND FACILITIES | 4 |
| 7.7 | EXPECTED CLINICAL BENEFIT | 4 |
| 7.8 | UNDESIRABLE SIDE EFFECTS | 4 |
| 7.9 | SUMMARY OF SAFETY AND CLINICAL PERFORMANCE | 5 |
| 7.10 | ELECTRONIC INSTRUCTIONS FOR USE | 5 |
| 8 | PERFORMANCE CHARACTERISTICS | 5 |
| 8.1 | DEVICE ATTRIBUTES | 5 |
| 9 | COMPATIBLE DEVICES/ACCESSORIES | 5 |
| 9.1 | RULER | 5 |
| 10 | INSTRUCTIONS FOR USE | 5 |
| 11 | MRI SAFETY INFORMATION | 5 |
| 12 | PACKAGING, STERILIZATION, AND REUSE | 6 |
| 12.1 | STERILE/NON-STERILE OBTURATORS | 6 |
| 12.2 | REUSE..... | 6 |
| 13 | STORAGE CONDITIONS | 6 |
| 14 | DISPOSAL | 6 |
| 15 | INCIDENT REPORTING | 6 |
| 15.1 | EUROPEAN UNION | 6 |
| 15.2 | CANADA | 6 |
| 15.3 | UNITED STATES..... | 6 |
| 16 | ADMINISTRATIVE | 36 |

1 SYMBOL GLOSSARY

| Symbol | Description |
|---|--|
|  | Manufacturer/Hersteller/Fabricante/Fabricant/ Fabbricante/Fabricante/Üretici |
|  | Authorized representative in the European Community/European Union/Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/ Representante autorizado en la Comunidad Europea/Unión Europea/Mandataire dans la Communauté européenne/Union européenne/ Mandatario nella Comunità Europea/Unione Europea/Mandatário na Comunidade Europeia/ União Europeia/Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi |
|  | Date of manufacture/Herstellungsdatum/Fecha de fabricación/Date de fabrication/Data di fabbricazione/ Data de fabrico/Üretim tarihi |
|  | Use-by date/Haltbarkeitsdatum/Fecha de caducidad/ Date limite d'utilisation/Data di scadenza/Data de validade/Son kullanma tarihi |
|  | Batch code/Chargenbezeichnung/Código de lote/ Code de lot/Codice lotto/Código de lote/Parti kodu |
|  | Catalog number/Bestellnummer/Número de catálogo/ Número de référence/Numero di catalogo/Número de catálogo/Katalog numarası |
|  | Non-sterile. Sterilize before use./Unsteril. Vor Gebrauch sterilisieren./No estéril. Esterilizar antes del uso./Non stérile. Stériliser avant utilisation./Non sterile. Sterilizzare prima dell'uso./Não estéril. Esterilizar antes de utilizar./Steril değildir. Kullanmadan önce sterilize edin. |
|  | Sterilized using ethylene oxide/Sterilisiert mit Ethylenoxid/Esterilizado con óxido de etileno/Stérilisé à l'oxyde d'éthylène/Sterilizzato con ossido di etilene/Esterilizado por óxido de etileno/Etilen oksit kullanılarak sterilize edilmiştir |
|  | Do not use if the package is damaged and consult instructions for use/Nicht verwenden, wenn Verpackung beschädigt ist und Gebrauchsanweisung beachten/No utilizar si el paquete está dañado y consultar las instrucciones de uso/Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation/Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso/Não utilizar se existirem danos na embalagem e consultar as instruções de utilização/Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun |
|  | Double sterile barrier system/Doppel-Sterilbarriere- system/Sistema de barrera estéril doble/Système de double barrière stérile/Sistema a doppia barriera sterile/Sistema de barreira estéril dupla/Çift steril bariyer sistemi |
|  | Fragile, handle with care/Zerbrechlich, sorgfältig handhaben/Frágil, manipular con cuidado/Fragile, manipuler avec soin/Fragile, maneggiare con cura/ Frágil, manusear com cuidado/Kırılgan, dikkatli taşıyın |
|  | Keep dry/Trocken aufbewahren/Mantener seco/Tenir au sec/Mantener asciutto/Manter seco/Kuru tutun |

Instructions for Use Obturators

| Symbol | Description |
|---|---|
|  | <p>Temperature limit Upper = 30 °C, Lower = 15 °C Temperaturgrenze obere = 30 °C, untere = 15 °C Límites de temperatura Superior = 30 °C, Inferior = 15 °C Limites de température Supérieure = 30 °C, inférieure = 15 °C Limite di temperatura Superiore = 30 °C, inferiore = 15 °C Limite de temperatura Máximo = 30 °C, Mínimo = 15 °C Sıcaklık sınırı Üst=30 °C, Alt=15 °C</p> |
|  | <p>Humidity limitation Upper = 60%, Lower = 20% Feuchtigkeitsbeschränkung obere = 60 %, untere = 20 % Límites de humedad Superior = 60 %, Inferior = 20 % Limites d'humidité Supérieure = 60 %, inférieure = 20 % Limitazione di umidità Superiore = 60%, inferiore = 20% Limite de humidade Máximo = 60%, Mínimo = 20% Nem sınırları Üst=%60, Alt=%20</p> |
|  | <p>Consult instructions for use/Gebrauchsanweisung beachten/Consultar las instrucciones de uso/Consultar la notice d'utilisation/Consultare le istruzioni per l'uso/ Consultar as instruções de utilização/Kullanma talimatına başvurun</p> |
|  | <p>Caution/Vorsicht/Atención/Mise en garde/Attenzione/Atenção/Dikkat</p> |
|  | <p>Medical device/Medizinprodukt/Producto sanitario/ Dispositif médical/Dispositivo medico/Dispositivo médico/Tıbbi cihaz</p> |
|  | <p>Unique device identifier/Einmalige Produktkennung (UDI)/Identificador único de dispositivo/Identifiant unique du dispositif/Identificativo univoco del dispositivo/Identificação única do dispositivo/ Benzersiz cihaz tanımlayıcı</p> |
|  | <p>A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment./Ein Medizinprodukt, das inakzeptable Risiken für den Patienten, medizinisches Personal oder andere Personen in der MR-Umgebung darstellt./Un producto sanitario que plantea riesgos inaceptables para el paciente, el personal médico u otras personas dentro del entorno de RM./Un dispositif médical présentant des risques inacceptables pour le patient, le personnel médical, ou toute autre personne à l'intérieur de l'environnement de RM./Un dispositivo medico che comporta rischi inaccettabili per il paziente, gli operatori o qualsiasi altro individuo in ambiente RM./Um dispositivo médico que represente riscos inaceitáveis para o doente, equipa médica ou outras pessoas dentro do ambiente de RM./MR ortamında hasta, sağlık çalışanları veya diğer kişiler için kabul edilemeyecek riskler taşıyan bir tıbbi cihaz.</p> |

| Symbol | Description |
|---|--|
|  | <p>Federal Law (USA) restricts this device to sale by or on the order of a physician./In den USA darf dieses Produkt laut Bundesgesetz nur durch einen Arzt oder auf ärztliche Verordnung verkauft werden./La legislación federal (EE. UU.) restringe la venta de este dispositivo a médicos o por prescripción facultativa./La loi fédérale (États-Unis) limite la vente de ce dispositif par un médecin ou sur ordonnance médicale./Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo a un medico o su prescrizione di un medico./A lei federal (EUA) restringe a venda deste dispositivo a médicos ou mediante prescrição médica./Federal Yasalar (ABD), bu cihazın satışını hekim tarafından veya talimatıyla yapılacak şekilde sınırlamaktadır.</p> |

2 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Obturators – Non-sterile
- AD-TECH® Obturators – Sterile

3 IMPORTANT NOTES

| | |
|---|---|
|  | <p>Before using these products, carefully review these Instructions for Use. Keep it readily available for reference.</p> |
|  | <p>Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.</p> |

4 DESCRIPTION

The Obturators are optional accessories to support the placement of the following Depth Electrodes:

- AD-TECH® Spencer Probe Depth Electrodes, Styles LD, RD, and SD
- AD-TECH® Behnke Fried Depth Electrode (BF)
- AD-TECH® Macro-Micro Depth Electrode (MM)

An Obturator is a solid stainless steel rod with a rounded tip and a Delrin® hub.

The Obturator is used in conjunction with an Anchor Bolt. The Obturator is advanced through the Anchor Bolt along the planned placement trajectory to create a pathway in the brain.


The Obturator is withdrawn, and the Depth Electrode is advanced in the same manner along the path made by the Obturator until the Depth Electrode reaches the target location.

Obturator selection is dependent on the type of Depth Electrode.

The Obturators are available either sterile or non-sterile.

The Obturators are reusable up to twenty-five times.

See Figure 1 for details of the available Obturators.

| Figure 1 – Available Obturators | | |
|--|------------|-----------------------|
| REF | OB-20-190N | Contains: 1 Obturator |
| REF | OB-17-240N | Contains: 1 Obturator |
|  | | |

Instructions for Use Obturators

5 WARNINGS

| | |
|---|---|
|  | <p><i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories.</i></p> <p><i>Incorrect use of this device may result in severe patient harm.</i></p> |
|  | <p><i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i></p> |
|  | <p><i>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</i></p> |
|  | <p><i>The Obturator is for surgical use only.</i></p> <p><i>Placement and extraction must be in a sterile surgical environment.</i></p> <p><i>Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i></p> |
|  | <p><i>Do not use the Obturator if the packaging is damaged or unintentionally opened before use.</i></p> <p><i>Using a device from damaged packing may result in contamination and severe patient harm.</i></p> |
|  | <p><i>Clean and sterilize the Obturator as directed.</i></p> <p><i>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</i></p> |
|  | <p><i>The Depth Electrode should be introduced easily and NOT FORCED.</i></p> <p><i>Forcing the Depth Electrode into position may result in severe patient harm.</i></p> |
|  | <p><i>The Obturators are MRI unsafe.</i></p> <p><i>Using the Obturators in an MRI environment may severely harm the user and patient</i></p> |

6 PRECAUTIONS

| | |
|---|---|
|  | <p><i>Handle the Obturators with care to prevent damage.</i></p> |
|  | <p><i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i></p> |
|  | <p><i>Do not use the Obturators if the third-party interface (e.g., stereotactic/head frame) is damaged or flawed.</i></p> |
|  | <p><i>Discontinue using the Obturator if you experience bending during the placement procedure.</i></p> |
|  | <p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p> |
|  | <p><i>Dispose of the Obturators following the healthcare facility's biohazardous and pathological waste procedures.</i></p> |

| | |
|---|--|
|  | <p><i>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</i></p> |
|  | <p><i>Federal Law (USA) restricts this device to sale by or on the order of a physician.</i></p> |

7 DEVICE APPLICATION

7.1 INTENDED USE/PURPOSE

The Obturators are optional accessories for use with Depth Electrodes. The Obturator may be applied before the placement of the Depth Electrodes. The Obturator is applied to create a path within the brain to facilitate the placement of the Depth Electrode.

7.2 INDICATIONS FOR USE

The Depth Electrodes are indicated for patients undergoing diagnostic electrical brain mapping as part of neurological investigations and treatment planning.

7.3 CONTRAINDICATIONS

| | |
|---|--|
|  | <p><i>The Obturators are not for use on any patient whom the physician/surgeon considers at risk of infection.</i></p> |
|  | <p><i>The Obturator must not be used for the application of thermal energy.</i></p> |

7.4 PATIENT TARGET GROUP

The Obturators are used in pediatric and adult populations, regardless of gender.

7.5 INTENDED USER

The Obturators are to be used by physicians/surgeons experienced in Depth Electrodes and the associated Placement Accessories.

7.6 ENVIRONMENT AND FACILITIES

Obturators are for surgical use only.

Placement and removal must be in a sterile surgical environment.

The Obturators are placed at the target location using an appropriate stereotactic guidance system and supporting equipment.

7.7 EXPECTED CLINICAL BENEFIT

The Obturator creates a path through the brain to receive the Depth Electrode, reducing tissue trauma and supporting accurate placement.

7.8 UNDESIRABLE SIDE EFFECTS

Use of the Depth Electrodes and the Obturators carries the following clinical risks:

- Abscess
- Bleeding
- CSF leakage/drainage/pseudomeningocele
- Death
- Epidural hemorrhage

Instructions for Use Obturators

- Intracerebral hemorrhage
- Medical complications including deep vein thrombosis, pulmonary embolisms, allergic reactions, and psychiatric changes
- Meningitis
- Permanent neurologic deficit
- Subarachnoid hemorrhage
- Subdural hemorrhage
- Superficial wound infection
- Transient neurologic deficit
- Worsening seizures

For further details, refer to the Summary of Safety and Clinical Performance (SSCP)

7.9 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) Unique ID: **SSCP-00002**

The SSCP is available from:

- EUDAMED under SRN: US-MF-000004795 and the Unique ID
- Email request to sales@adtechmedical.com, quoting the SSCP Unique ID

7.10 ELECTRONIC INSTRUCTIONS FOR USE

An electronic copy of this instructions for use is available from:

- Ad-Tech's website www.adtechmedical.com
- Email request to sales@adtechmedical.com, quoting the IFU document number AT10070-B

8 PERFORMANCE CHARACTERISTICS

8.1 DEVICE ATTRIBUTES

The following Obturators are available:

| Catalog Number | Compatible Electrodes | Length (L) | Diameter | Supplied Sterile |
|----------------|-----------------------|------------|----------|------------------|
| OB-20-190X | LD, RD, SD, BF, MM | 190 mm | 0.86 mm | Yes |
| OB-20-190N | LD, RD, SD, BF, MM | 190 mm | 0.86 mm | No |
| OB-17-240X | LD, SD, BF, MM | 240 mm | 1.17 mm | Yes |
| OB-17-240N | LD, SD, BF, MM | 240 mm | 1.17 mm | No |

9 COMPATIBLE DEVICES/ACCESSORIES

9.1 RULER

The AD-TECH® Ruler is an optional accessory to the Depth Electrodes and Obturator. The Ruler is used to make critical measurements when placing the Depth Electrodes.

| Catalog # | Description |
|-----------|---|
| RULER-20N | 20cm Channeled Ruler for Drill Stop & Depth Electrode measurement. Non-Sterile, Reusable. |

10 INSTRUCTIONS FOR USE

| | |
|---|--|
|  | <i>This product should only be used by a physician/surgeon experienced in the use of Depth Electrodes and the associated Placement Accessories. Incorrect use of this device may result in severe patient harm.</i> |
|  | <i>Placement of the Depth Electrodes without the correct pre-operative planning may result in severe patient harm.</i> |
|  | <i>Placement of the Depth Electrodes without the correct stereotactic guidance may result in severe patient harm.</i> |
|  | <i>The Obturator is for surgical use only. Placement and extraction must be in a sterile surgical environment. Failure to observe surgical asepsis and sterile techniques may result in infection and severe patient harm.</i> |
|  | <i>The Depth Electrode should be introduced easily and NOT FORCED. Forcing the Depth Electrode into position may result in severe patient harm.</i> |
|  | <i>Handle the Obturators with care to prevent damage.</i> |
|  | <i>Do not use the Obturators if the third-party interface (e.g., stereotactic/head frame) is damaged or flawed.</i> |
|  | <i>Discontinue using the Obturator if you experience bending during the placement procedure.</i> |

The Obturator is used in conjunction with an Anchor Bolt. For placement of the Anchor Bolt, follow the Anchor Bolt instructions for use.

Regarding Figure 2, the steps for deploying the Obturator are as follows:

Step 1: Prepare the Depth Electrode for insertion, measuring and marking the insertion depth on the body of the electrode using the top of the Anchor Bolt as the reference point.

Step 2: Prepare the Obturator for insertion, measuring and marking the insertion depth on the body of the Obturator using the top of the Anchor Bolt as the reference point.

Step 3: Slide the Obturator into the Anchor Bolt and push it to the desired depth.

Step 4: Remove the Obturator and discard. The Anchor Bolt is now ready to receive the Depth Electrode.

Refer to the Depth Electrode IFUs for further information on their placement procedure.

11 MRI SAFETY INFORMATION

| | |
|---|---|
|  | <i>The Obturators are MRI unsafe. Using the Obturators in an MRI environment may lead to severe patient harm.</i> |
|---|---|

AD-TECH® has not verified the safety and performance of the Obturators in the MRI environment.

AD-TECH® has classified the Obturators as MRI Unsafe and labelled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices And Other Items For Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

Instructions for Use Obturators

AD-TECH® has labelled the Obturators as MRI Unsafe using the symbol:

12 PACKAGING, STERILIZATION, AND REUSE

| | |
|--|---|
| | <p><i>Clean and sterilize the Obturators as directed.</i></p> <p><i>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</i></p> |
|--|---|

12.1 STERILE/NON-STERILE OBTURATORS

Observe the following packaging, sterilization, and reuse instructions:

| | |
|--|---|
| | AD-TECH® sterilizes the Obturators using ethylene oxide. |
| | <p><i>Do not use the Obturator if the packaging is damaged or unintentionally opened before use.</i></p> <p><i>Using a device from damaged packing may result in contamination and severe patient harm.</i></p> |
| | <i>The Obturator is packaged in a double sterile barrier to support surgical asepsis procedures.</i> |

| | |
|--|------------------------------------|
| | Non-sterile. Sterilize before use. |
|--|------------------------------------|

Before and after use, clean and sterilize the Obturator as per AT10054-B AD-TECH® CLEANING AND STERILIZATION INSTRUCTIONS PLACEMENT TOOLS AND ACCESSORIES.

NOTE: If the device was provided sterile by the manufacturer, the product has already undergone 1 cycle.

12.2 REUSE

| | |
|--|--|
| | <i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i> |
|--|--|

The Obturators can be cleaned, re-sterilized, and reused twenty-five times.

13 STORAGE CONDITIONS

| | |
|--|---|
| | <p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p> |
|--|---|

Store the Obturators in their outer packaging, observing the handling and storage conditions detailed below.

| | |
|--|----------------------------|
| | Fragile, handle with care. |
| | Keep away from sunlight. |

| | |
|--|--|
| | Keep dry |
| | <p>Temperature limit</p> <p>Upper = 30°C, Lower = 15°C</p> |
| | <p>Humidity limitation</p> <p>Upper = 60%, Lower = 20%</p> |

14 DISPOSAL

| | |
|--|---|
| | <i>Dispose of the Obturators following the healthcare facility's biohazardous and pathological waste procedures.</i> |
| | <i>Dispose of any devices/accessories that pose the risk of needle sticks, cuts, or puncture injuries in a sharps disposal container.</i> |

15 INCIDENT REPORTING

15.1 EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

15.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

15.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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