INSTRUCTIONS FOR USE

Slotted Cannula Sleeve Guides

SYMBOL GLOSSARY

Symbol Símbolo Symbole Simbolo Símbolo Sembol	Description/ Beschreibung/ Descripción/ Description/ Descrizione/ Descrição/ Açıklama
	Manufacturer/ Hersteller/ Fabricante/ Fabricant/ Fabbicante/ Fabricante/ Üretici
	Authorized representative in the European Community/European Union/ Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/ Representante autorizado en la Comunidad Europea/Unión Europea/ Mandataire dans la Communauté européenne/Union européenne/ Mandatario nella Comunità Europea/Unione Europea/ Mandatário na Comunidade Europeia/União Europeia/ Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi
	Use-by date/ Haltbarkeitsdatum/ Fecha de caducidad/ Date limite d'utilisation/ Data di scadenza/ Data de validade/ Son kullanma tarihi
	Batch code/ Chargenbezeichnung/ Código de lote/ Code de lot/ Codice lotto/ Código de lote/ Parti kodu
	Catalog number/ Bestellnummer/ Número de catálogo/ Numéro de référence/ Numero di catalogo/ Número de catálogo/ Katalog numarası
	Non-sterile. Sterilize before use./ Unsteril. Vor Gebrauch sterilisieren./ No estéril. Esterilizar antes del uso./ Non stérile. Stériliser avant utilisation./ Non sterile. Sterilizzare prima dell'uso./ Não estéril. Esterilizar antes de utilizar./ Steril değildir. Kullanmadan önce sterilize edin.
	Sterilized using ethylene oxide/ Sterilisiert mit Ethylenoxid/ Esterilizado con óxido de etileno/ Stérilisé à l'oxyde d'éthylène/ Sterilizzato con ossido di etilene/ Esterilizado por óxido de etileno/ Etilen oksit kullanılarak sterilize edilmiştir
	Do not use if the package is damaged and consult instructions for use/ Nicht verwenden, wenn Verpackung beschädigt ist und Gebrauchsanweisung beachten/ No utilizar si el envase está dañado y consultar instrucciones de uso/ Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation/ Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso/ Não utilizar se existirem danos na embalagem e consultar as instruções de utilização/ Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
	Do not re-sterilize/ Nicht reesterilisieren/ No reesterilizar/ Ne pas resteriliser/ Non risterilizzare/ Não reesterilizar/ Tekrar sterilize etmeyin
	Do not reuse/ Nicht wiederverwenden/ No reutilizar/ Ne pas réutiliser/ Non riutilizzare/ Não reutilizar/ Tekrar kullanmayın
	Double sterile barrier system/ Doppel Sterilbarriere-System/ Sistema de barrera estéril doble/ Système de double barrière stérile/ Sistema a doppia barriera sterile/ Sistema de barreira estéril dupla/ Çift steril bariyer sistemi
	Fragile, handle with care/ Zerbrechlich, sorgfältig handhaben/ Frágil, manipular con cuidado/ Fragile, manipuler avec soin/ Fragile, maneggiare con cura/ Frágil, manusear com cuidado/ Kırılgan, dikkatli taşıyın
	Keep dry / Trocken aufbewahren/ Mantener seco/ Tenir au sec/ Mantenere asciutto/ Manter seco/ Kuru tutun

Table of Contents

1	SCOPE	3
2	IMPORTANT NOTES	3

	Temperature limit/ Temperaturgrenze/ Limite de temperatura/ Limites de température/ Limite di temperatura/ Limite de temperatura/ Sıcaklık sınırı Upper/ obere / Superior/ Supérieure / Superiore/ Máximo/ Üst = 30°C, Lower/ untere/ Inferior/ inférieure/ inferiore/ Mínimo/ Alt = 15°C
	Humidity limitation/ Feuchtigkeitsbeschränkung/ Limite de humedad/ Limites d'humidité/ Limitazione di umidità/ Limite de humidade/ Nem sınırları Upper/ obere/ Superior/ Supérieure / Superiore/ Máximo/ Üst = 60%, Lower/ untere / Inferior/ inférieure/ inferior/ Mínimo/ Alt = 20%
	Consult instructions for use / Gebrauchsanweisung beachten/ Consultar las instrucciones de uso/ Consulter la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consultar as instruções de utilização/ Kullanma talimatına başvurun
	Caution/ Vorsicht/ Atención/ Mise en garde/ Attenzione/ Atenção/ Dikkat
	Medical device/ Medizinprodukt/ Producto sanitario/ Dispositif médical/ Dispositivo medico/ Dispositivo médico/ Tıbbi cihaz
	Unique device identifier/ Einmalige Produktkennung (UDI)/ Identificador único de dispositivo/ Identifiant unique du dispositif/ Identificativo univoco del dispositivo/ Identificação única do dispositivo/ Benzersiz cihaz tanımlayıcı
	A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment./ Ein Medizinprodukt, das inakzeptable Risiken für den Patienten, medizinisches Personal oder andere Personen in der MR- Umgebung darstellt./ Un producto sanitario que plantea riesgos inaceptables para el paciente, el personal médico u otras personas dentro del entorno de RM./ Un dispositif médical présentant des risques inacceptables pour le patient, le personnel médical, ou toute autre personne à l'intérieur de l'environnement de RM./ Un dispositivo medico che comporta rischi inaccettabili per il paziente, gli operatori o qualsiasi altro individuo in ambiente RM./ Um dispositivo médico que represente riscos inaceitáveis para o doente, equipa médica ou outras pessoas dentro do ambiente de RM./ MR ortamında hasta, sağlık çalışanları veya diğer kişiler için kabul edilemeyecek riskler taşıyan bir tıbbi cihaz.
	Federal Law (USA) restricts this device to sale by or on the order of a physician. / In den USA darf dieses Produkt laut Bundesgesetz nur durch einen Arzt oder auf ärztliche Verordnung verkauft werden./ La legislación federal (EE. UU.) restringe la venta de este dispositivo a médicos o por prescripción facultativa./ La loi fédérale (États-Unis) limite la vente de ce dispositif par un médecin ou sur ordonnance médicale./ Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo a un medico o su prescrizione di un medico./ A lei federal (EUA) restringe a venda deste dispositivo a médicos ou mediante prescrição médica./ Federal Yasalar (ABD), bu cihazın satışını hekim tarafından veya talimatıyla yapılacak şekilde sınırlamaktadır.

3	DESCRIPTION	3
4	WARNINGS	3
5	PRECAUTIONS	3
6	DEVICE APPLICATION	3
6.1	INTENDED USE/PURPOSE	3
6.2	INDICATIONS FOR USE	3
6.3	CONTRAINDICATIONS	3
6.4	PATIENT TARGET GROUP	3
6.5	INTENDED USER	3
6.6	ENVIRONMENT AND FACILITIES	3
6.7	EXPECTED CLINICAL BENEFIT	3

6.8	UNDESIRABLE SIDE EFFECTS	3
7	PERFORMANCE CHARACTERISTICS.....	4
7.1	DEVICE ATTRIBUTES.....	4
8	INSTRUCTIONS FOR USE.....	4
9	MRI SAFETY INFORMATION	4
10.1	MRI SAFETY, Slotted Cannula Sleeve Guides (SCSG)	4
10	PACKAGING, STERILIZATION, AND REUSE.....	4
11.1	NON-STERILE Slotted Cannula Sleeve Guides (SCSG)	4
11.2	REUSE	4
11	STORAGE CONDITIONS	4
12	DISPOSAL.....	5
13	INCIDENT REPORTING	5
13.1	EUROPEAN UNION.....	5
13.2	CANADA.....	5
13.3	UNITED STATES	5

1 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Slotted Cannula Sleeve Guide

2 IMPORTANT NOTES

	Before using these products, carefully review these Instructions for Use. Keep it readily available for reference.
	Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.


3 DESCRIPTION

The AD-TECH® Slotted Cannula Sleeve Guide (SCSG) is an accessory to the Slotted Cannula. The SCSG fits compatible stereotactic head frames to provide trajectory alignment and set the penetration depth for the Slotted Cannula.

The SCSG is available non-sterile. End-user sterilized by steam (Auto-clave), STERRAD® (Vaporized Hydrogen Peroxide Gas Plasma Sterilization), or Ethylene Oxide (EO).

The SCSGs are reusable up to twenty-five times.

See Figure 1 for details of the available Slotted Cannula Sleeve Guide (SCSG).

4 WARNINGS

No known warnings for the Slotted Cannula Sleeve Guide.

5 PRECAUTIONS

	Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.
	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
	Federal Law (USA) restricts this device to sale by or on the order of a physician.

6 DEVICE APPLICATION

6.1 INTENDED USE/PURPOSE

The Slotted Cannula Sleeve Guides are accessories to the Ad-Tech Slotted Cannulas.

The Slotted Cannulas are optional accessories for use with Depth Electrodes. The Slotted Cannulas may be applied when it is desired to place Ad-Tech's Depth Electrodes through a craniotomy or burr hole and without the use of an anchor bolt. The Slotted Cannula is applied to create a path within the brain to facilitate the placement of Ad-Tech's Depth Electrodes.

6.2 INDICATIONS FOR USE

There are no specific indications; this is an accessory to the Slotted Cannulas.

6.3 CONTRAINDICATIONS

None known.

6.4 PATIENT TARGET GROUP

Applicable to the Depth Electrodes and the Placement Accessories.

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

6.5 INTENDED USER

The SCSGs are to be used by physicians/surgeons experienced in the use of Slotted Cannula Sleeve Guides.

6.6 ENVIRONMENT AND FACILITIES

SCSGs are for surgical use only.

Use of the device must be in a sterile surgical environment.

6.7 EXPECTED CLINICAL BENEFIT

The SCSG is an accessory to the Slotted Cannula. The SCSG fits compatible stereotactic head frames to provide trajectory alignment and set the penetration depth for the Slotted Cannula.

6.8 UNDESIRABLE SIDE EFFECTS

No clinical risks with the use of the SCSG.

7 PERFORMANCE CHARACTERISTICS

7.1 DEVICE ATTRIBUTES

The SCSG is an accessory to the Slotted Cannula. The SCSG fits compatible stereotactic head frames to provide trajectory alignment and set the penetration depth for the Slotted Cannula.

Catalog Number	Name	Length (YY)	Inner Diameter (A.AA)	Slot Width	Compatible Electrodes	Supplied Sterile Yes/No
SCSG-90-2.18N	Slotted Cannula Sleeve Guide	90 mm	2.18mm	1.19 mm	SD, RD, SR	No
SCSG2-90-2.18N	Slotted Cannula Sleeve Guide	90 mm	2.18mm	1.73 mm	BF, MM	No

8 INSTRUCTIONS FOR USE

Slotted Cannulas are placed at the target location using an appropriate stereotactic guidance system and supporting equipment, including the Slotted Cannula Sleeve Guide (SCSG).

The steps for placing the Depth Electrode using the SCSG are as follows (see Figure 2):

Step 1a: Slide the assembled Slotted Cannula into the SCSG, through the burr hole, and into the patient's brain.


Step 2: Remove the Slotted Cannula's Inner Obturator, keeping the Outer Sheath in position.

Step 3: Prepare the Depth Electrode for insertion, measuring and marking the insertion depth on the body of the electrode using the top of the Slotted Cannula as the reference point.

Step 3a: Insert the Depth Electrode into the top of the Outer Sheath and push it into the patient's brain, stopping at the calculated depth. The Depth Electrode should be introduced easily and NOT FORCED.

Step 3a (cont'd): Withdraw the Outer Sheath of the Slotted Cannula, guiding the Depth Electrode tail through the slot in the Outer Sheath body.

Figure 2:

9 MRI SAFETY INFORMATION

	The SCSGs are MRI unsafe. Using the SCSGs in an MRI environment may lead to severe patient harm.
--	--

10.1 MRI SAFETY, Slotted Cannula Sleeve Guides (SCSG)

AD-TECH® has not verified the safety and performance of the SCSGs in the MRI environment.

AD-TECH® has classified these sleeve guides as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices and Other Items for Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

AD-TECH® has labeled the SCSGs as MRI Unsafe using the symbol:

10 PACKAGING, STERILIZATION, AND REUSE

	Clean and sterilize the SCSGs as directed. Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.
--	--

11.1 NON-STERILE Slotted Cannula Sleeve Guides (SCSG)

	Non-sterile. Sterilize before use.
--	---

Before and after use, clean and sterilize the SCSGs as per AT10054-B AD-TECH® CLEANING AND STERILIZATION INSTRUCTIONS PLACEMENT TOOLS AND ACCESSORIES.

11.2 REUSE

	Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.
--	--

The SCSGs can be cleaned and re-sterilized for use, up to a total of twenty-five times.

11 STORAGE CONDITIONS

	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
--	---

Store the SCSGs in their outer packaging, observing the handling and storage conditions detailed below.

	Fragile, handle with care.
	Keep dry
	Temperature limit Upper = 30°C, Lower = 15°C
	Humidity limitation Upper = 60%, Lower = 20%

12 DISPOSAL

	Dispose of the SCSGs following the healthcare facility's biohazardous and pathological waste procedures.
--	--

13 INCIDENT REPORTING**13.1 EUROPEAN UNION**

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

13.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

13.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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EC REP Limited 5 Fitzwilliam Square East Dublin 2, D02 R744, Irlanda Tel.: +353 1 2 544 944