INSTRUCTIONS FOR USE

Signum Connection Systems

SYMBOL GLOSSARY

Symbol	Description
	Manufacturer/Hersteller/Fabricante/Fabricant/Fabbricante/Fabricante/Üretici
	Authorized representative in the European Community/European Union/Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/Representante autorizado en la Comunidad Europea/Unión Europea/Mandataire dans la Communauté européenne/Unione europea/Mandatario nella Comunità Europea/Unione Europea/Mandatário na Comunidade Europeia/União Europeia/Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi
	Date of manufacture/Herstellungsdatum/Fecha de fabricación/Data de fabrication/Data di fabbricazione/Data de fabrico/Üretim tarihi
	Use-by date/Haltbarkeitsdatum/Fecha de caducidad/Date limite d'utilisation/Data di scadenza/Data de validade/Son kullanma tarihi
	Batch code/Chargenbezeichnung/Código de lote/Code de lot/Codice lotto/Código de lote/Parti kodu
	Catalog number/Bestellnummer/Número de catálogo/Número de reference/Numero di catalogo/Número de catálogo/Katalog numarası
	Non-sterile. Sterilize before use./Unsteril. Vor Gebrauch sterilisieren./No estéril. Esterilizar antes del uso./Non stérile. Stériliser avant utilisation./Non sterile. Sterilizzare prima dell'uso./Não esteril. Esterilizar antes de utilizar./Steril değildir. Kullanmadan önce sterilize edin.
	Do not use if the package is damaged and consult instructions for use./Nicht verwenden, wenn Verpackung beschädigt ist, und Gebrauchsanweisung beachten./No utilizar si el paquete está dañado y consultar las instrucciones de uso./Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation./Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso./Não utilizar se existirem danos na embalagem e consultar as instruções de utilização./Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
	Fragile, handle with care/Zerbrechlich, sorgfältig handhaben/Frágil, manipular con cuidado/Fragile, manipuler avec soin/Frágile, maneggiare con cura/Frágil, manusear com cuidado/Kırılğan, dikkatli taşıyın
	Keep dry/Trocken aufbewahren/Mantenere secco/Tenir au sec/Mantenere asciutto/Manter seco/Kurutu tutun
	Temperature limit Upper = 30 °C, Lower = 15 °C Temperaturgrenze Obere = 30 °C, untere = 15 °C Limites de temperatura Superior = 30 °C, Inferior = 15 °C Limites de température Supérieure = 30 °C, inférieure = 15 °C Limite di temperatura Superiore = 30 °C, inferiore = 15 °C Limite de temperatura Máximo = 30 °C, Mínimo = 15 °C Sıcaklık sınırı Üst = 30 °C, Alt = 15 °C
	Humidity limitation Upper = 60%, Lower = 20% Feuchtigkeitsbeschränkung Obere = 60 %, untere = 20 % Limites de humedad Superior = 60 %, Inferior = 20 % Limites d'humidité Supérieure = 60 %, inférieure = 20 % Limitazione di umidità Superiore = 60%, inferiore = 20% Limite de humidade Máximo = 60%, Mínimo = 20% Nem sınırları Üst = %60, Alt = %20

	Consult instructions for use or consult electronic instructions for use/Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten/Consulte las instrucciones de uso o las instrucciones de uso electrónicas/Consultar la notice d'utilisation ou la notice d'utilisation électronique/Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche/Consultar as instruções de utilização em papel ou eletrónicas/Kullanma talimatına veya elektronik kullanma talimatına başvurun
	Caution/Vorsicht/Atención/Mise en garde/Attenzione/Atenção/Dikkat
	Medical device/Medizinprodukt/Producto sanitario/Dispositif medical/Dispositivo medico/Dispositivo medico/Tıbbi cihaz
	Unique device identifier/Einmalige Produktkennung/Identificador único de dispositivo/Identifiant unique du dispositif/Identificativo univoco del dispositivo/Identificação única do dispositivo/Benzersiz cihaz tanımlayıcı
	Federal Law (USA) restricts this device to sale by or on the order of a physician./In den USA darf dieses Produkt laut Bundesgesetz nur durch einen Arzt oder auf ärztliche Verordnung verkauft werden./La legislación federal (EE. UU.) restringe la venta de este dispositivo a médicos o por prescripción facultativa./La loi fédérale (Etats-Unis) limite la vente de ce dispositif par un médecin ou sur ordonnance médicale./Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo a un medico o su prescrizione di un medico./A lei federal (EUA) restringe a venda deste dispositivo a médicos ou mediante prescrição médica./Federal Yasalar (ABD), bu cihazın satışını hekim tarafından veya talimatıyla yapılacak şekilde sınırlamaktadır.

Table of Contents

1	SCOPE	3
2	IMPORTANT NOTES	3
3	DESCRIPTION	3
4	WARNINGS	3
5	PRECAUTIONS	3
6	DEVICE APPLICATION	3
7	Performance Characteristics	3
8	COMPATIBLE DEVICES/ACCESSORIES	3
9	INSTRUCTIONS FOR USE	4
10	PACKAGING, STERILIZATION, AND REUSE	4
11	STORAGE CONDITIONS	5
12	DISPOSAL	5
13	Incident Reporting	5
14	ADMINISTRATIVE	28

1 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Signum Cable

2 IMPORTANT NOTES

	Before using these products, carefully review these Instructions for Use. Keep it readily available for reference.
	Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.

3 DESCRIPTION


The Signum Cable is designed to provide an easy method for making a secure connection between an Ad-Tech® electrode and standard clinic EEG monitoring and/or stimulation equipment.

The Signum Cable is available non-sterile and must be cleaned and sterilized prior to use. End-user sterilized by steam STERRAD® (Vaporized Hydrogen Peroxide Gas Plasma Sterilization).

The Signum Cables are reusable up to five times. Refer to cleaning and sterilization instructions in this document for proper cleaning and sterilization methods.

See Figure 1 for details of the available Signum Cable.

Figure 1 – Signum Cable

4 WARNINGS

	The Signum Cable is provided NON STERILE and must be cleaned and sterilized prior to use. Refer to the cleaning and sterilization instructions in this document.
	Disconnect from monitoring equipment during cardiac defibrillation.

5 PRECAUTIONS

	Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.
	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.

	A direct pull or stress on the electrode may cause damage to the electrode or a loss of contact recordings. If the patient will become mobile, place and secure a strain relief using the EEG cable.
	The Signum Cable may be used up to five (5) times. Refer to the cleaning and sterilization instructions in this document for proper cleaning and sterilization methods.
	Do not immerse the Signum Cable in fluids.
	Federal Law (USA) restricts this device to sale by or on the order of a physician.

6 DEVICE APPLICATION

INTENDED USE/PURPOSE

The Signum Cable is designed to provide a safe, secure, and easy method to connect Ad-Tech® electrodes with typical EEG monitoring and/or stimulation equipment.

INDICATIONS FOR USE

The Signum Cable is designed to provide a safe, secure, and easy method to connect Ad-Tech® electrodes with typical EEG monitoring and/or stimulation equipment.

CONTRAINDICATIONS

None known.

PATIENT TARGET GROUP

This device is used in patients diagnosed with intractable epilepsy who are undergoing brain mapping procedures to define the location of epileptogenic foci.

INTENDED USER

This product should only be used by a healthcare professional trained in the use of the Signum connection system.

ENVIRONMENT AND FACILITIES

Signum Cables are used in facilities which can accommodate EEG monitoring / stimulation equipment

EXPECTED CLINICAL BENEFIT

The Signum Cable provides a secure connection between the Ad-Tech® electrode and the EEG recording equipment.

UNDESIRABLE SIDE EFFECTS

No clinical risks with the use of the Signum Cable.

ELECTRONIC INSTRUCTIONS FOR USE

An electronic copy of these instructions for use is available from:

- Ad-Tech's website www.adtechmedical.com
- Email request to sales@adtechmedical.com, quoting the IFU document number AT10083-B

7 PERFORMANCE CHARACTERISTICS

DEVICE ATTRIBUTES

The Signum Cable consists of a connection cable with an attached electrode connector block on one end and the other end terminating in a standard 1.5 mm safety socket DIN connector.

DESIGN VARIANTS

The most common variants involve the number of contacts. Ad Tech® electrodes can have up to a maximum of sixteen (16) contacts per tail and the cables are often built to match. Variations to the cable style, cable length, and connector style are available to accommodate the procedural technique the physician defines as appropriate for the specific patient. These variants do not alter the indications or the intended use of the Signum Cable.

Catalog Number	Description	Cable Configuration (y)
L-JBL-y/DIN-aa*	Light Weight Signum Cables	2-64

*Note: "aa" denotes the standard configurations available from AD-TECH®. Contact an AD-TECH® sales representative for support.

8 COMPATIBLE DEVICES/ACCESSORIES

AD-TECH ELECTRODES

Electrode Family	Style
Depth Electrodes	Spencer Probe
	Macro-Micro
	Behnke Fried (Macro) Wire Bundle (Micro)

Subdural Electrodes	Grid Subdural
	Strip
	Dual-Sided Interhemispheric
	Multi-Strip and Split Grid

9 INSTRUCTIONS FOR USE

1. Inspect the Signum Cable for wear or damage. A cable that is or is suspected to be worn or damaged should be discarded.
2. Ensure that the connection block is in the open position. See **Figure 1**.


Fig. 1: Block in the Open Position

3. Immediately upon completion of electrode placement, the tail end contacts (1) need to be attached to the cable via the connector block. Carefully insert the tail of the electrode (1) into the connector block and ensure that the contacts on the electrode tail align with the pins in the connection block. See **Figure 2** for correct alignment.


Figure 2. Electrode Tail Fully Inserted into Block with Contacts Aligned

- 3.1. Ensure the electrode tail (1) is inserted fully by checking that the last contact inserted is **past** the protrusion (2) on the mating surface (See figure 2)
- 3.2. For electrodes with fewer contacts, insert the electrode until the last contact being inserted aligns with the far side of the first opening, and is past the protrusion on the mating surface. See Figure 3.


Figure 3. Correct positioning of electrode tail with fewer contacts


Figure 4. Incorrect positioning of electrode with fewer contacts.

- 3.3. Do not fully insert electrode with fewer contacts (shown in figure 4). The windows on the far side (away from where the electrode is inserted) should be left empty for electrodes with fewer contacts, as shown in figure 3.
- 3.4. If an electrode with fewer contacts is used incorrectly, contacting numbering of the DIN ends will not align with contact numbering of the electrodes.
4. Close the connection block using the flange handles to secure the electrode into the block. See **Figure 5**. At no time should the tail end contacts be exposed or unprotected.


Figure 5. Block Closed to Complete Connection

5. It is important to place and secure a strain relief using the EEG cable to prevent placing stress or a direct pull on the electrode which may cause loss of contact recordings and/or damage to the electrode.
6. Connect DIN end terminations to monitoring equipment.
7. If the patient will become mobile, the electrode must remain connected to the cable.

10 PACKAGING, STERILIZATION, AND REUSE

	Clean and sterilize the Signum Cable as directed. Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.
--	---

CLEANING AND STERILIZATION

	Non-sterile. Sterile before use.
--	----------------------------------

To help ensure patient safety, prevent cross contamination and provide maximum life from this product, proper reprocessing must be completed before each use.

Adhering to these validated reprocessing instructions, the cable may be used up to five (5) times.

It is the responsibility of the end user to ensure that the reprocessing is actually performed using appropriate equipment, materials, and personnel to achieve the desired result. Any deviation from these validated instructions should be evaluated for effectiveness and potential adverse consequences.

Manual Cleaning

The cable should be completely cleaned and free of visible soil, stains, and debris. If visible contaminants are present after cleaning, repeat the cleaning procedure. The cable can be cleaned using either Enzyme or Alkaline detergents.

Cleaning with Detergent

Enzymatic Detergent: Enzo® Enzymatic Detergent (Advanced Sterilization Products)

Alkaline Detergent: Prolystica® 2X Concentrate Alkaline Detergent (Steris Healthcare)

NOTE: At no point should the cable be immersed.

1. Soak a clean, soft cloth in lukewarm tap water. Use the cloth to wipe the cable and remove any heavy soil that may be present.
2. Prepare cleaning solution per manufacturer's recommendations.
3. Soak a clean, soft cloth in the prepared detergent. Use the cloth to wipe the cable and remove any visible soil, stains, and debris.
4. Use a soft bristled brush to remove soil from crevices. Dip brush in the prepared detergent prior to scrubbing the cable. Scrub all surfaces and hinged areas, reapplying detergent to the brush as necessary. Scrub for at least 2 minutes.
5. Re-wet the cloth and re-wipe the cable as needed to keep wet to mimic the manufacturer's suggested soak time.
6. Soak a clean, soft cloth in lukewarm tap water. Use the cloth to wipe the cable and remove the residual detergent. Re-wet the cloth and re-wipe the cable as needed.
7. Dry the cable with a clean, lint free cloth.
8. Visually inspect each cable for visible soil. If soil is still present on cable, repeat the cleaning procedure.

Sterilization

Ensure that the cable block is in the open position (See Figure 1), and not in the closed position (See Figure 6) when performing sterilization.

The cable is compatible with the following STERRAD® Systems (Advanced Sterilization Products):

- STERRAD® 100S
- STERRAD® NX™
- STERRAD® 100NX™

Follow the directions provided in the STERRAD® System User's Guides for the sterilization procedure.

STERRAD® Cycle Settings

Instrument	Cycle
STERRAD® 100S	Short
STERRAD® NX™	Standard
STERRAD® 100NX™	Standard

The cables are also compatible with the following STERIS® V-PRO Systems:

STERIS V-PRO® validated low temperature sterilization systems / cycle

Sterilizer / Cycle	Standard / Lumen	Non-Lumen	Flexible	Fast Non-Lumen	Fast
V-PRO® maX	✓	✓	✓		
V-PRO® maX 2	✓	✓	✓	✓	
V-PRO® 60	✓	✓	✓		
V-PRO® s2	✓	✓	✓		✓

REUSE

	Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles
--	--

The Signum Cable can be cleaned, re-sterilized, and reused 5 times.

11 STORAGE CONDITIONS

	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
--	---

Store the Signum Cables in their outer packaging, observing the handling and storage conditions detailed below.

	Fragile, handle with care.
	Keep away from sunlight.

	Keep dry
	Temperature limit Upper = 30 °C, Lower = 15 °C
	Humidity limitation Upper = 60%, Lower = 20%

12 DISPOSAL

	This device does not require any special disposal procedures, however, if the device has been contaminated with potentially infectious substances of human origin (such as bodily fluids), it should be handled according to hospital protocol for infectious medical waste. Dispose of the Signum Cable following the healthcare facility's biohazardous and pathological waste procedures.
--	---

13 INCIDENT REPORTING

EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or another person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the following:

- Death of a patient, user, or other person,
- Serious deterioration in the state of health of a patient, user, or other person,
- Potential for death or serious deterioration in the health of a patient, user, or other person if the incident were to recur.

UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening,
- Results in permanent impairment of a body function or permanent damage to a body structure,
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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