**INSTRUCTIONS FOR USE**  
**Products to work with BrainLab**  
**Navigation System - NaviVision™**

**SYMBOL GLOSSARY**

Symbol Símbolo Symbole Símbolo Símbolo Sembol	Description/ Beschreibung/ Descripción/ Description/ Descrizione/ Descrição/ Açıklama
	Manufacturer/ Hersteller/ Fabricante/ Fabricant/ Fabbicante/ Fabricante/ Üretici
	Authorized representative in the European Community/European Union/ Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/ Representante autorizado en la Comunidad Europea/Unión Europea/ Mandataire dans la Communauté européenne/Union européenne/ Mandatario nella Comunità Europea/Unione Europea/ Mandatário na Comunidade Europeia/União Europeia/ Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi
	Use-by date/ Haltbarkeitsdatum/ Fecha de caducidad/ Date limite d'utilisation/ Data di scadenza/ Data de validade/ Son kullanma tarihi
	Batch code/ Chargenbezeichnung/ Código de lote/ Code de lot/ Codice lotto/ Código de lote/ Parti kodu
	Catalog number/ Bestellnummer/ Número de catálogo/ Numéro de référence/ Numero di catalogo/ Número de catálogo/ Katalog numarası
	Non-sterile. Sterilize before use./ Unsteril. Vor Gebrauch sterilisieren./ No estéril. Esterilizar antes del uso./ Non stérile. Stériliser avant utilisation./ Non sterile. Sterilizzare prima dell'uso./ Não estéril. Esterilizar antes de utilizar./ Steril değildir. Kullanmadan önce sterilize edin.
	Sterilized using ethylene oxide/ Sterilisiert mit Ethylenoxid/ Esterilizado con óxido de etileno/ Stérilisé à l'oxyde d'éthylène/ Sterilizzato con ossido di etilene/ Esterilizado por óxido de etileno/ Etilen oksit kullanılarak sterilize edilmiştir
	Do not use if the package is damaged and consult instructions for use/ Nicht verwenden, wenn Verpackung beschädigt ist und Gebrauchsanweisung beachten/ No utilizar si el envase está dañado y consultar instrucciones de uso/ Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation/ Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso/ Não utilizar se existirem danos na embalagem e consultar as instruções de utilização/ Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
	Do not re-sterilize/ Nicht resterilisieren/ No reesterilizar/ Ne pas restériliser/ Non risterilizzare/ Não reesterilizar/ Tekrar sterilize etmeyin
	Do not reuse/ Nicht wiederverwenden/ No reutilizar/ Ne pas réutiliser/ Non riutilizzare/ Não reutilizar/ Tekrar kullanmayın
	Double sterile barrier system/ Doppel Sterilbarriere-System/ Sistema de barrera estéril doble/ Système de double barrière stérile/ Sistema a doppia barriera sterile/ Sistema de barreira estéril dupla/ Çift steril bariyer sistemi
	Fragile, handle with care/ Zerbrechlich, sorgfältig handhaben/ Frágil, manipular con cuidado/ Fragile, manipuler avec soin/ Fragile, maneggiare con cura/ Frágil, manusear com cuidado/ Kırılgan, dikkatli taşıyın
	Keep dry / Trocken aufbewahren/ Mantener seco/ Tenir au sec/ Mantenerse asciutto/ Manter seco/ Kuru tutun


	Temperature limit/ Temperaturgrenze/ Limite de temperatura/ Limites de température/ Limite di temperatura/ Limite de temperatura/ Sıcaklık sınırı  Upper/ obere / Superior/ Supérieure / Superiore/ Máximo/ Üst = 30°C,  Lower/ untere/ Inferior/ inférieure/ inferiore/ Mínimo/ Alt = 15°C
	Humidity limitation/ Feuchtigkeitsbeschränkung/ Limite de humedad/ Limites d'humidité/ Limitazione di umidità/ Limite de humidade/ Nem sınırları  Upper/ obere/ Superior/ Supérieure / Superiore/ Máximo/ Üst = 60%,  Lower/ untere / Inferior/ inférieure/ inferior/ Mínimo/ Alt = 20%
	Consult instructions for use / Gebrauchsanweisung beachten/ Consultar las instrucciones de uso/ Consulter la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consultar as instruções de utilização/ Kullanma talimatına başvurun
	Caution/ Vorsicht/ Atención/ Mise en garde/ Attenzione/ Atenção/ Dikkat
	Medical device/ Medizinprodukt/ Producto sanitario/ Dispositif médical/ Dispositivo medico/ Dispositivo médico/ Tıbbi cihaz
	Unique device identifier/ Einmalige Produktkennung (UDI)/ Identificador único de dispositivo/ Identifiant unique du dispositif/ Identificativo univoco del dispositivo/ Identificação única do dispositivo/ Bzersiz cihaz tanımlayıcı
	A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment./ Ein Medizinprodukt, das inakzeptable Risiken für den Patienten, medizinisches Personal oder andere Personen in der MR- Umgebung darstellt./ Un producto sanitario que plantea riesgos inaceptables para el paciente, el personal médico u otras personas dentro del entorno de RM./ Un dispositif médical présentant des risques inacceptables pour le patient, le personnel médical, ou toute autre personne à l'intérieur de l'environnement de RM./ Un dispositivo medico che comporta rischi inaccettabili per il paziente, gli operatori o qualsiasi altro individuo in ambiente RM./ Um dispositivo médico que represente riscos inaceitáveis para o doente, equipa médica ou outras pessoas dentro do ambiente de RM./ MR ortamında hasta, sağlık çalışanları veya diğer kişiler için kabul edilemeyecek riskler taşıyan bir tıbbi cihaz.
	Federal Law (USA) restricts this device to sale by or on the order of a physician. / In den USA darf dieses Produkt laut Bundesgesetz nur durch einen Arzt oder auf ärztliche Verordnung verkauft werden./ La legislación federal (EE. UU.) restringe la venta de este dispositivo a médicos o por prescripción facultativa./ La loi fédérale (États-Unis) limite la vente de ce dispositif par un médecin ou sur ordonnance médicale./ Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo a un medico o su prescrizione di un medico./ A lei federal (EUA) restringe a venda deste dispositivo a médicos ou mediante prescrição médica./ Federal Yasalar (ABD), bu cihazın satışını hekim tarafından veya talimatıyla yapılacak şekilde sınırlamaktadır.

Table of Contents **en**

**1 SCOPE** ..... 3

**2 IMPORTANT NOTES** ..... 3

**3 DESCRIPTION** ..... 3

**4 WARNINGS** ..... 3

**5 PRECAUTIONS** ..... 3

**6 DEVICE APPLICATION** ..... 3

6.1 INTENDED USE/PURPOSE ..... 3

6.2 INDICATIONS FOR USE ..... 3

6.3 CONTRAINDICATIONS ..... 3

6.4 PATIENT TARGET GROUP ..... 3

6.5 INTENDED USER ..... 3

6.6 ENVIRONMENT AND FACILITIES ..... 4

6.7 EXPECTED CLINICAL BENEFIT ..... 4

6.8 UNDESIRABLE SIDE EFFECTS ..... 4

**7 PERFORMANCE CHARACTERISTICS** ..... 4

7.1 DEVICE ATTRIBUTES ..... 4

**8 INSTRUCTIONS FOR USE** ..... 4

**9 MRI SAFETY INFORMATION** ..... 4

**10 PACKAGING, STERILIZATION, AND REUSE** ..... 5

10.1 STERILE BrainLab Navigation Guide ..... 5

10.2 REUSE ..... 5

**11 STORAGE CONDITIONS** ..... 5

**12 DISPOSAL** ..... 5


**13 INCIDENT REPORTING** ..... 5

13.1 EUROPEAN UNION ..... 5

13.2 CANADA ..... 5

13.3 UNITED STATES ..... 5

See Figure 1 for details of the available BrainLab Navigation Guide.

**4 WARNINGS**

No known warnings for the BrainLab Navigation Guide (BLNG).

**5 PRECAUTIONS**

	Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.
	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
	Federal Law (USA) restricts this device to sale by or on the order of a physician.

**6 DEVICE APPLICATION**

**6.1 INTENDED USE/PURPOSE**

The AD-TECH® BrainLab Navigation Guide (BLNG) is a depth placement device used as an accessory item for placing or using Ad-Tech's electrodes, electrode systems, or Biopsy Needles.

**6.2 INDICATIONS FOR USE**

There are no specific indications; this is an accessory to the Depth Electrodes

**6.3 CONTRAINDICATIONS**

None known.

**6.4 PATIENT TARGET GROUP**

Applicable to the Depth Electrodes and the Placement Accessories.

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

Applicable to the Brain Biopsy Needle and the Placement Accessories:

The Brain Biopsy Needles are used in pediatric and adult populations, regardless of gender.

**6.5 INTENDED USER**

The BrainLab Navigation Guide (BLNG) is to be used by physicians/surgeons experienced in using Ad-Tech's electrodes, electrode systems, or Biopsy Needles.

**1 SCOPE**

These instructions for use apply to the following device:

- BLNG-2X BrainLab Navigation Guide

**2 IMPORTANT NOTES**

	Before using these products, carefully review these Instructions for Use. Keep it readily available for reference.
	Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.

**3 DESCRIPTION**

The AD-TECH® BrainLab Navigation Guide (BLNG) is a depth placement device used as an accessory item for placing or using Ad-Tech's electrodes, electrode systems, or Biopsy Needles.

The BrainLab Navigation Guide is available sterile. End-user sterilized by steam (Auto-clave), STERRAD® (Vaporized Hydrogen Peroxide Gas Plasma Sterilization), or Ethylene Oxide (EO).

The BrainLab Navigation Guides are reusable up to twenty-five times. *Note: As the device was provided sterile by the manufacturer, the product has already undergone 1 cycle.*

## 6.6 ENVIRONMENT AND FACILITIES

BrainLab Navigation Guides (BLNG) are for surgical use only.  
Placement and removal must be in a sterile surgical environment.

## 6.7 EXPECTED CLINICAL BENEFIT

The AD-TECH® BrainLab Navigation Guide (BLNG) is a depth placement device used as an accessory item for placing or using Ad-Tech's electrodes, electrode systems, or Biopsy Needles.

## 6.8 UNDESIRABLE SIDE EFFECTS

No clinical risks with the use of the BrainLab Navigation Guide (BLNG).

## 7 PERFORMANCE CHARACTERISTICS

### 7.1 DEVICE ATTRIBUTES

The AD-TECH® BrainLab Navigation Guide (BLNG) is a depth placement device used as an accessory item for placing or using Ad-Tech's electrodes, electrode systems, or Biopsy Needles.

Catalog #	Description
BLNG-2X	BrainLab Navigation Guide

## 8 INSTRUCTIONS FOR USE

Use of the BLNG device with interfacing Two-piece Slotted Cannula is illustrated on the following page (Figure 2).


Descriptions of each step are defined below:

- 1a. Insert BLNG into BrainLab navigation system, tightening the lower hub as pictured.
- 1b. Insert Two-piece slotted cannula into the inner lumen of the BLNG device (ensure set screw is not blocking insertion path of interfacing device).
- 2a. Tighten set screw to secure Slotted Cannula device (do not overtighten).
- 3a. Remove inner obturator from slotted cannula device.
- 4a. Insert depth electrode.
- 4b. Remove electrode stylet.
- 4c. Loosen set screw securing slotted cannula device.
- 4d. Remove outer sheath of slotted cannula device.
- 4e. Loosen lower BLNG hub.
- 4f. Remove BLNG device, ensuring electrode position is maintained during removal.

Use of the BLNG device with interfacing Biopsy Needle is defined below (relevant illustrations from Figure 2 are referenced).

- 1a. Insert BLNG into BrainLab navigation system, tightening the lower hub as pictured.
- 1b. Insert Biopsy Needle into the inner lumen of the BLNG device (ensure set screw is not blocking insertion path of interfacing device).
2. Tighten set screw to secure Biopsy Needle device (do not overtighten) (shown in 2a).
3. Use Biopsy Needle Device in accordance with the associated Biopsy Needle Directions for Use.
4. Loosen set screw securing Biopsy Needle device (4c).
5. Remove Biopsy Needle from BLNG device.
6. Loosen lower BLNG hub (4e).
7. Remove BLNG device (4f).

Figure 2:

## 9 MRI SAFETY INFORMATION

The BrainLab Navigation Guides are MRI-unsafe.

Using the BrainLab Navigation Guides in an MRI environment may severely harm the user and patient.

AD-TECH® has not verified the safety and performance of the BrainLab Navigation Guides in the MRI environment.

AD-TECH® has classified the BrainLab Navigation Guides as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices And Other Items For Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

AD-TECH® has labeled the BrainLab Navigation Guides as MRI Unsafe using the symbol:

**10 PACKAGING, STERILIZATION, AND REUSE**

	<p><i>Clean and sterilize the BrainLab Navigation Guide as directed.</i></p> <p><i>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</i></p>
--	--

**10.1 STERILE BrainLab Navigation Guide**

	<p>Sterile.</p>
--	-----------------

After use, clean and sterilize the BrainLab Navigation Guide as per AT10054-B AD-TECH® CLEANING AND STERILIZATION INSTRUCTIONS PLACEMENT TOOLS AND ACCESSORIES.

**10.2 REUSE**

	<p><i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i></p>
--	---

The BrainLab Navigation Guide can be cleaned and re-sterilized for use, up to a total of twenty-five times.

*NOTE: If the device was provided sterile by the manufacturer, the product has already undergone 1 cycle.*

**11 STORAGE CONDITIONS**

	<p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p>
--	---

Store the BrainLab Navigation Guide in their outer packaging, observing the handling and storage conditions detailed below.

	<p>Fragile, handle with care.</p>
	<p>Keep dry</p>
	<p>Temperature limit Upper = 30°C, Lower = 15°C</p>
	<p>Humidity limitation Upper = 60%, Lower = 20%</p>

**12 DISPOSAL**

	<p><i>Dispose of the BrainLab Navigation Guide following the healthcare facility's biohazardous and pathological waste procedures.</i></p>
--	--

**13 INCIDENT REPORTING**

**13.1 EUROPEAN UNION**

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

**13.2 CANADA**

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

**13.3 UNITED STATES**

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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