INSTRUCTIONS FOR USE

Disposable Drill Stop and Wrench

SYMBOL GLOSSARY

Symbol Símbolo Symbole Símbolo Símbolo Símbolo	Description/ Beschreibung/ Descripción/ Description/ Descrizione/ Descrição/ Açıklama
	Manufacturer/ Hersteller/ Fabricante/ Fabricant/ Fabbricante/ Fabricante/ Üretici
	Authorized representative in the European Community/European Union/ Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union/ Representante autorizado en la Comunidad Europea/Unión Europea/ Mandataire dans la Communauté européenne/Union européenne/ Mandatario nella Comunità Europea/Unione Europea/ Mandatário na Comunidade Europeia/União Europeia/ Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi
	Use-by date/ Haltbarkeitsdatum/ Fecha de caducidad/ Date limite d'utilisation/ Data di scadenza/ Data de validade/ Son kullanma tarihi
	Batch code/ Chargenbezeichnung/ Código de lote/ Code de lot/ Codice lotto/ Código de lote/ Parti kodu
	Catalog number/ Bestellnummer/ Número de catálogo/ Numéro de référence/ Numero di catalogo/ Número de catálogo/ Katalog numarası
	Non-sterile. Sterilize before use./ Unsteril. Vor Gebrauch sterilisieren./ No estéril. Esterilizar antes del uso./ Non stérile. Stériliser avant utilisation./ Non sterile. Sterilizzare prima dell'uso./ Não estéril. Esterilizar antes de utilizar./ Steril değildir. Kullanmadan önce sterilize edin.
	Sterilized using ethylene oxide/ Sterilisiert mit Ethylenoxid/ Esterilizado con óxido de etileno/ Stérilisé à l'oxyde d'éthylène/ Sterilizzato con ossido di etilene/ Esterilizado por óxido de etileno/ Etilen oksit kullanılarak sterilize edilmiştir
	Do not use if the package is damaged and consult instructions for use/ Nicht verwenden, wenn Verpackung beschädigt ist und Gebrauchsanweisung beachten/ No utilizar si el envase está dañado y consultar instrucciones de uso/ Ne pas utiliser si le conditionnement est endommagé et consulter la notice d'utilisation/ Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso/ Não utilizar se existirem danos na embalagem e consultar as instruções de utilização/ Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
	Do not re-sterilize/ Nicht reesterilisieren/ No reesterilizar/ Ne pas resteriliser/ Non risterilizzare/ Não reesterilizar/ Tekrar sterilize etmeyin
	Do not reuse/ Nicht wiederverwenden/ No reutilizar/ Ne pas réutiliser/ Non riutilizzare/ Não reutilizar/ Tekrar kullanmayın
	Double sterile barrier system/ Doppel Sterilbarriere-System/ Sistema de barrera estéril doble/ Système de double barrière stérile/ Sistema a doppia barriera sterile/ Sistema de barreira estéril dupla/ Çift steril bariyer sistemi
	Fragile, handle with care/ Zerbrechlich, sorgfältig handhaben/ Frágil, manipular con cuidado/ Fragile, manipuler avec soin/ Fragile, maneggiare con cura/ Frágil, manusear com cuidado/ Kırılgan, dikkatli taşıyın
	Keep dry / Trocken aufbewahren/ Mantener seco/ Tenir au sec/ Mantenerse asciutto/ Manter seco/ Kuru tutun

	Temperature limit/ Temperaturgrenze/ Limite de temperatura/ Limites de température/ Limite di temperatura/ Limite de temperatura/ Sıcaklık sınırı Upper/ obere / Superior/ Supérieure / Superiore/ Máximo/ Üst = 30°C, Lower /untere/ Inferior/ inférieure/ inferiore/ Mínimo/ Alt = 15°C
	Humidity limitation/ Feuchtigkeitsbeschränkung/ Limite de humedad/ Limites d'humidité/ Limitazione di umidità/ Limite de humidade/ Nem sınırları Upper/ obere/ Superior/ Supérieure / Superiore/ Máximo/ Üst = 60%, Lower/ untere / Inferior/ inférieure/ inferior/ Mínimo/ Alt = 20%
	Consult instructions for use / Gebrauchsanweisung beachten/ Consultar las instrucciones de uso/ Consulter la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consultar as instruções de utilização/ Kullanma talimatına başvurun
	Caution/ Vorsicht/ Atención/ Mise en garde/ Attenzione/ Atenção/ Dikkat
	Medical device/ Medizinprodukt/ Producto sanitario/ Dispositif médical/ Dispositivo medico/ Dispositivo médico/ Tıbbi cihaz
	Unique device identifier/ Einmalige Produktkennung (UDI)/ Identificador único de dispositivo/ Identifiant unique du dispositif/ Identificativo univoco del dispositivo/ Identificação única do dispositivo/ Benzersiz cihaz tanımlayıcı
	A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment./ Ein Medizinprodukt, das inakzeptable Risiken für den Patienten, medizinisches Personal oder andere Personen in der MR- Umgebung darstellt./ Un producto sanitario que plantea riesgos inaceptables para el paciente, el personal médico u otras personas dentro del entorno de RM./ Un dispositif médical présentant des risques inacceptables pour le patient, le personnel médical, ou toute autre personne à l'intérieur de l'environnement de RM./ Un dispositivo medico che comporta rischi inaccettabili per il paziente, gli operatori o qualsiasi altro individuo in ambiente RM./ Um dispositivo médico que represente riscos inaceitáveis para o doente, equipa médica ou outras pessoas dentro do ambiente de RM./ MR ortamında hasta, sağlık çalışanları veya diğer kişiler için kabul edilemeyecek riskler taşıyan bir tıbbi cihaz.
	Federal Law (USA) restricts this device to sale by or on the order of a physician. / In den USA darf dieses Produkt laut Bundesgesetz nur durch einen Arzt oder auf ärztliche Verordnung verkauft werden./ La legislación federal (EE. UU.) restringe la venta de este dispositivo a médicos o por prescripción facultativa./ La loi fédérale (États-Unis) limite la vente de ce dispositif par un médecin ou sur ordonnance médicale./ Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo a un medico o su prescrizione di un medico./ A lei federal (EUA) restringe a venda deste dispositivo a médicos ou mediante prescrição médica./ Federal Yasalar (ABD), bu cihazın satışını hekim tarafından veya talimatıyla yapılacak şekilde sınırlamaktadır.

Table of Contents **en**

1 SCOPE 3

2 IMPORTANT NOTES 3

3 DESCRIPTION 3

4 WARNINGS 3

5 PRECAUTIONS 3

6 DEVICE APPLICATION 3

6.1 INTENDED USE/PURPOSE 3

6.2 INDICATIONS FOR USE 3

6.3 CONTRAINDICATIONS 3

6.4 PATIENT TARGET GROUP 4

6.5 INTENDED USER 4

6.6 ENVIRONMENT AND FACILITIES 4

6.7 EXPECTED CLINICAL BENEFIT 4

6.8 UNDESIRABLE SIDE EFFECTS 4

7 PERFORMANCE CHARACTERISTICS 4

7.1 DEVICE ATTRIBUTES 4

8 INSTRUCTIONS FOR USE 5

9 MRI SAFETY INFORMATION 5

9.1 MRI SAFETY, Drill Stop and Drill Stop Wrench 5

10 PACKAGING, STERILIZATION, AND REUSE 5

10.1 NON-STERILE Drill Stop and Drill Stop Wrench 6

10.2 STERILE Drill Stop and Drill Stop Wrench (within the Disposable Drill Kit) 6

10.3 REUSE 6

11 STORAGE CONDITIONS 6

12 DISPOSAL 6

13 INCIDENT REPORTING 6

13.1 EUROPEAN UNION 6

13.2 CANADA 6

13.3 UNITED STATES 6

1 ADMINISTRATIVE/ ADMINISTRATIV/ INFORMACIÓN ADMINISTRATIVA/ ADMINISTRATIV/ SEZIONE AMMINISTRATIVA/ DADOS ADMINISTRATIVOS/ IDARI 31

1 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Drill Stop
- AD-TECH® Drill Stop Wrench

2 IMPORTANT NOTES

	Before using these products, carefully review these Instructions for Use. Keep it readily available for reference.
	Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.

3 DESCRIPTION

The AD-TECH® Drill Stop and Drill Stop Wrench are accessories to the AD-TECH® Cranial Drill Bits. The Drill Stop is attached to a Cranial Drill Bit by the Drill Stop Wrench as a depth guide to the drill bit.

The Cranial Drill Bits and Accessories are intended to be used for drilling holes in the skull for neurological procedures, such as brain biopsy, brain contacting electrode and electrode accessory device placement.

These accessories are included with the Disposable Drill Kit or sold separately.

The Drill Stop and Drill Stop Wrench are available non-sterile (if sold separately) or sterile by EO (if included with the Disposable Drill Kit).

Sterile product by Ethylene Oxide (EO). Non-sterile product is end-user sterilized by steam (Auto-clave), STERRAD® (Vaporized Hydrogen Peroxide Gas Plasma Sterilization), or EO.

See Figure 1 for details of the available Drill Stop and Drill Stop Wrench.

Drill Stop and Drill Stop Wrench		
	DDS-2.4N	Contains: 1 Drill Stop 
	DDS-2.8N	Contains: 1 Drill Stop  1 Wrench 
	DDSW-N	Contains: 1 Wrench 

4 WARNINGS

No known warnings for the Drill Stop and Drill Stop Wrench.

5 PRECAUTIONS

	Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.
	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
	Federal Law (USA) restricts this device to sale by or on the order of a physician.

6 DEVICE APPLICATION

6.1 INTENDED USE/PURPOSE

The Drill Stop and Drill Stop Wrench are accessories to the Ad-Tech Cranial Drill Bits. The Drill Stop is attached to a Cranial Drill Bit by the Drill Stop Wrench as a depth guide to the drill bit.

The Cranial Drill Bits and accessories are intended to be used for drilling holes in the skull for neurological procedures, such as brain biopsy, brain contacting electrode, and electrode accessory device placement.

6.2 INDICATIONS FOR USE

There are no specific indications; this is an accessory to the Cranial Drill Bits.

6.3 CONTRAINDICATIONS

	Disposable Drill Stops and Wrenches are to be used with Ad-Tech's Disposable Drill Kits and should not be used on any patient whom the physician/surgeon considers at risk for infection or on whom the use cannot be performed safely.
---	---

6.4 PATIENT TARGET GROUP

Applicable to the Depth Electrodes and the Placement Accessories.

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

6.5 INTENDED USER

The Drill Stop and Drill Stop Wrench are to be used by physicians/surgeons experienced in the use of Disposable Cranial Drill Kits.

6.6 ENVIRONMENT AND FACILITIES

Drill Stop and Drill Stop Wrench are for surgical use only.

The use of these devices must be in a sterile surgical environment.

6.7 EXPECTED CLINICAL BENEFIT

The Drill Stop and Drill Stop Wrench are accessories to the Ad-Tech Cranial Drill Bits. The Drill Stop is attached to a Cranial Drill Bit by the Drill Stop Wrench as a depth guide to the drill bit.

6.8 UNDESIRABLE SIDE EFFECTS

No clinical risks with the use of the Drill Stop and Drill Stop Wrench.

7 PERFORMANCE CHARACTERISTICS


7.1 DEVICE ATTRIBUTES

The AD-TECH® Drill Stop and Drill Stop Wrench are accessories to the AD-TECH® Cranial Drill Bits. The Drill Stop is attached to a Cranial Drill Bit by the Drill Stop Wrench as a depth guide to the drill bit.

These accessories are included with the Disposable Drill Kit or sold separately.

Disposable Drill Stop and Drill Wrench provided with the Disposable Drill Kit:

Catalog Number	Description	Drill			Drill Stop Quantity	Drill Wrench Quantity	Sterile Yes/No
		Quantity	Length (YY)	Diameter (a.a)			
DDK2-2.4-30N	Disposable Drill Kit – Non-sterile	2	30cm	2.4mm	2	2	No
DRL-2.4-22N	Disposable Drill Bit – Non-sterile	1	22cm	2.4mm	0	0	No
DDK2-2.4-16X	Disposable Drill Kit – Sterile	2	16cm	2.4mm	2	2	Yes
DDK2-2.4-30X	Disposable Drill Kit – Sterile	2	30cm	2.4mm	2	2	Yes
DDK2-2.8-30X	Disposable Drill Kit – Sterile	2	30cm	2.8mm	2	2	Yes
DDK2-3.2-30X	Disposable Drill Kit – Sterile	2	30cm	3.2mm	2	2	Yes


Sold separately:

Catalog Number	Name	Inner Diameter	Supplied Sterile Yes/No
DDS-2.4N	Drill Stop	2.4mm	No
DDS-2.8N	Drill Stop	2.8mm	No
DDSW-N	Drill Stop Wrench	N/A	No

8 INSTRUCTIONS FOR USE

	<p>Do not over-tighten the Drill Stop set screw when applying the drill bit.</p> <p>Overtightening may deform the drill bit.</p>
	<p>The Drill Stop will not stop the drill. The Drill Stop is designed only to provide the user with a marker for drilling depth.</p>

Regarding Figure 2, the steps for deploying the Cranial Drill Kit are as follows:
Step 1a: Slide the Drill Stop onto the distal end of the Cranial Drill Bit and set the position of the Drill Stop to achieve the required drilling depth.
Step 1b, c: Tighten the Drill Stop set screw using the Drill Wrench in a clockwise direction taking care not to overtighten the screw. The Cranial Drill Bit and Drill Stop are now ready for use.

Refer to the Cranial Drill Kit IFU for use.

9 MRI SAFETY INFORMATION

	<p>The Drill Stop and Drill Stop Wrench are MRI unsafe.</p> <p>Using the Drill Stop and Drill Stop Wrench in an MRI environment may lead to severe patient harm.</p>
--	--

9.1 MRI SAFETY, Drill Stop and Drill Stop Wrench

AD-TECH® has not verified the safety and performance of the Drill Stop and Drill Stop Wrench in the MRI environment.

AD-TECH® has classified these accessories as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices and Other Items for Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

AD-TECH® has labeled the Drill Stop and Drill Stop Wrench as MRI Unsafe using the symbol:

10 PACKAGING, STERILIZATION, AND REUSE

	<p>Clean and sterilize the Drill Stop and Drill Stop Wrench as directed.</p> <p>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</p>
--	---

10.1 NON-STERILE Drill Stop and Drill Stop Wrench

	Non-sterile. Sterilize before use.
--	------------------------------------

Before use, clean and sterilize the Drill Stop and Drill Stop Wrench as per AT10054-B AD-TECH® CLEANING AND STERILIZATION INSTRUCTIONS PLACEMENT TOOLS AND ACCESSORIES.

10.2 STERILE Drill Stop and Drill Stop Wrench (within the Disposable Drill Kit)

	AD-TECH® sterilizes the Cranial Drill Bits and Accessories using ethylene oxide.
	Do not use the Cranial Drill Bits and Accessories if the packaging is damaged or unintentionally opened before use. Using a device from damaged packing may result in contamination and severe patient harm.
	Do not re-sterilize the sterile Cranial Drill Bits and Accessories. Re-sterilization may result in device failure.
	The Cranial Drill Bits and Accessories are packaged in a double sterile barrier to support surgical asepsis procedures.

10.3 REUSE

	The Cranial Drill Bits and Accessories are for single use only. Do not reuse the devices. Reuse of the devices may result in contamination and severe patient harm.
---	--

The Drill Stop and Drill Stop Wrench can be cleaned and re-sterilized for use, up to a total of twenty-five times.

NOTE: If the devices were provided sterile by the manufacturer, the product has already undergone 1 cycle.

11 STORAGE CONDITIONS

	Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.
--	--

Store the Drill Stop and Drill Stop Wrench in their outer packaging, observing the handling and storage conditions detailed below.

	Fragile, handle with care.
	Keep dry
	Temperature limit Upper = 30°C, Lower = 15°C
	Humidity limitation Upper = 60%, Lower = 20%

12 DISPOSAL

	Dispose of the Drill Stop and Drill Stop Wrench following the healthcare facility's biohazardous and pathological waste procedures.
---	---

13 INCIDENT REPORTING

13.1 EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

13.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

13.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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