INSTRUCTIONS FOR USE

Leksell Guide Blocks

SYMBOL GLOSSARY

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10 MRI SAFETY INFORMATION 4

10.1 MRI SAFETY, LEKSELL GUIDE BLOCKS 4

11 PACKAGING, STERILIZATION, AND REUSE 5

11.1 NON-STERILE GUIDE BLOCKS 5

11.2 REUSE 5

12 STORAGE CONDITIONS 5

13 DISPOSAL 5

14 INCIDENT REPORTING 5

14.1 EUROPEAN UNION 5

14.2 CANADA 5

14.3 UNITED STATES 5

1 SCOPE

These instructions for use apply to the following devices:

- AD-TECH® Leksell Guide Blocks

2 IMPORTANT NOTES

	<i>Before using these products, carefully review these Instructions for Use. Keep it readily available for reference.</i>
	<i>Carefully read all the symbol-marked cautions and warnings. Incorrect use of the products may lead to severe patient harm.</i>

3 DESCRIPTION


The AD-TECH® Leksell Guide Block is an accessory to the Cranial Drill Bits. The Guide Blocks support the placement of the Drill Sleeve Guide on a stereotactic frame.

The Cranial Drill Bits and Accessories are intended to be used for drilling holes in the skull for neurological procedures, such as brain biopsy, brain contacting electrode and electrode accessory device placement.

The Guide Block is available non-sterile. End-user sterilized by steam (Autoclave), STERRAD® (Vaporized Hydrogen Peroxide Gas Plasma Sterilization), or Ethylene Oxide (EO).

The Guide Blocks are reusable up to twenty-five times.

See Figure 1 for details of the available Leksell Guide Block.

4 WARNINGS

No known warnings for the Guide Block.

5 PRECAUTIONS

	<i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i>
	<i>Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.</i>
	<i>Federal Law (USA) restricts this device to sale by or on the order of a physician.</i>

6 DEVICE APPLICATION

6.1 INTENDED USE/PURPOSE

The Guide Blocks are accessories to the Ad-Tech Cranial Drill Bits. The Guide Blocks support placement of the Ad-Tech Drill Sleeve Guide on a stereotactic frame.

6.2 INDICATIONS FOR USE

There are no specific indications; this is an accessory to the Cranial Drill Bits.

6.3 CONTRAINDICATIONS

	<i>Ad-Tech's Guide Blocks are an optional accessory to be used to aid in the placement of Ad-Tech's Depth Electrodes and Ad-Tech Cranial Drill Bits. The Guide Blocks should not be used on any patient whom the physician/surgeon considers at risk for infection or on whom the use cannot be performed safely.</i>
--	---

6.4 PATIENT TARGET GROUP

Applicable to the Depth Electrodes and the Placement Accessories.

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

6.5 INTENDED USER

The Guide Blocks are to be used by physicians/surgeons experienced in the use of Disposable Cranial Drill Kits.

6.6 ENVIRONMENT AND FACILITIES

Guide Blocks are for surgical use only.

Placement and removal must be in a sterile surgical environment.

6.7 EXPECTED CLINICAL BENEFIT

The Guide Blocks are accessories to the Ad-Tech Cranial Drill Bits. The Guide Blocks support placement of the Ad-Tech Drill Sleeve Guide on a stereotactic frame.

6.8 UNDESIRABLE SIDE EFFECTS

No clinical risks with the use of the Guide Block.

7 PERFORMANCE CHARACTERISTICS

7.1 DEVICE ATTRIBUTES

The AD-TECH® Leksell Guide Block is an accessory to the Cranial Drill Bits. The Guide Blocks support the placement of the Drill Sleeve Guide on a stereotactic frame.

Catalog #	Name	Description	Supplied Sterile Yes/No
LK-GB-N-13	Guide Blocks	Leksell Coordinate G-Frame, Model 1002248 Stereotactic System®.	No

8 COMPATIBLE DEVICES/ACCESSORIES

8.1 DRILL SLEEVE GUIDE

The Drill Sleeve Guide is an accessory to the Ad-Tech Cranial Drill Bits. It is a long, stainless steel tube with a through-hole inner diameter sized to interface with the Cranial Drill Bit, which maintains a specified trajectory along the skull when used with stereotactic equipment. The Guide Blocks support the placement of the Drill Sleeve Guide on a stereotactic frame.

The following Drill Sleeve Guide is available:

Catalog Number	Name	Length (YYY)	Inner Diameter (ID)	Outer Diameter (OD)	Supplied Sterile Yes/No
DSG-6.3-90-2.4N	Drill Sleeve Guide	90 mm	2.4 mm	6.3 mm	No

9 INSTRUCTIONS FOR USE

	<i>If the Guide Blocks do not securely fit within the noted stereotactic frame, do not force connection.</i>
	<i>Do not apply to a stereotactic frame that is damaged.</i>

The Guide Blocks support the placement of the Drill Sleeve Guide on the indicated Stereotactic Frame.

Catalog Number	Frame
LK-GB-N-13	Leksell Coordinate G-Frame, Model 1002248 Stereotactic System®.

The steps for placing the Leksell Guide Blocks are as follows:

- Step 1:** Once the applicable stereotactic frame listed above is in place for the trajectory, securely fit the Guide Blocks in the frame.
- Step 2:** Insert the Drill Sleeve Guide in the Guide Blocks.
- Step 3:** Place the Drill Bit through the Drill Sleeve Guide and drill hole through the skull for the Anchor Bolt.
- Step 4:** Once the hole is made for the Anchor Bolt, remove the Drill Bit and Drill Sleeve Guide.
- Step 5:** Remove Guide Blocks.

10 MRI SAFETY INFORMATION

	<i>The Leksell Guide Blocks are MRI unsafe. Using the Guide Blocks in an MRI environment may lead to severe patient harm.</i>
---	---

10.1 MRI SAFETY, LEKSELL GUIDE BLOCKS

AD-TECH® has not verified the safety and performance of the Leksell Guide Blocks in the MRI environment.

AD-TECH® has classified these Guide Blocks as MRI Unsafe and labeled them appropriately.

ASTM F2503-20 Standard Practice For Marking Medical Devices and Other Items for Safety In The Magnetic Resonance Environment defines MRI Unsafe as:

A medical device that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

AD-TECH® has labeled the Leksell Guide Blocks as MRI Unsafe using the symbol:

11 PACKAGING, STERILIZATION, AND REUSE

	<p><i>Clean and sterilize the Guide Blocks as directed.</i></p> <p><i>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</i></p>
--	---

11.1 NON-STERILE GUIDE BLOCKS

	<p>Non-sterile. Sterilize before use.</p>
--	---

Before and after use, clean and sterilize the Guide Blocks as per AT10054-B AD-TECH® CLEANING AND STERILIZATION INSTRUCTIONS PLACEMENT TOOLS AND ACCESSORIES.

11.2 REUSE

	<p><i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i></p>
--	---

The Guide Blocks can be cleaned and re-sterilized for use, up to a total of twenty-five times.

12 STORAGE CONDITIONS

	<p><i>Observe the recommended storage conditions.</i></p> <p><i>Failure to observe storage conditions may result in device failure.</i></p>
--	---

Store the Guide Blocks in their outer packaging, observing the handling and storage conditions detailed below.

	<p>Fragile, handle with care.</p>
	<p>Keep dry</p>
	<p>Temperature limit Upper = 30°C, Lower = 15°C</p>
	<p>Humidity limitation Upper = 60%, Lower = 20%</p>

13 DISPOSAL

	<p><i>Dispose of the Guide Blocks following the healthcare facility's biohazardous and pathological waste procedures.</i></p>
--	---

14 INCIDENT REPORTING

14.1 EUROPEAN UNION

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

14.2 CANADA

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

14.3 UNITED STATES

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or
- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

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