# **INSTRUCTIONS FOR USE**

## **LSB Anchor Bolt Placement Wrenches**

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
7.1 DEVICE ATTRIBUTES ..... 4

8 INSTRUCTIONS FOR USE ..... 4

9 PACKAGING, STERILIZATION, AND REUSE ..... 4

10.1 STERILE/NON-STERILE Anchor Bolt Placement Wrenches ..... 4

10.2 REUSE ..... 4

10 STORAGE CONDITIONS ..... 5

11 DISPOSAL ..... 5

12 INCIDENT REPORTING ..... 5

13.1 EUROPEAN UNION ..... 5

13.2 CANADA ..... 5

13.3 UNITED STATES ..... 5

1 ADMINISTRATIVE/ ADMINISTRATIV/ INFORMACIÓN ADMINISTRATIVA/ ADMINISTRATIF/ SEZIONE AMMINISTRATIVA/ DADOS ADMINISTRATIVOS/ IDARI ..... 24

**1 SCOPE**

These instructions for use apply to the following devices:

- Ad-Tech® LSB Anchor Bolt Placement Wrenches – Non-Sterile
- Ad-Tech® LSB Anchor Bolt Placement Wrenches - Sterile

**2 IMPORTANT NOTES**

	<p>Before using these products, carefully review these Instructions for Use.</p> <p>Keep it readily available for reference.</p>
	<p>Carefully read all the symbol-marked cautions and warnings.</p> <p>Incorrect use of the products may lead to severe patient harm.</p>

**3 DESCRIPTION**


The AD-TECH® LSB Anchor Bolt Wrench is used to screw the Anchor Bolt into the patient’s skull and remove it at the end of the monitoring period.

The Anchor Bolt Placement wrenches are available either sterile or non-sterile.

The Anchor Bolt Placement wrenches are reusable up to twenty-five times.

See Figure 1 for details of the available Placement wrenches.

Figure 1:

**4 WARNINGS**

No known warnings for the Anchor Bolt Placement Wrenches.

**5 PRECAUTIONS**

	<p>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</p>
	<p>Observe the recommended storage conditions.</p> <p>Failure to observe storage conditions may result in device failure.</p>
	<p>Federal Law (USA) restricts this device to sale by or on the order of a physician.</p>
	<p>Discontinue Anchor Bolt placement if you experience any instrument (e.g., the Placement/Removal Wrench, Cranial Drill, etc.) bending during the placement procedure.</p>

**6 DEVICE APPLICATION**

**6.1 INTENDED USE/PURPOSE**

The Anchor Bolt Placement Wrench is an accessory to the Ad-Tech Anchor Bolts. The Wrench is used to screw the Anchor Bolt into the patient’s skull and remove it at the end of the monitoring period.

**6.2 INDICATIONS FOR USE**

There are no specific indications; this is an accessory to the Anchor Bolts.

**6.3 CONTRAINDICATIONS**

None known.

**6.4 PATIENT TARGET GROUP**

Applicable to the Depth Electrodes and the Placement Accessories.

The Depth Electrodes are used in pediatric and adult populations, regardless of gender.

**6.5 INTENDED USER**

The Anchor Bolt Placement Wrench is to be used by physicians/surgeons experienced in the placement of Ad-Tech’s Anchor Bolts.

**6.6 ENVIRONMENT AND FACILITIES**

The Ad-Tech Anchor Bolt Placement Wrenches are for surgical use only.

Placement and removal must be in a sterile surgical environment.

The surgeon must use the AD-TECH® LSB Anchor Bolt Wrench to place and remove the Anchor Bolts.

**6.7 EXPECTED CLINICAL BENEFIT**

Ad-Tech's Placement Wrenches are a placement tool to aid in the placement of Ad-Tech's Anchor Bolts, and should not be used on any patient whom the physician/surgeon considers at risk for infection or on whom the use cannot be performed safely.

**6.8 UNDESIRABLE SIDE EFFECTS**

No clinical risks with the use of the Anchor Bolt Placement Wrench.

**7 PERFORMANCE CHARACTERISTICS**

**7.1 DEVICE ATTRIBUTES**

The following Anchor Bolt Placement Wrenches are available.

Catalog #	Shaft Diameter (inches)	Length (inches)	Supplied Sterile Yes/No
LSB-PWL-2.4-6X	0.094 (2.4 mm)	6	Yes
LSB-PWL-2.4-10X	0.094 (2.4 mm)	10	Yes
LSB-PWL-2.4-6N	0.094 (2.4 mm)	6	No
LSB-PWL-2.4-10N	0.094 (2.4 mm)	10	No
LSB-PWL-2.4-11N	0.094 (2.4 mm)	11	No

**8 INSTRUCTIONS FOR USE**

	<i>Discontinue Anchor Bolt placement if you experience any instrument (e.g., the LSB Anchor Bolt Wrench, Cranial Drill, etc.) bending during the placement procedure.</i>
--	---

Anchor Bolts are placed at the target location using an appropriate stereotactic guidance system and supporting equipment.

Ad-Tech's Placement Wrenches are a placement tool to aid in the placement of Ad-Tech's Anchor Bolts. Refer to the Anchor Bolt IFU for details.

When used with the compatible stereotactic system, the AD-TECH® Drill Sleeve Guide supports the alignment of the Cranial Drill Bit and LSB Anchor Bolt Wrench. Refer to the Disposable Cranial Drill Bit Kit IFU for details.

Regarding Figure 2, the steps for placement and removal of the Anchor Bolt in the skull are as follows:


**Step 1: Placement of the Anchor Bolt.** Using the Cranial Drill Kit and a Drill Bit of the appropriate diameter, drill a hole in the patient's skull in the position and trajectory determined by pre-operative planning. Remove the Anchor Bolt Cap. Place the Anchor Bolt in the pre-drilled hole and use the LSB Anchor Bolt Wrench in a clockwise direction to screw it into place. The Anchor Bolt has a tapered thread and will naturally tighten when rotated. Do not over-tighten.

**Step 2: Removal of the Anchor Bolt.** Remove the Silicone Cap from the Anchor Bolt Cap and Depth Electrode. Untwist the Anchor Bolt Cap anticlockwise and remove it. Slowly withdraw the Depth Electrode from the Anchor Bolt. Use the LSB Anchor Bolt Wrench to unscrew the Anchor Bolt in an anticlockwise direction. Dispose of the Anchor Bolt as directed.

In the event of Anchor Bolt breakage:

1. Stop use immediately and monitor the patient.
2. Attempt to remove the remaining Anchor Bolt using the LSB Anchor Bolt Wrench.
3. Use other surgical equipment, e.g., a clamp, to grip and unscrew the retained portion.
4. Consider other surgical interventions as necessary.

Figure 2:

**9 PACKAGING, STERILIZATION, AND REUSE**

	<i>Clean and sterilize the Anchor Bolt Placement Wrench as directed.</i>  <i>Incorrect cleaning and sterilization of the device may result in contamination and severe patient harm.</i>
--	--

**10.1 STERILE/NON-STERILE Anchor Bolt Placement Wrenches**

Observe the following packaging, sterilization, and reuse instructions:

	<i>AD-TECH® sterilizes the Placement Wrenches using ethylene oxide.</i>
	<i>Do not use the Placement Wrench if the packaging is damaged or unintentionally opened before use.</i>  <i>Using a device from damaged packaging may result in contamination and severe patient harm.</i>
	<i>The Placement Wrench is packaged in a double sterile barrier to support surgical asepsis procedures.</i>

	<i>Non-sterile. Sterilize before use.</i>
--	---

After use, clean and sterilize the Placement Wrench as per AT10054-B AD-TECH® CLEANING AND STERILIZATION INSTRUCTIONS PLACEMENT TOOLS AND ACCESSORIES.

**10.2 REUSE**

	<i>Dispose of devices that show signs of wear or damage, irrespective of the number of reuse cycles.</i>
--	--

The Placement Wrench can be cleaned and re-sterilized for use, up to a total of twenty-five times.

*NOTE: If the device was provided sterile by the manufacturer, the product has already undergone 1 cycle.*

**10 STORAGE CONDITIONS**

- Necessitates medical or surgical intervention to preclude permanent damage or impairment.

	<p>Observe the recommended storage conditions. Failure to observe storage conditions may result in device failure.</p>
--	--

Store the Anchor Bolt Placement Wrench in their outer packaging, observing the handling and storage conditions detailed below.

	Fragile, handle with care.
	Keep dry
	<p>Temperature limit Upper = 30°C, Lower = 15°C</p>
	<p>Humidity limitation Upper = 60%, Lower = 20%</p>

**11 DISPOSAL**

	<p>Dispose of the Anchor Bolt Placement Wrench following the healthcare facility's biohazardous and pathological waste procedures.</p>
--	--

**12 INCIDENT REPORTING**

**13.1 EUROPEAN UNION**

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

**13.2 CANADA**

A medical device incident (MDI) that has occurred in relation to the device should be reported to the manufacturer and Health Canada.

A medical device incident (MDI) is defined as an incident related to the failure of a medical device, deterioration in its effectiveness, or inadequacy in its labeling or directions that led to the:

- Death of a patient, user or other person
- Serious deterioration in the state of health of a patient, user or other person
- Potential for death or serious deterioration in health of a patient, user or other person if the incident were to recur

**13.3 UNITED STATES**

Report a suspected medical device-related death to both the FDA and the manufacturer.

Report a medical device-related serious injury to the manufacturer or the FDA if the medical device manufacturer is unknown.

Serious injury means an injury or illness that is:

- Life-threatening;
- Results in permanent impairment of a body function or permanent damage to a body structure; or

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