Lightweight CABRIO Cable

Contents

en	Lightweight CABRIO Cable	1
cs	Kabel Lightweight CABRIO	5
da	Lightweight CABRIO-kabel.....	9
de	Lightweight CABRIO Kabel	13
el	Καλώδιο Lightweight CABRIO.....	17
es	Cable Lightweight CABRIO	21
fr	Câble CABRIO Lightweight	25
it	Cavo Lightweight CABRIO	29
nl	Lightweight CABRIO-kabel.....	33
pt	Cabo CABRIO Lightweight.....	37
tr	Lightweight CABRIO Kablosu.....	41

Lightweight CABRIO Cable


Symbol Definitions / Definice symbolů / **Symboldefinitioner / Symboldefinitionen /** Ορισμοί συμβόλων / **Definiciones de símbolos / Définitions des symboles / Definizioni dei simboli / Definities van symbolen / Definições dos símbolos / Sembol Tanımları**

Manufacturer / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Fabricant / Fabbricante / Fabrikant / Fabricante / Üretici

Use-by date / Datum použitelnosti / Udløbsdato / Haltbarkeitsdatum / Ημερομηνία λήξης / Fecha de caducidad / Date limite d'utilisation / Data di scadenza / Uiterste gebruiksdatum / Data de validade / Son kullanma tarihi


Authorized representative in the European Community / Zplnomocněný zástupce v Evropském společenství / Autoriseret repræsentant i EU / Bevollmächtigter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Gemachtigde vertegenwoordiger in de Europese Gemeenschap / Mandatário na Comunidade Europeia / Avrupa Topluluğu yetkili temsilcisi

Do not use if package is damaged / Nepoužívejte, pokud je obal poškozený / Må ikke anvendes, hvis pakken er beskadiget / Inhalt bei beschädigter Verpackung nicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά / No utilice si el envase está dañado / Ne pas utiliser si le conditionnement est endommagé / Non usare se la confezione è danneggiata / Niet gebruiken als de verpakking is beschadigd / Não utilizar se a embalagem estiver danificada / Ambalaj hasarlıysa kullanmayın

Batch Code / Kód šarže / Batchkode / Chargenbezeichnung / Κωδικός παρτίδας / Código de lote / Code de lot / Codice lotto / Partijcode / Código de lote / Parti Kodu

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DESCRIPTION

The Lightweight CABRIO Cable is designed to provide an easy method for making a secure connection between an Ad-Tech® electrode and standard clinic EEG monitoring and/or stimulation equipment.


Figure 1 - Lightweight CABRIO Cable

The Lightweight CABRIO Cable consists of a connection cable with an attached electrode connector block on one end and the other end terminating in a standard 1.5 mm safety socket DIN connector.

Depth Electrodes	Spencer Probe
	Macro-Micro
Subdural Electrodes	Grid
	Strip
	Dual-Sided Interhemispheric
	Multi-Strip and Split Grid

Table 1 - Compatible Ad-Tech® Electrodes

Design Variants

The most common variants involve the number of contacts. Ad-Tech® electrodes can have up to a maximum of sixteen (16) contacts per tail and the cables are often built to match. Variations to the cable style, cable length, and connector style are available to accommodate the procedural technique the physician defines as appropriate for the specific patient. These variants do not alter the indications or the intended use of the Lightweight CABRIO Cable.

The CABRIO Cable is provided NON-STERILE and must be cleaned and sterilized prior to use. The cable may be used up to five (5) times.

INDICATIONS AND INTENDED USE

The Lightweight CABRIO Cable is designed to provide a safe, secure, and easy method to connect Ad-Tech® electrodes with typical EEG monitoring and/or stimulation equipment.

INTENDED USER

This product should only be used by a healthcare professional trained in the use of the CABRIO connection system.

PATIENT TARGET GROUP

This device is used in patients diagnosed with intractable epilepsy who are undergoing brain mapping procedures to define the location of epileptogenic foci.

EXPECTED CLINICAL BENEFIT

The Lightweight CABRIO Cable provides a secure connection between the Ad-Tech® electrode and the EEG recording equipment.

CONTRAINDICATIONS

None known.

WARNINGS

- The CABRIO Cable is provided NON-STERILE and must be cleaned and sterilized prior to use. Refer to the cleaning and sterilization instructions in this document.
- Disconnect from monitoring equipment during cardiac defibrillation.

PRECAUTIONS

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- A direct pull or stress on the electrode may cause damage to the electrode or a loss of contact recordings. If the patient will become mobile, place and secure a strain relief using the EEG cable.
- **The CABRIO Cable may be used up to five (5) times.** Refer to the cleaning and sterilization instructions in this document for proper cleaning and sterilization methods.
- Do not immerse the CABRIO cable in fluids.
- The cable will only perform reliably when maintained in its original condition. A cable

that is or is suspected of being worn or damaged should not be used.

ADVERSE EVENTS

None known.

INSTRUCTIONS FOR USE

1. Inspect the cable for wear or damage. A cable that is or is suspected to be worn or damaged should be discarded.
2. Ensure that the clear connection block is in the open position. See Figure 1. The connection block should not be closed without the tail of an electrode inserted.


Figure 2. Block in the Open Position

3. Immediately upon completion of electrode placement, the tail end contacts need to be attached to the cable via the connector block.
4. Carefully insert the tail of the electrode into the connector block and ensure that the contacts on the electrode tail align with the pins in the connection block. See Figure 2.


Figure 3. Electrode Tail Inserted into Block

5. Close the connection block to secure the electrode into the block. See Figure 3. At no time should the tail end contacts be exposed or unprotected.

Figure 4. Block Closed to Complete Connection

6. It is important to place and secure a strain relief using the EEG cable to prevent placing stress or a direct pull on the electrode which may cause loss of contact recordings and/or damage to the electrode.
7. If the patient will become mobile, the electrode must remain connected to the cable.

CLEANING AND STERILIZATION

To help ensure patient safety, prevent cross contamination and provide maximum life from this product, proper reprocessing must be completed before each use.

Adhering to these validated reprocessing instructions, the cable may be used up to five (5) times.

It is the responsibility of the end user to ensure that the reprocessing is actually performed using appropriate equipment, materials, and personnel to achieve the desired result. Any deviation from these validated instructions should be evaluated for effectiveness and potential adverse consequences.

Manual Cleaning

The cable should be completely cleaned and free of visible soil, stains, and debris. If visible contaminants are present after cleaning, repeat the cleaning procedure. The cable can be cleaned using either Enzyme or Alkaline detergents or CaviWipes™.

Cleaning with Detergent

Enzymatic Detergent: Enzol® Enzymatic Detergent (Advanced Sterilization Products)

Alkaline Detergent: Polystica® 2X Concentrate Alkaline Detergent (Steris Healthcare)

NOTE: At no point should the cable be immersed.

1. Soak a clean, soft cloth in lukewarm tap water. Use the cloth to wipe the cable and remove any heavy soil that may be present.

2. Prepare cleaning solution per manufacturer's recommendations.
3. Soak a clean, soft cloth in the prepared detergent. Use the cloth to wipe the cable and remove any visible soil, stains, and debris. If necessary, a soft bristled brush can be used to remove soil from crevices. Dip brush in the prepared detergent prior to scrubbing the cable.
4. Re-wet the cloth and re-wipe the cable as needed to keep wet to mimic the manufacturer's suggested soak time.
5. Soak a clean, soft cloth in lukewarm tap water. Use the cloth to wipe the cable and remove the residual detergent. Re-wet the cloth and re-wipe the cable as needed.
6. Dry the cable with a clean, lint free cloth.
7. Visually inspect each cable for visible soil. If soil is still present on cable, repeat the cleaning procedure.

- STERRAD® 100S
- STERRAD® NX™
- STERRAD® 100NX™

Follow the directions provided in the STERRAD® System User's Guides for the sterilization procedure.

Instrument	Cycle
STERRAD® 100S	Short
STERRAD® NX™	Standard
STERRAD® 100NX™	Standard

Table 2 - STERRAD® Cycle Settings

The cables are also compatible with the following STERIS® V-PRO Systems:

Sterilizer/ Cycle	Lumen	Non - Lumen	Flexible	Fast Non - Lumen	Fast
V-PRO® maX	✓	✓	✓		
V-PRO® maX 2	✓	✓	✓	✓	
V-PRO® 60	✓	✓	✓		
V-PRO® s2	✓	✓	✓		✓

Table 3 - STERIS V-PRO® validated low temperature sterilization systems / cycle **STORAGE**

Store and use the CABRIO Cable in a dry, clean location maintained at typical hospital/office ambient temperature and humidity conditions.

DEVICE DISPOSAL

This device does not require any special disposal procedures, however, if the device has been contaminated with potentially infectious substances of human origin (such as bodily fluids), it should be handled according to hospital protocol for infectious medical waste.

INCIDENT REPORTING

Any serious incident that has occurred in relation to

Cleaning with CaviWipes™

CaviWipes™ (Metrex Infection Prevention)

CaviCide™ (Metrex Infection Prevention)

1. Use the CaviWipes™ to completely remove any visible soil, stains, and debris from the cable. Use as many wipes as needed. If necessary, a soft bristled brush sprayed with CaviCide™ can be used to remove soil from crevices.
2. Using a new CaviWipe™, keep the cable wet for the length of time listed on the product label. Re-wipe the cable with new wipes as needed to keep the cable wet.
3. Allow the cable to air dry.
4. Visually inspect each cable for visible soil. If soil is still present on cable, repeat the cleaning procedure.

Sterilization

The cable is compatible with the following STERRAD® Systems (Advanced Sterilization Products):

the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

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© 2023 AD-TECH
AD-TECH® Medical Instrument Corporation
E C REP Limited